CAST[®] Allergens

for

CAST[®] ELISA
Flow CAST[®]

Drug Allergens

A Commitment to Diagnostics

Introduction	4	BAG2-CERY	Erythromycin	28
Literature	15	BAG2-CRIF	Rifampicin	28
ACE inhibitors	16	Antiseptics		29
BAG2-CRAM	16	BAG2-CCHX	Chlorhexidine	29
Analgesics	17	Beta Blockers		30
BAG2-C51	17	BAG2-CBIS	Bisoprolol	30
BAG2-C52	17	BAG2-CMET	Metoprolol	30
BAG2-C53	17	Contrast Media		31
BAG2-C54	18	BAG2-CIOBI	Iobitridol	31
BAG2-C55	18	BAG2-CIODI	Iodixanol	31
BAG2-C56	18	BAG2-CIOHE	IoHexol	31
BAG2-C57	19	BAG2-CIOME	Iomeprol	32
BAG2-C58	19	BAG2-CIOPA	Iopamidol	32
BAG2-C59	19	BAG2-CIOPR	Iopromide	32
BAG2-CNAP	20	BAG2-CIOXA	Ioxaglate	33
BAG2-C114	20	Local Anesthetics		34
BAG2-CTRA	20	BAG2-CART	Articaine	34
Antibiotics	21	BAG2-CBUP	Bupivacaine	34
BAG2-C1	21	BAG2-CMEP	Mepivacaine	34
BAG2-C11	21	BAG2-CLID	Lidocaine	35
BAG2-C12	21	Myorelaxants		36
BAG2-C2	22	BAG2-CATR	Atracurium	36
BAG2-C3	22	BAG2-CCAT	Cis-Atracurium	36
BAG2-C31	22	BAG2-CMIV	Mivacurium	36
BAG2-C32	23	BAG2-CPAN	Pancuronium	37
BAG2-C33	23	BAG2-CPRO	Propofol	37
BAG2-C34	23	BAG2-CROC	Rocuronium	37
BAG2-C35	24	BAG2-CSUX	Suxamethonium	38
BAG2-C36	24	BAG2-CVEC	Vecuronium	38
BAG2-C61	24	BAG2-CNBX	Neuromuscular Blockers Mix	39
BAG2-C62	25	Proton Pump Inhibitors		40
BAG2-C75	25	BAG2-COME	Omeprazole	40
BAG2-C76	25	BAG2-CESO	Esomeprazole	40
BAG2-C81	26	BAG2-CPANT	Pantoprazole	40
BAG2-C82	26	BAG2-CLAN	Lansoprazole	41
BAG2-C203	26			
BAG2-C204	27			
BAG2-CCLA	27			
BAG2-CCLAX	27			
BAG2-CCLIN	28			

Introduction

Intended Use

BÜHLMANN CAST® Allergens are optimized for *in vitro* Basophil Activation using BÜHLMANN CAST® assays: CAST® ELISA (EK-CAST) and Flow CAST® (FK-CCR).

CAST® Allergens are not included in BÜHLMANN CAST® assays. They are provided upon request. Refer to the allergen list.

Instruction for Use

The following instructions are intended for the following groups of drug allergens and mixes of them (refer to TOC page 3).

- Adjuvants
- Analgesics
- Antibiotics
- Antiseptics
- Contrast Media
- Local Anesthetics
- Myorelaxants
- Proton Pump Inhibitors (PPI)

CAST® drug allergens are delivered lyophilized. Before usage in CAST® assays, they are to be reconstituted with test specific Stimulation Buffer: CAST® ELISA: B-CAST-STB and Flow CAST®: B-CCR-STB.

- Add 250 µL Stimulation Buffer into the vial and vortex until the lyophilized material has completely dissolved.
- For some allergens we recommend an additional dilution with Stimulation Buffer (refer to the paragraph „Additional dilutions“ in the allergen specifications).
- E.g. an „additional dilution“ of 1:5 is carried out by adding 160 µL Stimulation Buffer to 40 µL allergen solution.
- For stimulation of cells, use these dilutions according to the respective CAST® IFU.

Remark: Some individuals react positive either at high or at low allergen concentrations, whereas others react positive over a broad range of concentration.

Storage and Stability

Unopened CAST® Allergens are to be stored at the temperature specified on the label until the expiration date.

Important: Reconstituted or diluted allergens must not be stored and used again!

Optimal time period for stimulation tests

For optimal results, CAST® assays should be run between 3 and 12 weeks after an allergic reaction to the presumed allergen [Lit.8], although specific basophil reaction is often being persistent over a longer time period.

Blood samples must be drawn for CAST® assays before carrying out a skin test or *in vivo* provocation with the presumed allergen as *in vivo* exposure to an allergen can cause an activation of the leucocytes.

Limitations

A negative CAST® assay result for a specific allergen does not exclude potential (even serious) clinical symptoms. Individuals who developed „allergic“ reactions towards protein or drug allergens in the past and whose CAST® assay was negative, should thus be verified with additional methods like the *in vivo* provocation or skin test before a drug or allergen is administered.

Concentrations and cut-off

Concentrations and cut-off of these allergen groups were optimized regarding maximum specificity and sensitivity. Concentrations and cut-offs of individual allergens are reported in the specifications of the allergens (refer to TOC page 3).

For Flow CAST®, BÜHLMANN recommends the application of the stimulation index (SI) as an additional diagnostic criterion. The SI is defined as the ratio of allergen specific basophil activation and background activation.

Remarks:

- In order to apply allergen specific cut-offs, QC criteria published in the CAST® instructions must be fulfilled.
- Cut-offs should serve only as recommendations. Clinically defined cut-offs should be established by each laboratory or through clinical studies.

Recommendations how to use CAST® assays in the diagnosis of drug allergies

The diagnosis of drug allergies is a complex matter including heterogeneous symptoms and a diverse pathophysiology. Thus the following recommendations should be considered when using CAST® assays in order to diagnose drug allergies. These recommendations should be considered only as general guidelines.

Clinical pictures, justifying the application of CAST® assays

Direct "allergic" reactions *in vivo* after exposure to drugs or chemicals can either be IgE or non-IgE mediated (pseudo-allergic reactions) stimulation of effector cells (e.g. basophils, mast cells, and eosinophils). CAST® assays may serve as a model for the above mentioned reactions by which the individuals' leucocytes are stimulated *in vitro*. As for CAST® assays pure allergens and a defined but artificial buffer system are applied, the method can only serve as a simplified model of the real *in vivo* conditions.

CAST® assays were optimized for the detection of immediate type allergic reactions (Type I) and pseudo allergic reactions. CAST® assays are not applicable for the diagnosis of delayed type allergies (e.g. Type IV reactions). Thus, these assays should be applied only for immediate type allergic reactions, especially if one of the following symptoms is detected:

- Anaphylactic or anaphylactoid reactions
- Rhino conjunctivitis
- Asthma bronchiale
- Urticaria / Angioedema
- Gastrointestinal reactions

Clinical pictures that do not justify the application of CAST® assays

CAST® assays should not be used for the diagnosis of T-cell mediated reactions. Based on the pathogenetical mechanism it is not recommended using CAST® assays routinely for the following patient groups:

- Maculopapular and pustular exanthema
- Vasculitis
- Allergic contact dermatitis
- Bullous exanthemas such as Stevens-Johnson Syndrome and Lyell's disease

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1. Eberlein B, León Suárez I, Darsow U, Ruëff F, Behrendt H, Ring J: A new basophil activation test using CD63 and CCR3 in allergy to antibiotics. *Clin & Exp Allergy* 2010, 40: 411-8.
2. Garcia-Ortega P, Marin A, Quilez E, Rodriguez-Leor O: New tools for allergy diagnosis and clinical advice after immediate reactions to ICM. *Br J Radiol* 2010, 83: 632-33.
3. De Weck AL, Sanz ML, Gamboa PM, Jermann JM, Kowalski M, Medrala W, Sainte-Laudy J, Schneider MS, Weber JM, Wolanczyk-Medrala A: Nonsteroidal anti-inflammatory drug hypersensitivity syndrome. A multicenter study II. Basophil activation by nonsteroidal anti-inflammatory drugs and its impact on pathogenesis. *J Investig Allergol Clin Immunol* 2010, 20: 39-57.
4. Hausmann OV, Gentinetta T, Bridts CH, Ebo DG: The basophil activation test in immediate-type drug allergy. *Immunol Allergy Clin N Am* 2009, 29: 555-66.
5. Dewachter P, Nicaise-Roland P, Kalaboka S, Lefèvre J, Chollet-Martin S: Anaphylaxis to amidotrizoate proved by skin testing and flow cytometry-based basophil activation test. *Allergy* 2009, 64: 501-2.
6. De Weck AL, Sanz ML, Gamboa PM, Aberer W, Sturm G, Bilo MB, Montroni M, Blanca M, Torres MJ, Mayorga L, Campi P, Manfredi M, Drouet M, Staint-Laudy J, Romano A, Merk H, Weber JM, Jermann TM, members of ENDA (European Network for Drug Allergy): Diagnosis of immediate-type β -Lactam allergy in vitro by flow-cytometric basophil activation test and sulfidoleukotriene production: a multicenter study. *J Investig Allergol Clin Immunol* 2009, 19: 91-109.
7. Steiner U, Gentinetta T, Hausmann O: IgE-mediated anaphylaxis to intraarticular glucocorticoid preparations. *AJR Am J Roentgenol* 2009, 193:156-7.
8. Gómez E, Blanca-Lopez N, Torresz MJ, Requenaz G, Rondoz C, Canto G, Blanca M, Mayorga C: Immunoglobulin E-mediated immediate allergic reactions to dipyrone: value of basophil activation test in the identification of patients. *Clinical & Experimental Allergy* 2009, 39: 1217-1224.
9. Kvedariene V, Kamey S, Ryckwaert Y, Rongier M, Bousquet J, Demoly P, Arnoux B: Diagnosis of neuromuscular blocking agent hypersensitivity reactions using cytofluorimetric analysis of basophils. *Allergy* 2006, 61: 311-5.
10. Scherer K, Studer W, Figueiredo V, Bircher A: Anaphylaxis to isosulfan blue and cross-reactivity to patent blue V: Case report and review of the nomenclature of vital blue dyes. *Ann Allergy Asthma Immunol* 2006, 96: 497-0.
11. Erdmann SM, Ventocilla S, Moll-Slodowy S, Sauer I, Merk HF: Basophilenaktivierungstests in der Diagnostik von Arzneimittelreaktionen. *Hautarzt* 2005, 56: 38-43.
12. Sanz ML, Gamboa P, de Weck AL: A new combined test with flowcytometric basophil activation and determination of sulfidoleukotrienes is useful for in vitro diagnosis of Hypersensitivity to Aspirin and other nonsteroidal anti-inflammatory drugs. *Int Arch Allergy Immunol* 2005, 136: 58-72.
13. Sudherr PS, Haal JE, Read GF, Rowbottom AW, Williams PE: Flow cytometric investigation of peri-anaesthetic anaphylaxis using CD63 and CD203c. *Anaesthesia* 2005, 60: 251-56.
14. Gamboa P, Sanz ML, Caballero MR, Urrutia I, Antepará I, Esparza R, de Weck AL: The flowcytometric determination of basophil activation induced aspirin and other non-steroidal anti-inflammatory drugs (NSAIDs) is useful for in vitro diagnosis of the NSAID hypersensitivity syndrome. *Clin Exp Allergy* 2004, 34: 1448-57.
15. Torres MJ, Padial A, Mayorga C, Fernandez T, Sanchez-Sabate E, Cornejo-Garcia JA, Antunez C, Blanca M: The diagnostic interpretation of basophil activation test in immediate allergic reactions to betalactams. *Clin Exp Allergy* 2004, 34: 1768-75.
16. Gamboa PM, Sanz ML, Caballero MR, Antepará I, Urrutia I, Jáuregui I, González G, Diéguez I, De Weck AL: Use of CD63 expression as a marker of in vitro basophil activation and leukotrienes determination in metamizole allergic patients. *Allergy* 2003, 58: 312-7.
17. Sanz ML, Gamboa PM, De Weck AL: Clinical evaluation of in vitro tests in Diagnosis of immediate allergic reactions to β -lactam antibiotics. *ACI International* 2002, 14: 185-93.
18. Monneret G, Benoit Y, Debard AL, Gutowski MC, Topenot I, Bienvenu J: Monitoring of basophil activation using CD63 and CCR3 in allergy to muscle relaxant drugs. *Clin Immunol* 2002, 102: 192-9.
19. Sanz ML, Gamboa PM, Antepará I, Uasuf C, Vila L, Garcia-Aviles C, Chazot M, De Weck AL: Flow cytometric basophil activation test by detection of CD63 expression in patients with immediate-type reactions to betalactam antibiotics. *Clin Exp Allergy* 2002, 32: 277-86.
20. Monneret G, Benoit Y, Gutowski MC, Bienvenu J: Detection of basophil activation by flow cytometry in patients with allergy to muscle-relaxant drugs. *Anesthesiology* 2000, 92: 275-7.
21. Abuaf N, Rajoely B, Ghazouani E, Levy DA, Pecquet C, Chabane H, Leynadier F: Validation of a flow cytometric assay detecting in vitro basophil activation for the diagnosis of muscle relaxant allergy. *J Allergy Clin Immunol* 1999, 104: 411-8.

Order code	Name	$C_{23}H_{32}N_2O_5$
BAG2-CRAM	Ramipril	416.51

Trade Names Delix; Hypren; Lannapril; Ramicard; Ramiclaire; Ramilich; Triatec; Tritace; Vesdil

Content per vial	12.5 µg	Conc after reconst.	50 µg/mL
------------------	---------	---------------------	----------

Additional dilutions 1:5

Assay	CAST® ELISA	Flow CAST®
-------	--------------------	-------------------

Conc. in stimulation	10 µg/mL	11.4 µg/mL
----------------------	----------	------------

Cut-off	≥40 pg/mL	≥6% CD63; SI*≥2
---------	-----------	-----------------

Analgesics

Order code	Name
BAG2-C51	Lys-Aspirin

$C_{15}H_{22}N_2O_6$
MW 326.4

Trade Names Acetophen; Acetosal; Acylpyrin; Alcacyl; Asatard; Aspidol; Aspro; Caprin; Colfarit; Contrheuma; Delgesic; Duramax; Ecotrin; Endydol; Entrophen; Rhodine; Salacatin; Salcetogen; Saletin; Solpyron; Venopirin

Content per vial	625 µg	Conc after reconst.	2.50 mg/mL
------------------	--------	---------------------	------------

Additional dilutions 1:5

Assay	CAST® ELISA	Flow CAST®
-------	--------------------	-------------------

Conc. in stimulation	500 µg/mL	570 µg/mL
----------------------	-----------	-----------

Cut-off	≥90 pg/mL	≥5% CD63; SI*≥2
---------	-----------	-----------------

Order code	Name
BAG2-C52	Diclofenac

$C_{14}H_{10}Cl_2NNaO_2$
MW 318.1

Trade Names Allvoran; Benfofen; Cataflam; Diclo-Phlogont; Diclo-Puren; Diclorem; Dolobasan; Duravolten; Ecofenac; Effekton; Novapirina; Primofenac; Prophenatin; Rhumalgen; Valetan; Voldal; Voltaren; Voltarol; Xenid

Content per vial	6.25 µg	Conc after reconst.	25 µg/mL
------------------	---------	---------------------	----------

Additional dilutions 1:5

Assay	CAST® ELISA	Flow CAST®
-------	--------------------	-------------------

Conc. in stimulation	5 µg/mL	5.7 µg/mL
----------------------	---------	-----------

Cut-off	≥40 pg/mL	≥5% CD63; SI*≥2
---------	-----------	-----------------

Order code	Name
BAG2-C53	Ibuprofen

$C_{13}H_{18}O_2$
MW 206.3

Trade Names Alges; Algifor; Brufen; Dismenol; Dentigoa; Dolgirid; Dolgit; Dolocyl; Epobron; Femadon; Grefen; Ifufen; Ibumetin; Ibutad; Ibutop; Nurofen; Optifen; Saridon; Sonotryl; Spedifen; Tabalon; Urem

Content per vial	25 µg	Conc after reconst.	100 µg/mL
------------------	-------	---------------------	-----------

Additional dilutions 1:5

Assay	CAST® ELISA	Flow CAST®
-------	--------------------	-------------------

Conc. in stimulation	20 µg/mL	22.7 µg/mL
----------------------	----------	------------

Cut-off	≥50 pg/mL	≥5% CD63; SI*≥2
---------	-----------	-----------------

Analgesics

Order code	Name
BAG2-C54	Indomethacin

$C_{19}H_{16}ClNO_4$
MW 357.8

Trade Names Elmetacin; Indocin; Indometacin; Indophtal; Indo-TOP; Indocolir; Indomet; Inflam; Protaxon; Reusin; Rheublamin; Sigadoc; Sjoebo

Content per vial 2.5 µg Conc after reconst. 10 µg/mL

Additional dilutions 1:5

Assay **CAST® ELISA** **Flow CAST®**

Conc. in stimulation 2.0 µg/mL 2.27 µg/mL

Cut-off ≥40 pg/mL ≥5% CD63; SI*≥2

Order code	Name
BAG2-C55	Acetaminophen

$C_8H_9NO_2$
MW 151.2

Trade Names Acetalgin; Azur; Ben-u-ron; Boxacin; Boxonal; Buscopan; Captin; Cervopyrin; Coffalon; Dafalgan; Dolereduct; Dolevar; Enelfa; Epsimol; Exedrin; Febrifant; Fibrex; Gelonida; Panadol; Perfalgan; Treupel; Zolben

Content per vial 2.5 µg Conc after reconst. 10 µg/mL

Additional dilutions 1:5

Assay **CAST® ELISA** **Flow CAST®**

Conc. in stimulation 2.0 µg/mL 2.27 µg/mL

Cut-off ≥60 pg/mL ≥5% CD63; SI*≥2

Order code	Name
BAG2-C56	Mefenamic Acid

$C_{15}H_{15}NO_2$
MW 241.3

Trade Names Lysalgo; Mefenacid; Mefenamin; Mephadolol; Ponstan; Spiralgin

Content per vial 6.25 µg Conc after reconst. 25 µg/mL

Additional dilutions 1:5

Assay **CAST® ELISA** **Flow CAST®**

Conc. in stimulation 5 µg/mL 5.7 µg/mL

Cut-off ≥60 pg/mL ≥5% CD63; SI*≥2

Analgesics

Order code	Name
BAG2-C57	Phenylbutazone

$C_{19}H_{20}N_2O_2$
MW 308.4
sodium salt

Trade Names Ambene; Artrizin; Azolid; Bizolin; Butacote; Butadion; Butapirazol; Butadiona; Butatron; Butoz; Butazolidin; Buzon; Ecobutazone; Equipalazone; Exrheudon N; Fenibutol; Intrabutazone; Intrazone; Mepha-Butazon; Phenyzene

Content per vial	25 µg	Conc after reconst.	100 µg/mL
------------------	-------	---------------------	-----------

Additional dilutions 1:5

Assay	CAST® ELISA	Flow CAST®
Conc. in stimulation	20 µg/mL	22.7 µg/mL
Cut-off	≥80 pg/mL	≥5% CD63; SI*≥2

Order code	Name
BAG2-C58	Propyphenazone

$C_{14}H_{18}N_2O$
MW 230.3

Trade Names Amtalidin; Budiro; Causyth; Cistalgan; Eufibron; Isoprochin; Optalidon Saridon Spasmo

Content per vial	12.5 µg	Conc after reconst.	50 µg/mL
------------------	---------	---------------------	----------

Additional dilutions 1:5

Assay	CAST® ELISA	Flow CAST®
Conc. in stimulation	10 µg/mL	11.5 µg/mL
Cut-off	≥40 pg/mL	≥5% CD63; SI*≥2

Order code	Name
BAG2-C59	Dipyrone / Metamizole

$C_{13}H_{16}N_3NaO_4S \cdot H_2O$
MW 351.4
sodium salt

Trade Names Analgina; Analgin; Berlosin; Inalgon; Metalgin; Minalgin; Nopain; Novalgin; Novalgina; Novaminsulfon; Nolotil; Sulfonovin; Vetalgin

Content per vial	12.5 µg	Conc after reconst.	50 µg/mL
------------------	---------	---------------------	----------

Additional dilutions 1:5

Assay	CAST® ELISA	Flow CAST®
Conc. in stimulation	10 µg/mL	11.5 µg/mL
Cut-off	≥50 pg/mL	≥5% CD63; SI*≥2

Analgesics

Order code	Name
BAG2-CNAP	Naproxen

$C_{14}H_{13}NaO_3$
MW 252.2
sodium salt

Trade Names Alacetan; Aleve; Anaprox; Apranax; Dolormin; Dysmenalgit; Mobilat; Momendol; Naprobene; Proxen

Content per vial	100 µg	Conc after reconst.	400 µg/mL
Additional dilutions	1:10		
Assay	CAST® ELISA	Flow CAST®	
Conc. in stimulation	80 µg/mL	91 µg/mL	
Cut-off	≥200 pg/mL	≥5% CD63; SI*≥2	

Order code	Name
BAG2-C114	Salicylic acid

$C_7H_5NaO_3$
MW 160.1
sodium salt

Trade Names Metabolite of Aspirin. Alysine; Idocyl novum; Enterosalicyl; Magsalyl; Parbocyl-Rev; Sodium 2-hydroxybenzoate; 2-Hydroxybenzoic acid

Content per vial	250 µg	Conc after reconst.	1 mg/mL
Additional dilutions	1:5		
Assay	CAST® ELISA	Flow CAST®	
Conc. in stimulation	200 µg/mL	227 µg/mL	
Cut-off	≥120 pg/mL	≥5% CD63; SI*≥2	

Order code	Name
BAG2-CTRA	Tramadol

$C_{16}H_{25}NO_2$
MW 263.4
hydrochloride

Trade Names ConZip; Rybix ODT; Ryzolt; Ultram; Ultram ER

Content per vial	440 µg	Conc after reconst.	1.76 mg/mL
Additional dilutions	1:5		
Assay	CAST® ELISA	Flow CAST®	
Conc. in stimulation	353 µg/mL	400 µg/mL	
Cut-off	≥50 pg/mL	≥5% CD63; SI*≥2	

Antibiotics

Order code	Name
BAG2-C1	Penicillin G

 $C_{16}H_{17}N_2NaO_4S$

MW 356.4

sodium salt

Trade Names Benzylpenicillin Sodium; Sodium Penicillin II; Sodium Benzylpenicillin; Pencillinic Acid Sodium Salt; Penicillin: American Penicillin; Nalpen G; Novocillin; Penilaryn; Pen-A-Brasive; Veticillin

Content per vial	1 mg	Conc after reconst.	4 mg/mL
------------------	------	---------------------	---------

Additional dilutions	1:5
----------------------	-----

Assay	CAST® ELISA	Flow CAST®
-------	--------------------	-------------------

Conc. in stimulation	800 µg/mL	910 µg/mL
----------------------	-----------	-----------

Cut-off	≥50 pg/mL	≥5% CD63; SI*≥2
---------	-----------	-----------------

Order code	Name
BAG2-C11	Benzylpenicilloyl-Polylysine (PPL)

 $C_{28}H_{46}N_6O_8S$

MW 626.8

Trade Names Penicilloypolylysine; PPL; Cilligen; Pre-Pen; Testarpen. The BPO (Benzylpenicilloyl)-group is the major determinant of Typ I allergies against Benzylpenicillin and related antibiotics

Content per vial	12.5 µg	Conc after reconst.	50 µg/mL
------------------	---------	---------------------	----------

Additional dilutions	1:5
----------------------	-----

Assay	CAST® ELISA	Flow CAST®
-------	--------------------	-------------------

Conc. in stimulation	10 µg/mL	11.5 µg/mL
----------------------	----------	------------

Cut-off	≥50 pg/mL	≥5% CD63; SI*≥2
---------	-----------	-----------------

Order code	Name
BAG2-C12	Minor Determinant Mix (MDM)

Trade Names MDM; Mixture of Sodium-Benzylpenicillin and di-Sodium-Benzylpenicilloat

Content per vial	250 µg	Conc after reconst.	1 mg/mL
------------------	--------	---------------------	---------

Additional dilutions	1:5
----------------------	-----

Assay	CAST® ELISA	Flow CAST®
-------	--------------------	-------------------

Conc. in stimulation	200 µg/mL	225 µg/mL
----------------------	-----------	-----------

Cut-off	≥40 pg/mL	≥5% CD63; SI*≥2
---------	-----------	-----------------

Antibiotics

Order code	Name
BAG2-C2	Penicillin V

$C_{16}H_{18}N_2O_5SK$
MW 388.5
potassium salt

Trade Names Acipen-V; Apopen; Distaquaine V; Eskacillin V; Fenacillin; Fenospin; Meropenin; Orocillin; Oratren; Oспен; Pen-Oral; Pen-Vee; Phenopenicillin; P-Mega-Tablinen; Stabicillin; V-Gil; V-Gillin; Vebecillin

Content per vial	1 mg	Conc after reconst.	4 mg/mL
------------------	------	---------------------	---------

Additional dilutions 1:5

Assay	CAST® ELISA	Flow CAST®
-------	--------------------	-------------------

Conc. in stimulation 800 µg/mL 910 µg/mL

Cut-off ≥40 pg/mL ≥5% CD63; SI*≥2

Order code	Name
BAG2-C3	Cephalosporin C

$C_{16}H_{19}N_3O_8SZn$
MW 478.8
zinc salt

Trade Names None. Cephalosporin C is the basic structure of many antibiotics such as Cefazolin and Cefamandole.

Content per vial	25 µg	Conc after reconst.	100 µg/mL
------------------	-------	---------------------	-----------

Additional dilutions 1:5

Assay	CAST® ELISA	Flow CAST®
-------	--------------------	-------------------

Conc. in stimulation 20 µg/mL 22.5 µg/mL

Cut-off ≥40 pg/mL ≥5% CD63; SI*≥2

Order code	Name
BAG2-C31	Cefamandole

$C_{19}H_{17}N_6O_5S_2$
MW 462.5
formate sodium salt

Trade Names Bergacef; Cedol; Cefam; Cefiran; Cemando; Cemandil; Fado; Kefadol; Kefandol; Lampomandol; Mandokef; Mandol; Mandolsan; Neocefal; Pavecef

Content per vial	625 µg	Conc after reconst.	2.5 mg/mL
------------------	--------	---------------------	-----------

Additional dilutions 1:5

Assay	CAST® ELISA	Flow CAST®
-------	--------------------	-------------------

Conc. in stimulation 500 µg/mL 570 µg/mL

Cut-off ≥80 pg/mL ≥5% CD63; SI*≥2

Antibiotics

Order code	Name
BAG2-C32	Cefazoline

$C_{14}H_{14}N_6O_4S_3$
MW 454.5
sodium salt

Trade Names Acef; Ancef; Atirin; Biazolina; Bor-Cefazol; Cefacidal; Cefamedin; Cefamezin; Cefazil; Cefazina; Elzogram; Firmacef; Gramaxin; Kefzol; Lampocef; Liviclina; Totacef; Zolicef

Content per vial	625 µg	Conc after reconst.	2.5 mg/mL
------------------	--------	---------------------	-----------

Additional dilutions 1:5

Assay	CAST® ELISA	Flow CAST®
-------	--------------------	-------------------

Conc. in stimulation 500 µg/mL 570 µg/mL

Cut-off ≥80 pg/mL ≥5% CD63; SI*≥2

Order code	Name
BAG2-C33	Cefuroxime

$C_{16}H_{16}N_4O_8S$
MW 424.4
sodium salt

Trade Names Anaptivan; Biociclin; Biofurex; Bioxima; Cefamar; Cefoprom; Ceftin; Cefumax; Cefurax; Cefurex; Cefurin; Cepacine; Curocef; Curoxim; Duxima; Elobact; Kefurox; Novocef; Spectazole; Zinacef; Zinat; Zinnat

Content per vial	625 µg	Conc after reconst.	2.5 mg/mL
------------------	--------	---------------------	-----------

Additional dilutions 1:5

Assay	CAST® ELISA	Flow CAST®
-------	--------------------	-------------------

Conc. in stimulation 500 µg/mL 570 µg/mL

Cut-off ≥40 pg/mL ≥5% CD63; SI*≥2

Order code	Name
BAG2-C34	Cefaclor

$C_{15}H_{14}ClN_3O_4S$
MW 367.8

Trade Names Alenfral; Alfacet; Alfatil; Ceclor; Distaclor; Kefral; Keflor; Panacef; Panorol; Raniclор

Content per vial	0.75 mg	Conc after reconst.	3 mg/mL
------------------	---------	---------------------	---------

Additional dilutions 1:10

Assay	CAST® ELISA	Flow CAST®
-------	--------------------	-------------------

Conc. in stimulation 600 µg/mL 682 µg/mL

Cut-off ≥40 pg/mL ≥5% CD63; SI*≥2

Antibiotics

Order code	Name
BAG2-C35	Ceftriaxone

 $C_{18}H_{16}N_8Na_2O_7S_3$

MW 598.6

sodium salt

Trade Names Amplospec; Antibacin; Cefotrix; Rocefin; Rocephin; Tercefon

Content per vial	1 mg	Conc after reconst.	4 mg/mL
Additional dilutions	1:10		
Assay	CAST® ELISA	Flow CAST®	
Conc. in stimulation	720 µg/mL	818 µg/mL	
Cut-off	≥40 pg/mL	≥5% CD63; SI*≥2	

Order code	Name
BAG2-C36	Levofloxacin

 $C_{18}H_{20}FN_3O_4$

MW 361.4

Trade Names Cravit; Elequine; Floxel; Iquix. Leroxacin; Lesacin; Levaquin; Levokacin; Levox; Levoxacin; Mosardal; Nofaxin; Oftaquix; Quixin; Reskuin; Tavanic; Volequin

Content per vial	200 µg	Conc after reconst.	800 µg/mL
Additional dilutions	1:10		
Assay	CAST® ELISA	Flow CAST®	
Conc. in stimulation	160 µg/mL	182 µg/mL	
Cut-off	≥40 pg/mL	≥5% CD63; SI*≥2	

Order code	Name
BAG2-C61	Sulfamethoxazole

 $C_{10}H_{11}N_3O_3S$

MW 253.3

Trade Names Gantanol; Sinomin. Mixture with Trimethoprim: Bactrim; Berlocid; Cotrim; Cotrimox; Cotrimstada; Cotrimoxazol; Drylin; Escoprim; Eusaprim; Kepinol; Lagatrim; Nopil; Sigaprim; TMS forte

Content per vial	25 µg	Conc after reconst.	100 µg/mL
Additional dilutions	1:5		
Assay	CAST® ELISA	Flow CAST®	
Conc. in stimulation	20 µg/mL	22.5 µg/mL	
Cut-off	≥50 pg/mL	≥5% CD63; SI*≥2	

Antibiotics

Order code	Name	$C_{14}H_{18}N_4O_3$
BAG2-C62	Trimethoprim	MW 290.3

Trade Names Infectotrimet; Motrim; Triprim. Mixture with Sulfamethoxazole: Bactrim; Berlocid; Cotrim; Cotrimox; Cotrimstada; Cotrimoxazol; Drylin; Escoprim; Eusaprim; Kepinol; Lagatrim; Nopil; Sigaprim; TMS forte

Content per vial	25 µg	Conc after reconst.	100 µg/mL
------------------	-------	---------------------	-----------

Additional dilutions	1:5
----------------------	-----

Assay	CAST® ELISA	Flow CAST®
Conc. in stimulation	20 µg/mL	22.7 µg/mL
Cut-off	≥40 µg/mL	≥5% CD63; SI*≥2

Order code	Name	$C_{22}H_{24}N_2O_8$
BAG2-C75	Tetracycline	MW 444.4

Trade Names Actisite-Dental; Fluorex Plus; Imex; Mysteclin; Sumycin; Tefilin

Content per vial	25 µg	Conc after reconst.	100 µg/mL
------------------	-------	---------------------	-----------

Additional dilutions	1:5
----------------------	-----

Assay	CAST® ELISA	Flow CAST®
Conc. in stimulation	20 µg/mL	22.7 µg/mL
Cut-off	≥90 µg/mL	≥5% CD63; SI*≥2

Order code	Name	$C_{22}H_{24}N_2O_8$
BAG2-C76	Doxycycline	MW 444.4
		hydrochloride

Trade Names Aknefug Doxy; Antodox; Atridox; Doxakne; Dotur; Doxyclean; Doxyderma; DoxyHexal; Doxylag; Doxymono; Doxysol; Oraycea; Periostat; Rudocyclin; Supracyclin; Tasmacyclin; Vibramycin; Vibravenös; Zadorin

Content per vial	15.6 µg	Conc after reconst.	62.5 µg/mL
------------------	---------	---------------------	------------

Additional dilutions	1:5
----------------------	-----

Assay	CAST® ELISA	Flow CAST®
Conc. in stimulation	12.5 µg/mL	14.2 µg/mL
Cut-off	≥65 µg/mL	≥7% CD63; SI*≥2

Antibiotics

Order code	Name
BAG2-C81	Ciprofloxacin

$C_{17}H_{18}FN_3O_3$
MW 331.4

Trade Names Baycip; Ciflox; Giloxan; Ciprinol; Cipro; Ciprobay; Ciproxan; Ciproxin; Flociprin; Septicide; Velmonit

Content per vial	25 µg	Conc after reconst.	100 µg/mL
Additional dilutions	1:5		
Assay	CAST® ELISA	Flow CAST®	
Conc. in stimulation	20 µg/mL	22.5 µg/mL	
Cut-off	≥90 pg/mL	≥5% CD63; SI*≥2	

Order code	Name
BAG2-C82	Moxifloxacin

$C_{21}H_{24}FN_3O_4$
MW 401.4
hydrochloride

Trade Names Actira; Avelon; Avalox; Avelox; Octegra; Vigamox

Content per vial	31.25 µg	Conc after reconst.	125 µg/mL
Additional dilutions	1:5		
Assay	CAST® ELISA	Flow CAST®	
Conc. in stimulation	25 µg/mL	28.4 µg/mL	
Cut-off	≥40 pg/mL	≥7% CD63; SI*≥2	

Order code	Name
BAG2-C203	Ampicillin

$C_{16}H_{19}N_3O_4S$
MW 388.5

Trade Names Albipen; Amblosin; Ampicin; Ampilar; Ampitab; Amplital; Austrapen; Britacil; Cetampin; Doktacillin; Omnipen; Penbristol; Penbritin; Pentrex; Pentrexyl; Princillin; Principen; Tokiocillin; Totacillin; Ukapen; Vidopen

Content per vial	2.5 mg	Conc after reconst.	10 mg/mL
Additional dilutions	1:5		
Assay	CAST® ELISA	Flow CAST®	
Conc. in stimulation	2 mg/mL	2.27 mg/mL	
Cut-off	≥70 pg/mL	≥5% CD63; SI*≥2	

Antibiotics

Order code	Name
BAG2-C204	Amoxicillin

$C_{16}H_{19}N_3O_5S$
MW 365.4

Trade Names Alfamox; Amocilline; Amoxi; Amoxidal; Ardine; Betamox; Clamoxyl; Dura X; Flemoxin; Hiconcil; Ibiamax; Moxal; Neamoxyl; Ospamox; Pasetocin; Penimox; Raylina; Sawacillin; Silamax; Trimox; Ultimox; Zamocillin

Content per vial	625 µg	Conc after reconst.	2.5 mg/mL
------------------	--------	---------------------	-----------

Additional dilutions	1:5
----------------------	-----

Assay	CAST® ELISA	Flow CAST®
Conc. in stimulation	500 µg/mL	570 µg/mL
Cut-off	≥100 pg/mL	≥5% CD63; SI*≥2

Order code	Name
BAG2-CCLA	Clarithromycin

$C_{38}H_{69}NO_{13}$
MW 748.0

Trade Names Biaxin; Clathromycin; Cyllind; Klacid; Klaricid; Macladin; Mavid; Naxy; Veclam; Zeclar

Content per vial	25 µg	Conc after reconst.	100 µg/mL
------------------	-------	---------------------	-----------

Additional dilutions	1:5
----------------------	-----

Assay	CAST® ELISA	Flow CAST®
Conc. in stimulation	20 µg/mL	22.7 µg/mL
Cut-off	≥70 pg/mL	≥5% CD63; SI*≥2

Order code	Name
BAG2-CCLAX	Clavulanic Acid Amox

Trade Names Combined with amoxicillin (co-amoxiclav, trade names: Augmentin, Synulox, and others)

Content per vial	0.750 mg	Conc after reconst.	3.0 mg/mL
------------------	----------	---------------------	-----------

Additional dilutions	1:5
----------------------	-----

Assay	CAST® ELISA	Flow CAST®
Conc. in stimulation	600 µg/mL	682 µg/mL
Cut-off	≥50 pg/mL	≥5% CD63; SI ≥2

Antibiotics

Order code	Name
BAG2-CCLIN	Clindamycin

$C_{18}H_{33}ClN_2O_5S$

MW 425.0

hydrochloride

Trade Names Basocin; Clindac; Clin-Sanorania; Dalacin; Dentomycin; Jutaclin; Lanacine; Sobelin; Turimycin; Zindaclin

Content per vial 62.5 µg Conc after reconst. 250 µg/mL

Additional dilutions 1:5

Assay **CAST® ELISA** **Flow CAST®**

Conc. in stimulation 50 µg/mL 57 µg/mL

Cut-off ≥40 pg/mL ≥6% CD63; SI*≥2

Order code	Name
BAG2-CERY	Erythromycin

$C_{37}H_{67}NO_{13}$

MW 733.9

Trade Names Aknederm Ery; Aknefug-EL; Aknemycin; Aknilox; Erios; Eryaknen; Eryfluid; Erythrocin; Inderm; Infectomycin; Meromycin; Paediathrocin; Sanasepton; Stiemycine

Content per vial 125 µg Conc after reconst. 500 µg/mL

Additional dilutions 1:5

Assay **CAST® ELISA** **Flow CAST®**

Conc. in stimulation 100 µg/mL 114 µg/mL

Cut-off ≥40 pg/mL ≥6% CD63; SI*≥2

Order code	Name
BAG2-CRIF	Rifampicin

$C_{43}H_{58}N_4O_{12}$

MW 823.0

Trade Names Rifampicin; Rifaldazine; Rifamycin AMP; RIAMP; Abrifam; Eremfat; Rifa; Rifadine; Rifaldine Rifaldin; Rifaprodin; Rifoldin; Rimactane

Content per vial 25 µg Conc after reconst. 100 µg/mL

Additional dilutions 1:5

Assay **CAST® ELISA** **Flow CAST®**

Conc. in stimulation 20 µg/mL 22.7 µg/mL

Cut-off ≥100 pg/mL ≥5% CD63; SI*≥2

Antiseptics

Order code **Name**

BAG2-CCHX **Chlorhexidine**

$C_{22}H_{30}Cl_2N_{10} \cdot 2(C_6H_{12}O_7)$

MW 897.8

Digluconate

Trade Names Atoseptal; Chlorhexamed; Corsodyl; Dentohexin; Dosisseptine; Hexamedal; Hibidil; Hibiscrub; Lifo-Scrub; Periochip

Content per vial 1.25 µg Conc after reconst. 5 µg/mL

Additional dilutions 1:5

Assay **CAST® ELISA** **Flow CAST®**

Conc. in stimulation 1 µg/mL 1.14 µg/mL

Cut-off ≥100 pg/mL ≥5% CD63; SI*≥2

Beta Blockers

Order code	Name
BAG2-CBIS	Bisoprolol

$C_{18}H_{31}NO_4 \cdot 1/2C_4H_4O_4$
MW 383.5
Hemifumarate

Trade Names Bilol; Bisacardiol; Bisocor; Conor; Jutabis; Rivacor

Content per vial 12.5 µg Conc after reconst. 50 µg/mL

Additional dilutions 1:5

Assay **CAST® ELISA** **Flow CAST®**

Conc. in stimulation 10 µg/mL 11.4 µg/mL

Cut-off ≥50 pg/mL ≥6% CD63; SI*≥2

Order code	Name
BAG2-CMET	Metoprolol

$C_{15}H_{25}NO_3$
MW 267.4
Tartrate

Trade Names Toprol-XL; Lopressor; Metoprolol Tartrate

Content per vial 495 µg Conc after reconst. 1.98 mg/mL

Additional dilutions 1:5

Assay **CAST® ELISA** **Flow CAST®**

Conc. in stimulation 396 µg/mL 450 µg/mL

Cut-off ≥70 pg/mL ≥6% CD63; SI*≥2

Contrast Media

Order code	Name
BAG2-CIOME	Iomeprol

$C_{17}H_{22}I_3N_3O_8$
MW 777.1

Trade Names Iomeron

Content per vial 1.25 mg Conc after reconst. 5 mg/mL

Additional dilutions

Assay	CAST® ELISA	Flow CAST®
Conc. in stimulation	1 mg/mL	1.1 mg/mL
Cut-off	≥50 pg/mL	≥5% CD63; SI*≥2

Order code	Name
BAG2-CIOPA	Iopamidol

$C_{17}H_{22}I_3N_3O_8$
MW 777.1

Trade Names Iopamiro; Isovue; Iopamiron; Niopam; Scanlux

Content per vial 1.25 mg Conc after reconst. 5 mg/mL

Additional dilutions

Assay	CAST® ELISA	Flow CAST®
Conc. in stimulation	1 mg/mL	1.1 mg/mL
Cut-off	≥160 pg/mL	≥5% CD63; SI*≥2

Order code	Name
BAG2-CIOPR	Iopromide

$C_{18}H_{24}I_3N_3O_8$
MW 791.1

Trade Names Ultravist

Content per vial 1.25 mg Conc after reconst. 5 mg/mL

Additional dilutions

Assay	CAST® ELISA	Flow CAST®
Conc. in stimulation	1 mg/mL	1.1 mg/mL
Cut-off	≥200 pg/mL	≥5% CD63; SI*≥2

Contrast Media

Order code	Name
BAG2-CIOXA	Ioxaglate

$C_{24}H_{21}I_6N_5O_8$
MW 1268.9

Trade Names Hexabrix; Ioxaglic acid

Content per vial	1.25 mg	Conc after reconst.	5 mg/mL
------------------	---------	---------------------	---------

Additional dilutions

Assay	CAST® ELISA	Flow CAST®
Conc. in stimulation	1 mg/mL	1.1 mg/mL
Cut-off	≥50 pg/mL	≥5% CD63; SI*≥2

Local Anesthetics

Order code	Name
BAG2-CART	Articaine

$C_{13}H_{20}N_2O_3S \times HCl$
MW 320.4
hydrochloride

Trade Names Astracaine; Rudocain; Septanest; Septocaine; Ubistesin; Ultracaine

Content per vial	100 µg	Conc after reconst.	400 µg/mL
Additional dilutions	1:5		
Assay	CAST® ELISA	Flow CAST®	
Conc. in stimulation	80 µg/mL	91 µg/mL	
Cut-off	≥40 pg/mL	≥5% CD63; SI*≥2	

Order code	Name
BAG2-CBUP	Bupivacaine

$C_{18}H_{28}N_2O \times HCl$
MW 324.9
hydrochloride

Trade Names Carbostesin; Marcain; Marcaine; Sensorcaine; Vivacaine

Content per vial	100 µg	Conc after reconst.	400 µg/mL
Additional dilutions	1:10		
Assay	CAST® ELISA	Flow CAST®	
Conc. in stimulation	80 µg/mL	91 µg/mL	
Cut-off	≥40 pg/mL	≥5% CD63; SI*≥2	

Order code	Name
BAG2-CMEP	Mepivacaine

$C_{15}H_{22}N_2O \times HCl$
MW 282.8
hydrochloride

Trade Names Carbocaina; Isogaine; Meaverin; Mecain; Mepibil; Mepicaton; Mepident; Mepidont; Mepiforan; Mepihexal; Mepimynol; Mepivacaina; Mepivamol; Mepivastesin; Scandicain; Scandicaine; Tevacaine

Content per vial	100 µg	Conc after reconst.	400 µg/mL
Additional dilutions	1:10		
Assay	CAST® ELISA	Flow CAST®	
Conc. in stimulation	80 µg/mL	91 µg/mL	
Cut-off	≥40 pg/mL	≥5% CD63; SI*≥2	

Local Anesthetics

Order code	Name
BAG2-CLID	Lidocaine

$C_{14}H_{20}N_2O Na_2$

MW 278.3

sodium salt

Trade Names Cuivasil; Duncaine; Leostesin; Lidesthesin; Lidostesin; Lignavet; Odontalg; Rucaina; Sedagul;
Xylocaine; Xylocard; Xylocitin; Xyloneural; Xylotox;

Content per vial	62.5 µg	Conc after reconst.	250 µg/mL
------------------	---------	---------------------	-----------

Additional dilutions 1:5

Assay	CAST® ELISA	Flow CAST®
-------	--------------------	-------------------

Conc. in stimulation	50 µg/mL	57 µg/mL
----------------------	----------	----------

Cut-off	≥40 pg/mL	≥5% CD63; SI*≥2
---------	-----------	-----------------

Myorelaxants

Order code	Name
BAG2-CATR	Atracurium

 $C_{65}H_{82}N_2O_{18}S_2$

MW 929.1

besylate salt

Trade Names Acuron; Aratan; Atratan; Atrelax; Atrican; Tracrium; Trarium; Wellcome 33-A-74

Content per vial	1.25 µg	Conc after reconst.	5 µg/mL
Additional dilutions	1:5		
Assay	CAST® ELISA	Flow CAST®	
Conc. in stimulation	1 µg/mL	1.2 µg/mL	
Cut-off	≥50 pg/mL	≥5% CD63; SI*≥2	

Order code	Name
BAG2-CCAT	Cis-Atracurium

 $C_{53}H_{72}N_2O_{12} \times 2C_6H_5SO_3$

MW 929.1

Besylate

Trade Names Nimbex

Content per vial	100 µg	Conc after reconst.	400 µg/mL
Additional dilutions	1:5		
Assay	CAST® ELISA	Flow CAST®	
Conc. in stimulation	80 µg/mL	91 µg/mL	
Cut-off	≥40 pg/mL	≥5% CD63; SI*≥2	

Order code	Name
BAG2-CMIV	Mivacurium

 $C_{58}H_{80}Cl_2N_2O_{14}$

MW 1100.2

Trade Names Mivacron

Content per vial	250 µg	Conc after reconst.	1 mg/mL
Additional dilutions	1:5		
Assay	CAST® ELISA	Flow CAST®	
Conc. in stimulation	200 µg/mL	227 µg/mL	
Cut-off	≥40 pg/mL	≥5% CD63; SI*≥2	

Myorelaxants

Order code **Name**
BAG2-CSUX **Suxamethonium**

$C_{14}H_{30}Cl_2N_2O_4$
MW361

chloride salt

Trade Names Lysthenon; Midarine; Pantolax; Succinolin

Content per vial 2.50 mg Conc after reconst. 10 mg/mL

Additional dilutions 1:5

Assay **CAST® ELISA** **Flow CAST®**

Conc. in stimulation 2 mg/mL 2.27 mg/mL

Cut-off ≥40 pg/mL ≥5% CD63; SI*≥2

Order code **Name**
BAG2-CVEC **Vecuronium**

$C_{34}H_{57}BrN_2O_4$
MW 637.7

bromid salt

Trade Names Musculax; Norcuron

Content per vial 62.5 µg Conc after reconst. 250 µg/mL

Additional dilutions 1:5

Assay **CAST® ELISA** **Flow CAST®**

Conc. in stimulation 50 µg/mL 51 µg/mL

Cut-off ≥40 pg/mL ≥5% CD63; SI*≥2

Myorelaxants

Order code	Name
BAG2-CN BX	Neuromuscular Blockers Mix

Ingredients	Content per vial µg	Con. after reconst. µg/mL
Suxamethonium	795	3180
Recuronium	43	171
Vecuronium	96	383
Pancuronium	95	381
Atracurium	12	50
Mivacurium	10	41
Cisatracurium	20	82

Content per vial	1.072 mg	Conc after reconst.	4.3 mg/mL
------------------	----------	---------------------	-----------

Additional dilutions 1:5

Assay	CAST® ELISA	Flow CAST®
-------	--------------------	-------------------

Conc. in stimulation	0.86 mg/mL	0.97 mg/mL
----------------------	------------	------------

Cut-off	≥200 pg/mL	≥5% CD63; SI*≥2
---------	------------	-----------------

Proton Pump Inhibitors

Order code **Name**

BAG2-CLAN **Lansoprazole**

$C_{16}H_{14}F_2N_3O_2S$
MW 396.4

Trade Names Agopton; Lanciprol; Lansobene; Lansoloc; Lansoprazol; Lansox; Lanton; Lanzo; Lanzoprial; Lanzor; Lanzostad; Limpidex; Monolitum; Prevacid; Zoton; Zomel

Content per vial 50 µg Conc after reconst. 200 µg/mL

Additional dilutions 1:5

Assay **CAST® ELISA** **Flow CAST®**

Conc. in stimulation 40 µg/mL 45.5 µg/mL

Cut-off ≥200 pg/mL ≥5% CD63; SI*≥2

BÜHLMANN Laboratories AG
Baselstrasse 55
4124 Schönenbuch
Switzerland

Phone +41 61 487 12 12
Fax orders +41 61 487 12 99
info@buhlmannlabs.ch
www.buhlmannlabs.ch

Ordering Code:

Flow CAST®	FK-CCR	100 tests
CAST® ELISA	EK-CAST	96 wells
	EK-CAST5	480 wells

CAST® is a registered trademark of
BÜHLMANN Laboratories AG