# Direct Saliva MELATONIN ELISA

EK-DSM 96 tests

Release Date: 2019-01-14  
Version A1

## ENGLISH

### INTENDED USE

The BÜHLMANN Direct Saliva Melatonin ELISA (EK-DSM) is intended for highly sensitive, quantitative determination of melatonin in human saliva (1-4).

### PRINCIPLE OF THE ASSAY

The BÜHLMANN Direct Saliva Melatonin ELISA is a competitive immunoassay using a capture antibody (Ab) technique. The polyclonal Kennaway G280 anti-melatonin antibody (5, 6) has been coated onto the microtiter plate, provided in the kit. After the first 16-20 hour overnight incubation, melatonin present in the pre-treated saliva and controls as well as in the calibrators, compete with biotinylated melatonin during a 3 hour incubation for the binding sites of this highly specific antibody. After washing, the enzyme label streptavidin conjugated to horseradish peroxidase (HRP) is added, which binds to the melatonin-biotin-antibody complexes captured on the coated wells during a 60 minute incubation step. Unbound enzyme label is then removed by a washing step and TMB substrate (tetramethylbenzidine) is added to the wells. In a further 30 minute incubation step, a chromophore is formed in inverse proportion to the amount of melatonin present in the sample. The color turns from blue to yellow after the addition of an acidic stop solution and can be measured at 450 nm.

### REAGENTS SUPPLIED AND PREPARATION

Reagents	Quantity	Code	Reconstitution
<b>Pretreatment Solution</b>	1 vial 5 mL	B-EKDSM-PRS	Ready to use <b>corrosive agent</b>
<b>Neutralizing Solution</b>	1 vial 5 mL	B-EKDSM-NS	Ready to use <b>irritant</b>
<b>Microtiter Plate</b> precoated with G280 anti-melatonin Ab	12x8 wells	B-EKDSM-MP	Wash 2x before use
<b>Plate Sealer</b>	3 pieces	-	Ready to use
<b>Wash Buffer Concentrate (10x)</b> with preservatives	1 bottle 100 mL	B-EKDSM-WB	Dilute with 900 mL deionized water
<b>Blanking Reagent<sup>1)</sup></b> lyophilized	1 vial lyoph	B-EKDSM-BR	Reconstitute in 1mL Incubation Buffer
<b>Incubation Buffer</b> (Zero Calibrator) melatonin-free buffer	1 vial 12 mL	B-EKDSM-IB	Ready to use
<b>Calibrators<sup>2)</sup></b> lyophilized; do not pretreat	5 vials lyoph	B-EKDSM-CASET	Reconstitute each with 1 mL Incubation Buffer
<b>Control low / high<sup>3)</sup></b> for pretreatment see page 3	2 vials lyoph	B-EKDSM-CONSET	Reconstitute each with 1 mL Incubation Buffer
<b>Biotin Conjugate</b>	1 vial 5.5 mL	B-EKDSM-BC	Ready to use
<b>Enzyme Label</b> Streptavidin conjugated to HRP	1 vial 11 mL	B-EKDSM-EL	Ready to use
<b>TMB Substrate</b> buffered with citrate	1 vial 11 mL	B-TMB	Ready to use

Reagents	Quantity	Code	Reconstitution
<b>Stop Solution</b> 0.25 M sulfuric acid (H <sub>2</sub> SO <sub>4</sub> )	1 vial 11 mL	B-ST5	Ready to use <b>corrosive agent</b>

Table 1

- The Blanking reagent contains a saturated melatonin solution. Prevent any contamination of other kit reagents.
- The Calibrators A, B, C, D and E contain the following melatonin concentration: 0.48, 1.2, 3.2, 8.0 and 20 pg/mL which are corrected for the 20% sample dilution during pretreatment and therefore, labeled with 0.6, 1.5, 4.0, 10, and 25 pg/mL of melatonin, respectively.
- Lot specific amount of melatonin see data sheet added to the kit.

### STORAGE AND SHELF LIFE OF REAGENTS

Sealed / Unopened Reagents	
Store at 2-8°C until expiration date. Do not use past expiration date.	
Opened / Reconstituted Reagents	
Microtiter Plate	Return unused strips immediately to the plastic pouch containing the desiccant pack and reseal along the entire edge of zip-seal. Store for up to 2 months at 2-8°C
Pretreatment Reagent	Store at 2-8°C until expiration date printed on the labels.
Neutralizing Solution	
Incubation Buffer	
Wash Buffer diluted	Store at 2-8°C up to 6 months
Blanking Reagent	Stable at 2-8°C up to 4 months.
Calibrators	
Controls	
Biotin Conjugate	Store at 2-8 °C until expiration date printed on the labels.
Enzyme Label	
TMB Substrate	
Stop Solution	

Table 2

### MATERIALS REQUIRED BUT NOT PROVIDED

- Precision pipettes with disposable tips: 5 µL, 50 µL, 100 µL and 1000 µL pipettes. Repeater or multichannel pipette for 50 µL and 100 µL.
- Disposable polystyrene or polypropylene tubes for the preparation of sample dilutions.
- 1000 mL cylinder for the dilution of the wash buffer.
- Microtiter plate washer or squeeze bottle for wash buffer.
- Blotting paper.
- Refrigerator
- Microtiter plate orbital shaker.
- Microtiter plate reader for measurement of absorbance at 450 nm.
- BÜHLMANN Saliva Collection Devices, B-SLEEP-CHECK16 or B-SVC (50).

### PRECAUTIONS

#### Safety precautions

- This test is for in vitro use only, and must be handled by qualified personnel, in accordance with good laboratory practices (GLP).
- Pretreatment/ neutralizing solution, substrate and stop solution: The pretreatment solution (B-EKDSM-PRS)

contains sodium hydroxide (NaOH) and the neutralizing solution (B-EKDSM-NS) contains hydrochloric acid (HCl). Substrate and stop solution: The substrate TMB (B-TMB) contains Tetramethylbenzidine. The stop solution (B-STs) contains sulfuric acid (0.25 M). Each of those reagents is irritant to eyes, skin and mucous membranes. Avoid contact with eyes, skin and clothes. Wear suitable protective clothing, gloves and eye protection. After contact with eyes or skin, wash immediately with plenty of water.

- Unused solution should be disposed of according to local state and federal regulations.

### **Technical precautions**

#### **Kit components**

- Read the instructions carefully before carrying out the test. Test performance will be adversely affected if reagents are incorrectly diluted, modified or stored under conditions other than those as detailed in this instruction for use.
- Residues in the microtiter plate wells result from the production process. They are removed in the washing step (assay procedure step 2) and do not affect the results.
- Components must not be used after the expiry date printed on the labels.
- Do not mix different lots of reagents.
- Every effort should be made to ensure that no cross contamination occurs between reagents, samples or between wells.
- Microwells cannot be re-used.
- Let the reagents adjust to reach room temperature. Mix the reagents well (vortex) before use.

#### **Assay procedure**

- The blanking reagent contains a saturated melatonin solution. Avoid any contamination of other reagents of this kit. Change disposable tips after each pipetting step.
- The assay procedure has been optimized for Sleep Check application. Therefore blank reagent and calibrators are assayed in duplicates, whereas controls and patient samples are measured in single determinations. This approach allows you to test 16 individual profiles (with 5 points each) per microtiter plate. For applications other than Sleep Check duplicate determinations are recommended.

---

### **SPECIMEN SHIPMENT AND STORAGE**

Shipment: The collected saliva samples must be shipped to the laboratory within two days. Collected saliva samples must be kept in the fridge at 2-8°C.

Samples should not be sent on Fridays, Saturdays or the day before holiday.

Storage: The saliva samples absorbed in the cotton swab may be stored in the saliva collection device for up to 7 days at 2-8°C. If not assayed within one week after collection, samples should be frozen and may be stored for at least 6 months at ≤ -20°C. Repeated freeze-thaw cycles should be avoided.

---

### **SPECIMEN COLLECTION**

Collect saliva using the BÜHLMANN Saliva Collection Devices. The devices can absorb up to 3 mL of saliva. The procedure calls for 0.2 mL of saliva.

- Do not use cotton swabs containing citric acid.
- Do not stimulate saliva flow by chewing gums or eating lemons.
- Patients should perform the collection on an evening without sporting activities and any intense efforts.
- When collecting saliva at night, a dim flash light or a ≤100 lux yellow light should be used in order to avoid a possible light influence on the melatonin profile.
- Nothing should be eaten during the collection time. The last meal must be taken at least 30 minutes before starting the collection. Bananas and chocolate should not be eaten during the entire day before the collection. Rinse the mouth with water 15 minutes before the collection.
- Drinks containing artificial colorants, caffeine (coffee, black or green tea, iced tea, cola) or alcohol are to be avoided on the evening of the collection.
- Patients should avoid brushing their teeth, with or without toothpaste, during sampling periods. It is likely that patients with gingivitis will contaminate the saliva with plasma or even whole blood leading to unknown consequences.
- On the collection day, if possible, no aspirin and medicines that contain ibuprofen (Brufen®, Algifor® Dismenol®, Dolocyl®, Ecoprofen®) should be taken. If your sleep or sleep-wake rhythm is treated with melatonin, this must be discontinued at least one week before the collection.

---

### **SAMPLE PRETREATMENT (LABORATORY)**

#### **Sample recovery from saliva collection devices**

Centrifuge the collection devices sent by the patient for around 5 min at 3000 rpm (~1500x g). Discard the suspended insert with the swab and store the tube at 2-8°C or -20°C.

#### **Pre-treatment of saliva samples and controls**

- Pipet 200 µL of controls and saliva samples, respectively, into correspondingly marked polypropylene tubes.
- Add 25 µL of pretreatment solution to each tube using a multichannel or repeater pipette.
- Vortex for 5 seconds and leave the tubes for 10 minutes at 18-28 °C.
- Add 25 µL of neutralizing solution to each tube using a multichannel or repeater pipette. Vortex for 5 seconds.
- Centrifuge the pre-treated samples for 5 minutes at 10'000 rpm. Proceed to the ELISA procedure.

---

### **ASSAY PROCEDURE**

1. Use a plate with enough 8-well strips to test the desired number of blanks, calibrators, controls and samples. Remove excess strips from the holder and re-seal them

in the foil pouch together with the two desiccant bags **without delay**. Store refrigerated.

2. Wash the coated strips twice using at least 300  $\mu\text{L}$  of Wash Buffer per well. Empty the wells and strike the plate firmly onto blotting paper.
  - 3a. Pipet 100  $\mu\text{L}$  of blanking reagent (blank) in duplicate into wells A1+A2.
  - 3b. Pipet 100  $\mu\text{L}$  of incubation buffer (zero calibrator) in duplicate into wells B1+B2.
  - 3c. Pipet 100  $\mu\text{L}$  of calibrator A in duplicate into wells C1+C2  
Pipet 100  $\mu\text{L}$  of calibrator B in duplicate into wells D1+D2  
Pipet 100  $\mu\text{L}$  of calibrator C in duplicate into wells E1+E2  
Pipet 100  $\mu\text{L}$  of calibrator D in duplicate into wells F1+F2  
Pipet 100  $\mu\text{L}$  of calibrator E in duplicate into wells G1+G2
  - 3d. Pipet 100  $\mu\text{L}$  of pretreated low control into well H1  
Pipet 100  $\mu\text{L}$  of pretreated high control into well H2
  - 3e. Pipet 100  $\mu\text{L}$  of each pretreated sample (single) into the subsequent wells.
  4. Cover the plate with a plate sealer and incubate for 16-20 hours at 2-8  $^{\circ}\text{C}$ .
  5. Remove and discard plate sealer. Add 50  $\mu\text{L}$  of biotin conjugate (blue solution) to each well. Cover the plate with a plate sealer and place it for 1 min on a plate orbital shaker set at 600 rpm.
  6. Incubate for 3 hours ( $\pm 5$  minutes) at 2-8  $^{\circ}\text{C}$ .
  7. Remove and discard the plate sealer. Aspirate or invert the plate to empty the solution from each well and wash four times using at least 300  $\mu\text{L}$  of wash buffer per well. Empty the wells and strike the plate firmly onto blotting paper.
  8. Add 100  $\mu\text{L}$  of enzyme label (yellow solution) to all wells.
  9. Cover the plate with a new plate sealer, place the plate on a plate orbital shaker set at 600 rpm and incubate for 60 minutes ( $\pm 5$  minutes) at 18-28  $^{\circ}\text{C}$ .
  10. Remove and discard the plate sealer. Aspirate or invert the plate to empty the solution from each well and wash four times using at least 300  $\mu\text{L}$  of wash buffer per well. Empty the wells and strike the plate firmly onto blotting paper.
- Important:** allow the TMB substrate to equilibrate to 18-28 $^{\circ}\text{C}$  prior to use.
11. Add 100  $\mu\text{L}$  of TMB substrate to all wells.
  12. Cover the plate, place it on a plate orbital shaker set at 600 rpm, protect the plate from direct light and incubate for 30 $\pm$ 5 minutes at 18-28 $^{\circ}\text{C}$ .
  13. Add 100  $\mu\text{L}$  of stop solution to all wells. Remove air bubbles by pricking them with a pipette tip. Proceed to step 14 within 30 minutes.
  14. Read the absorbance at 450 nm in a microtiter plate reader.

## QUALITY CONTROL

A thorough understanding of this instruction for use is necessary for the successful use of the product. Reliable results will be obtained only by precise laboratory techniques (current GLP guidelines) and accurately following this instruction for use.

Since there are no controls for saliva melatonin commercially available, we recommend using saliva pools containing different levels of melatonin for internal quality

controls. The reproducibility of standard curve parameters and control values should be within established limits of laboratory acceptability. The confidence limits for the controls are lot-specific and printed on the additional QC data sheet.

If the performance of the assay does not meet the established limits and repetition has excluded errors in technique, check the following issues: i) pipeting, temperature controlling and timing devices ii) ELISA reader settings iii) expiration dates of reagents iv) storage and incubation conditions v) TMB substrate solution should be colorless and vi) purity of water.

## STANDARDIZATION

Direct Saliva Melatonin ELISA is calibrated with UV/VIS:  $\epsilon_{278} = 6300 \text{ M}^{-1}\text{cm}^{-1}$  in methanol.

## RESULTS

### Standard Curve

Record the absorbance at 450 nm for each calibrator, incubation buffer and blank well. Average the duplicate values, subtract the average of the blank wells and record averages (=corrected average absorbance). Calculate the binding (B) of each pair of calibrator wells as a percent of incubation buffer ( $B_0$ ), with the blank-corrected absorbance of the incubation buffer taken as 100 %.

$$B/B_0 (\%) = \text{percent bound} = \frac{\text{net absorbance}}{\text{net absorbance of incubation buffer}} \times 100$$

Plot the percent bound (vertical axis) versus the concentration of melatonin in pg/mL (horizontal axis) using a lin/log graph paper. Draw the best fitting curve or calculate the standard curve using a four parameter algorithm.

### Samples and controls

- Record the absorbance at 450 nm for each sample, each control and well. Subtract the average of the blank wells and record the absorbance (=corrected average absorbance). Calculate, as described above, the binding of each pair of sample wells as a percent of incubation buffer ( $B_0$ ), with the blank-corrected absorbance of incubation buffer taken as 100%.
- Locate the  $B/B_0$  value of the samples on the vertical axis, draw a horizontal line intersecting the standard curve and read the melatonin concentration (pg/mL) from the horizontal axis.

See table 3 and figure 1 for examples of results and standard curves. *These results and standard curves are for demonstration purposes only. A standard curve must be generated for each set of samples to be assayed.*

## LIMITATIONS

Melatonin results should be interpreted in conjunction with information available from clinical assessment of the patient and other diagnostic procedures.

## PERFORMANCE CHARACTERISTICS

**Intra-Assay Precision: 12.6 %.** The intra-assay precision was calculated from the results of four different saliva samples within the standard range, measured 10 times in

duplicate in a single run. The results are presented in table 4.

**Inter-Assay Precision: 22.9 %.** The inter-assay precision was calculated from the results of 17 independent runs with 5 samples within the standard range. The results are presented in table 5.

**Detection Limit (LoB): 0.5 pg/mL.** 32 wells of incubation buffer (zero calibrator) were assayed in two independent runs. The minimum detectable concentration in 0.1 mL of incubation buffer was calculated by subtracting two standard deviations of averaged refer values from the OD of zero calibrator and intersecting the value with the standard curve obtained in the same run.

**Detection Limit (Limit of Quantification - LoQ): 1.6 – 20.5 pg/mL.** The limit of quantification of this assay is the melatonin concentration in saliva that can be measured with an inter-assay coefficient of variation (CV) of less than 30 %. The LoQ was determined from 7 different samples from 1.3 – 47.3 pg/mL each sample measured 17 times in duplicate in independent runs.

**Dilution Linearity/Parallelism: 92.2 %.** Three saliva samples with high amount of melatonin were sequentially diluted with incubation buffer and assayed according to the assay procedure. The results are presented in table 6.

Due to the complex matrix of saliva samples dilution with incubation buffer higher than 1:8 will cause a decreased linearity. Therefore sample dilution with incubation buffer higher than 1:4 is not recommended.

**Spiking Recovery: 97.9%.** Two saliva samples from the same donor, one collected during daytime and one during night time were titrated against each other and assayed according the assay procedure twice, independently. The results are presented in table 7.

Due to the complex and individual nature of the saliva matrix direct spiking of saliva with melatonin can lead to decreased recovery rates.

**Specificity:** The 50% binding (cross-reactivity) of the melatonin antiserum with different compounds were tested in the Direct Saliva Melatonin Radioimmunoassay (RK-DSM) from BÜHLMANN and are presented in table 8.

---

## METHOD COMPARISON

The comparison was done with 78 saliva samples from 10 different donors collected at different daytimes. The samples were analyzed using the presented EK-DSM assay as well as the Direct Saliva Melatonin Radioimmunoassay (RK-DSM) from BÜHLMANN. The subsequent linear regression analysis resulted in a correlation factor of  $R^2 = 0.84$ , an intercept of 0.77 pg/mL and a slope of 1.21. The correlation is presented in figure 3.

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
Figure 2

### Intra-Assay Precision

Sample	Mean [pg/mL]	SD [pg/mL]	CV [%]
S01	1.7	0.19	11.2
S04	5.2	1.20	22.9
S03	13.7	1.48	10.8
S05	15.9	0.84	5.3
Mean			12.6

Table 4

### Inter-Assay Precision

Sample	Mean [pg/mL]	SD [pg/mL]	CV [%]
Ctrl low	2.6	0.59	23.8
S11	2.4	0.41	17.2
S12	4.6	1.32	28.8
Ctrl high	5.2	1.20	23.2
S13	13.7	3.05	22.3
Mean			22.9

Table 5

### Dilution Linearity

Sample	Dilution Factor	Observed [pg/mL]	Expected [pg/mL]	Recovery O/E [%]
S06	1:1	14.9	--	--
	1:2	7.8	7.5	104.7
	1:4	3.0	3.7	80.5
	[1:8]	[1.1]	[1.9]	[59.1]
S07	1:1	22.7	--	--
	1:2	11.6	11.4	102.2
	1:4	5.0	5.7	88.1
	[1:8]	[1.8]	[2.8]	[63.4]
S08	1:1	20.4	--	--
	1:2	13.1	13.6	96.3
	1:4	6.4	6.8	94.1
	[1:8]	[2.7]	[3.4]	[79.4]
Mean				92.2

Table 6

**APPENDIX I**

**TABLES AND FIGURES / TABELLEN UND ABBILDUNGEN / TABLES ET FIGURES / TABELLE E FIGURE / TABLAS Y FIGURAS**

**Spiking Recovery**

Sample	Titration Ratio S5/S8	Expected [pg/mL]	Observed [pg/mL]	Recovery O/E [%]
S5/S8	5/0	1.2	1.2	--
	4/1	4.3	4.4	102.8
	3/2	7.4	5.3	72.0
	2/3	10.4	9.8	93.9
	1/4	13.5	13.9	102.8
	0/5	16.6	16.6	--
S5/S8	5/0	0.9	0.9	--
	4/1	3.7	3.6	96.8
	3/2	6.5	6.4	97.9
	2/3	9.4	10.5	112.2
	1/4	12.2	12.8	105.1
	0/5	15.0	15.0	--
Mean				97.9

Table 7

**Specificity**

Compound	Crossreactivity [%]
melatonin	100
serotonin	< 0.001
6-sulfatoxymelatonin	< 0.001
N-acetylserotonin	0.045
5-hydroxy-indole acetic acid	< 0.001
5-methoxytryptamine	0.007
5-methoxytryptophane	< 0.001
2-methyl-5-hydroxytryptamine	< 0.001
5-methoxyorsoralen	< 0.001
5-methoxytryptophol	0.002
chloramelatonin	1.3
caffeine	< 0.001
caffeine acid	<0.001
soluble coffee	<0.001
soluble coffee decaffeinated	<0.001

Table 8

**Correlation EK-DSM/ RK-DSM**


Figure 3

---

## APPENDIX II

### REFERENCES / LITERATURREFERENZEN / RÉFÉRENCES / RIFERIMENTI / REFERENCIAS

1. Koorengevel, KM, *et al.* *A forced desynchrony study of circadian pacemaker characteristics in seasonal affective disorder.* *J Biol Rhythms* **17**, 463-75 (2002).
2. Wirz-Justice, A, *et al.* *No evidence for a phase delay in human circadian rhythms after a single morning melatonin administration.* *J Pineal Res* **32**, 1-5 (2002).
3. Graw, P, *et al.* *Early morning melatonin administration impairs psychomotor vigilance.* *Behav Brain Res* **121**, 167-72 (2001).
4. Danilenko, KV, *et al.* *Is sleep per se a Zeitgeber in humans.* *J Biol Rhythms* **18**, 170-8 (2003).
5. Vaughan G M: *New sensitive serum melatonin radioimmunoassay employing the Kennaway G280 antibody: Syrian hamster morning adrenergic response.* *J Pineal Res* **15**, 88-103 (1993).
6. Voultsios, A, *et al.* *Salivary melatonin as a circadian phase marker: validation and comparison to plasma melatonin.* *J Biol Rhythms* **12**, 457-466 (1997).

**Direct Saliva Melatonin  
Sample Pretreatment (Saliva& Controls)**

Clean polypropylene tube


200 µl Saliva Sample or Control  
25 µl Pretreatment Solution

*Vortex, 5 sec*

*Incubate, 10 min, 18-28°C*

25 µl Neutralization Solution


*Vortex, 5 sec*

*Centrifuge at 10'000 rpm for 5 min*

Proceed to ELISA procedure

**Direct Saliva Melatonin  
ELISA Procedure**

Precoated Microtiter Plate

*Wash 2x with ≥300µL wash buffer*

100 µl Calibrators, Pretreated Controls  
or Samples

*Incubate 16-20 hours at 2-8°C*

add 50 µl Melatonin-Biotin-Conjugate

*1 minute on a plate orbital shaker*

*3 hours at 2-8°C*

*Wash 4x with ≥300µL wash buffer*


add 100 µl Enzyme Label

*60 minutes at 18-28°C  
on a plate orbital shaker*

*Wash 4x*

add 100 µl TMB Substrate

*Incubate 30 minutes at 18-28°C  
on a plate orbital shaker*

add 100 µl Stop Solution

➔ Read absorbance at 450 nm (within 30 minutes)

---

**APPENDIX IV**

**NOTES / NOTIZEN / NOTES / NOTE / NOTAS**

---

**APPENDIX IV**

**NOTES / NOTIZEN / NOTES / NOTE / NOTAS**

APPENDIX V

SYMBOLS / SYMBOLE / SYMBOLES / SIMBOLI / SIMBOLOS

Symbol	Explanation
	Use By Verwendbar bis Utiliser jusqu'au Utilizzare entro Fecha de caducidad
<b>REF</b>	Order Code Bestellnummer Code Codice Código
<b>LOT</b>	Batch Code Lotbezeichnung Code du lot Codice del lotto Codigo de lote
<b>IVD</b>	<i>In Vitro</i> Diagnostic Medical Device <i>In Vitro</i> Diagnostikum Dispositif médical de diagnostic <i>in vitro</i> Dispositivo medico-diagnostico <i>in vitro</i> Producto sanitario para diagnóstico <i>in vitro</i>
	Contains sufficient for <n> tests Ausreichend für „n“ Ansätze Contenu suffisant pour „n“ tests Contenuto sufficiente per „n“ saggi Contenido suficiente para <n> ensayos
	Consult Instructions for Use- Gebrauchsanweisung beachten Consulter le mode d'emploi Consultare le istruzioni per l'uso Consulte las instrucciones de uso
	Temperature Limitation Zulässiger Temperaturbereich Limites de température Limiti di temperatura Limite de temperatura
<b>REAG PRE</b>	Pretreatment Reagent Vorbehandlungs Reagenz réactif de prétraitement reagente di pretrattamento reactivo del tratamiento previo
<b>SOLN NEUT</b>	Neutralizing Solution Neutralisierungs-Lösung Solution neutralisante Soluzione neutralizzante Solución neutraliza
<b>MP</b>	Microtiter plate Mikrotiterplatte Microplaque Micropiastra Placa de microtitulación

Symbol	Explanation
<b>BUF WASH 10X</b>	Wash Buffer Concentrate (10x) Wasch-Puffer Konzentrat (10x) Concentré de tampon de lavage (10x) Tamponi di lavaggio concentrato (10x) Tampón de lavado concentrado (x10)
<b>REAG BLANK</b>	Blanking Reagent Nullwert-Reagenz Réactif blanc Reagente bianco Reactivo blanco
<b>BUF INC</b>	Incubation Buffer Inkubations-Puffer Tampon d'incubation Tamponi di incubazione Tampón de incubación
<b>CAL A - CAL E</b>	Calibrator A - E Kalibrator A - E Calibrateur A - E Calibratore A - E Calibrador A - E
<b>CONTROL L</b>	Control Low Kontrolle tief Contrôle bas Controllo basso Control bajo
<b>CONTROL H</b>	Control High Kontrolle hoch Contrôle élevé Controllo alto Control alto
<b>BC</b>	Biotin Conjugate Biotin-Konjugat Conjugué Biotine Coniugato biotinilato Conjugado de Biotina
<b>EL</b>	Enzyme Label Enzym-Marker Marqueur enzymatique Marcato enzimatico Marcador enzimático
<b>SUBS TMB</b>	TMB Substrate TMB-Substrat Substrat TMB Substrato di TMB Substrato de TMB
<b>SOLN STOP</b>	Stop Solution Stopp-Lösung Solution stop Soluzione stoppante Solución de interrupción