Quantum Blue[®] fCAL extended

Quantitative Lateral Flow Assay

LF-CALE25 25 tests

Revision date: 2018-02-15

ENGLISH

INTENDED USE

The BÜHLMANN Quantum Blue® fCAL extended is an *in vitro* diagnostic test for the quantitative determination of calprotectin in human stool specimens intended as an aid in the assessment of intestinal mucosal inflammation. The assay results can be used as an aid to diagnosis in distinguishing organic, inflammatory disease of the gastrointestinal tract (inflammatory bowel disease, IBD, e.g. Crohn's disease or ulcerative colitis, UC) from functional disease (irritable bowel syndrome, IBS) (ref. 1-7), in patients with chronic abdominal pain, above the age of four (ref. 8-9), and as an aid to IBD disease monitoring (ref. 10-20).

For professional use only*.

*Canada, Taiwan: For laboratory use only.

PRINCIPLE OF THE ASSAY

The test is designed for the selective measurement of calprotectin antigen by sandwich immunoassay. A monoclonal capture antibody (mAb) being highly specific for calprotectin is coated onto the test membrane. A second monoclonal detection antibody conjugated to gold colloids is deposited onto the conjugate release pad and released into the reaction system after addition of the extracted and diluted stool sample. The calprotectin / anti-calprotectin gold conjugate binds to the anti-calprotectin antibody coated on the test membrane (Test Line) and the remaining free anti-calprotectin gold conjugate binds to the goat anti-mouse antibody coated on the test membrane (Control Line). The signal intensities of the Test Line and the Control Line are measured quantitatively by the BÜHLMANN Quantum Blue® Reader.

REAGENTS SUPPLIED AND PREPARATION

Reagents	Quantity	Code	Comments
Test Cassette	25 pieces	B-LFCALUS-TC	vacuum-sealed in a foil pouch
Extraction Buffer	1 bottle 125 mL	B-CAL-EX	Ready to use
Controls Low* / High*	2 vials 0.5 mL	B-CALE-CONSET	Ready to use
RFID Chip Card	1 piece	B-CALE-RCC	White plastic card
RFID Chip Card	1 piece	B-CALE-RCC720	Green plastic card

Table 1

* The controls contain lot specific amounts of native human calprotectin. Refer to the additional QC data sheet for actual concentrations.

STORAGE AND SHELF LIFE OF REAGENTS

All kit components are stable at 2-8 °C until the expiration date printed on the labels.

REAGENTS & MATERIAL SUPPLIED SUPPLEMENTARY

Fecal extraction devices

Fecal extraction devices described below are not delivered with the kit and either of them has to be ordered with the kit.

Extraction devices	Quantity	Code
CALEX® Cap device	50, 200 or 500 tubes filled with 5 mL extraction buffer, ready to use	B-CALEX-C50 B-CALEX-C200 B-CALEX-C500
Smart-Prep	50 tubes consisting of spatulas and base caps	B-CAL-RD
ScheBo® Quick-Prep™	50 tubes consisting of tube, cone & dosing tip. Filled with 1.3 mL extraction buffer, ready to use	B-CAL-SO50

Table 2

MATERIALS REQUIRED BUT NOT PROVIDED

- Vortex mixer for stool extraction
- Precision pipettes with disposable tips: 10-100 µL, 100-1000 µL and 250-2500 µL
- Centrifuge
- 5 mL polypropylene or polystyrene tubes for dilution of the extracts
- Timer (optional)
- Quantum Blue® Reader available from BÜHLMANN (order code: BI-POCTR-ABS)
- Soft tissues or blotting paper

PRECAUTIONS

Safety precautions

- The controls of this test contain components of human origin. Although tested and found negative for HBV surface antigen, HCV and HIV1/2 antibodies, the reagents should be handled as if capable of transmitting infections and should be handled in accordance with Good Laboratory Practices (GLP) using appropriate precautions.
- Patient specimens should be handled as if capable of transmitting infections and should be handled in accordance with Good Laboratory Practice (GLP) using appropriate precautions.
- Reagents: Avoid contact of reagents with the skin, eyes, or mucous membranes. If contact does occur, immediately wash with generous amounts of water; otherwise, irritation can occur.
- Unused solution should be disposed of according to local state and federal regulations.

Technical precautions

Kit components

- All reagents and test samples must be equilibrated to room temperature (18-28 °C) before starting the assay.
- Mix well (vortex) the reagents before use.
- Components must not be used after the expiry date printed on the labels.
- Do not mix different lots of reagents.
- Test cassettes cannot be re-used.

Test procedure

- Read carefully the instructions prior to carrying out the assay. Assay performance will be adversely affected, if reagents are incorrectly diluted, handled or stored under conditions other than those as detailed in this instruction for use.
- The Quantum Blue® Reader must be switched on and programmed for the Quantum Blue® fCAL extended. B-CALE-RCC or B-CALE-RCC720 is to use before starting the assay (see Quantum Blue® Reader Manual).
- Use the RFID chip card in order to change lot-specific test parameters.
- Patient samples that are not properly handled may cause inaccurate results.
- In order to receive reliable and quantitative results it is important to homogenize the stool sample entirely in the extraction device.
- With BÜHLMANN Smart Prep and ScheBo® Quick-Prep™, it is important to centrifuge the extracts before storage. Centrifuge the tubes for 5 minutes at 3000 x g. After centrifugation the supernatant must be transferred into a fresh storage tube. With CALEX® Cap device you can store the extracts directly with the device. For re-use / re-measurement of the extracts see step 2 under the chapter assay procedure.

SPECIMEN COLLECTION AND STORAGE

If the extraction devices are used, less than 1 g of native stool specimen is needed for the extraction procedure.

Stool specimens should be collected into plain tubes.

Important: The specimen must be collected without any chemical or biological additions in the collection device.

Specimen transport

Stool specimens should be received by the laboratory within three days of collection. The specimens may be transported at room temperature (23 °C).

Specimen storage

Received stool specimens should be stored at 2-8 °C and extracted within three days.

Extract storage

Calprotectin in extracts obtained with the BÜHLMANN CALEX® Cap is stable at room temperature for three days, at 2-8 °C for six days and at -20 °C for 18 months.

Calprotectin in extracts obtained by manual weighing method, by BÜHLMANN Smart-Prep or by ScheBo® Quick-Prep™ is stable at 2-8 °C for six days or at -20 °C for 18 months.

ASSAY PROCEDURE

The assay procedure consists of three steps:

1. Extraction of stool samples


The extraction is described in the instruction for use delivered with the respective extraction devices.

2. Sample processing

- Smart Prep or ScheBo® Quick Prep™: Let the stool extract settle for 10 minutes after extraction. Dilute the supernatant 1:10 with extraction buffer (e.g. 50 µL stool extract and 450 µL extraction buffer) and mix well. Let the samples equilibrate for at least 5 minutes at 18-28 °C prior to proceeding to the next step (step no. 3).
- CALEX® Cap device: After extraction, let the stool extract settle for 10 minutes with the white head of the device down. Unscrew the blue cap. The supernatant can be used without further dilution in the lateral flow assay.


3. Lateral flow assay procedure and readout

Two alternative methods can be loaded from the respective RFID chip card: B-CALE-RCC or B-CALE-RCC720. Select one of the RFID chip cards before starting the experiments. Load the test method from the RFID chip card.


3.1. Method <B-CALE-RCC720> with internal timer

- Use the green plastic card.
- Load the test cassette onto the test cassette holder of the Reader.
- Add 60 µL of diluted stool extract onto the sample loading port of the test cassette.
- Close the test cassette holder and start the measurement by pressing the start button.
- The scan starts automatically after 12 minutes (720 seconds).
- For low / high controls: Repeat step 3.1 with 60 µL controls instead of diluted stool extract.

3.2. Method <B-CALE-RCC> without internal timer

- Use the white plastic card.
- Add 60 µL of diluted stool extract onto the sample loading port of the test cassette.
- Incubate for 12 minutes +/- 1 minute (set a timer manually).
- Load the test cassette onto the test cassette holder of the Reader.
- Scan the test cassette with the Quantum Blue® Reader by pressing the start button immediately.
- For low / high controls: Repeat step 3.2 with 60 µL controls instead of diluted stool extract.

Remark: Please refer to the Quantum Blue® Reader Manual to learn more about the basic functions and how to initialize and operate the Reader, especially how to select test methods, and how to load lot specific parameters from the RFID chip card in order to get the samples measured.

QUALITY CONTROL

- If the performance of the assay does not correlate with the established limits and repetition excludes errors in technique, check the following issues: i) pipetting, temperature controlling and timing ii) expiration dates of reagents and iii) storage and incubation conditions.
- Result of the self-test of the Quantum Blue® Reader performed at startup of the instrument has to be valid.

STANDARDIZATION

- The Quantum Blue® fCAL extended is standardized with the BÜHLMANN fCAL® ELISA (order code: EK-CAL).
- The BÜHLMANN Quantum Blue® Reader uses a lot-specific standard curve to calculate the calprotectin concentration. The assay range is between 30 and 1000 µg/g.
- To receive quantitative results for calprotectin concentration between 850-1800 µg/g, high samples reading above 850 µg/g can be re-tested with the BÜHLMANN Quantum Blue® high range test (order code: LF-CHR25).

VALIDATION OF RESULTS

- For a valid test result, the Control Line (C) must be visible in any case (see figures 1A and 1B). It is used as functional test control only and cannot be used for the interpretation of the Test Line (T). If the Test Line (T) is not detectable after 12 minutes of incubation time (figure 1A), the concentration of calprotectin present in the stool sample is below the detection limit. If a Test Line (T) is detectable after 12 minutes of incubation time (figure 1B), the calprotectin concentration present in the stool sample is calculated by the Quantum Blue® Reader.
- If only the Test Line (T) is detectable after 12 minutes of incubation time (figure 1C), the test result is invalid and the assay has to be repeated using another test cassette.
- If neither the Control Line (C) nor the Test Line (T) are detectable after 12 minutes of incubation time (figure 1D), the test result is invalid and the assay has to be repeated using another test cassette.
- As the Quantum Blue® Reader allows a quantitative evaluation of the Test (T) and Control (C) Lines, an additional validity check of the Control Line (C) is undertaken. If the signal intensity of the Control Line (C) is below a threshold after 12 minutes of incubation time, the test result is also invalid and the assay has to be repeated using another test cassette.

LIMITATIONS

- Reagents delivered with the BÜHLMANN Quantum Blue® fCAL extended kit are intended for the determination of calprotectin levels in human stool samples only.
- Fecal calprotectin values are intended as an aid to diagnosis in distinguishing organic disease from functional disease and as an aid to IBD monitoring. Results should always be interpreted in combination with other clinical and laboratory findings.
- For IBD disease monitoring, multiple fecal calprotectin measurements performed at up to 4 weeks intervals have been suggested to have best diagnostic accuracy in predicting clinical relapse in patients (ref. 21-22).
- In rare cases, when calprotectin levels are extremely high (above 4000 µg/g, ie. in acute UC), the test system may be prone to a high dose hook effect, that can result in values below the expected 850 µg/g upper limit of the linear range. It is advised to give particular attention to results above 300 µg/g when accompanied by strong symptoms, which could indicate acute inflammation. In this case retesting of a patient's stool sample by a laboratory in timely manner is recommended for confirmation.
- Patients who are taking NSAIDs regularly may have elevated fecal calprotectin levels.

INTERPRETATION OF RESULTS

DISTINGUISHING ORGANIC DISEASE FROM FUNCTIONAL GASTROINTESTINAL DISEASE

Determination of fecal calprotectin levels can be used as a reliable and simple aid in distinguishing organic from functional gastrointestinal diseases (ref. 1-7).

The following data were established with the BÜHLMANN fCAL® ELISA (order code: EK-CAL). In a clinical study, fecal calprotectin values of 344 symptomatic patients were compared with endoscopic findings. Endoscopy examination showed 264 patients with functional diseases, whereas 80 patients had various organic diseases (colitis, Crohn's, ulceration, diverticulitis, polyps, adenomas, cancer, or infectious diseases) (ref. 5).

ROC curve analysis (AUC: 0.962) resulted in an optimal clinical cut-off at 50 µg/g for the differentiation between organic and functional diseases, with a clinical sensitivity and specificity of 88.8 % and 93.6 %, respectively (refer to table 3). Fecal calprotectin levels from adults and children are comparable, whereas levels of newborns can be significantly increased (ref. 8).

These data support the following recommendation for interpretation of results:

Values below 50 µg/g:

Calprotectin values below 50 µg/g are not indicative of inflammation in the gastrointestinal tract. Patients with low calprotectin levels are not likely to be in need of invasive procedures to determine the inflammation cause (ref. 4).

Elevated values between 50 and 200 µg/g:

Calprotectin values between 50 and 200 µg/g can represent mild organic disease such as inflammation caused by NSAIDs, mild diverticulitis and IBD in remission phase. The low inflammatory response shown within this range may suggest repeating the measurement and performing further investigations.

Elevated values above 200 µg/g:

Calprotectin values above 200 µg/g are indicative of active organic disease with inflammation in the gastrointestinal tract. Appropriate further investigative procedures by specialists are suggested.

The cut-off level suggested for adults (50 µg/g) can be used for children aged from 4 to 17 years regardless of sex (ref. 8-9).

INTERPRETATION OF RESULTS

IBD MONITORING

The determination of fecal calprotectin is also a reliable and simple way to assist the monitoring of IBD patients (ref. 10-20).

Correlation of calprotectin levels and the inflammatory status of patient's intestinal mucosa, according to endoscopic evaluations, were determined in three independent studies using BÜHLMANN calprotectin tests (table 4). The diagnostic value of calprotectin in predicting clinical remission and relapse, according to patient's symptoms, clinical activity indices, unplanned need for therapy escalation, hospitalization or emergency was determined in three studies using BÜHLMANN calprotectin tests (table 5).

Condensed knowledge of published cut-offs and the above clinical performance studies support the following result categories:

Values below 100 µg/g:

Fecal calprotectin levels below 100 µg/g can reliably indicate patients, with low risk of clinical relapse, in endoscopic remission for whom invasive endoscopic procedures can be avoided (ref. 10-20).

Values between 100-300 µg/g:

Fecal calprotectin levels between 100-300 µg/g may indicate the necessity of tighter control in the following period to assess disease development tendencies.

Values above 300 µg/g:

Fecal calprotectin levels above 300 µg/g should be repeated and, if raised levels are confirmed, prompt further investigative procedures (ref. 10-20).

The above result categories are recommendations. It is advised that healthcare practitioners establish individual patient thresholds by determining the patient's baseline calprotectin level during disease remission.

A false negative result that is a calprotectin result below 100 µg/g category that should show a value above 300 µg/g, for a patient with endoscopic inflammation, although unlikely, may delay appropriate clinical decisions and patient treatment. Therefore, it is important that the patient remains under the care of a healthcare practitioner and reports any clinical symptoms.

Studies have shown that high calprotectin levels, above 300 µg/g, will not always indicate development of a clinical relapse (ref. 10-20). High calprotectin levels should be treated as a red flag signal and repeated. Confirmation of raised levels should prompt further investigative procedures.

PERFORMANCE CHARACTERISTICS

Method comparison

Bias at clinical decision points:

50 µg/g: 10.6 % (95 % CI: -2.1-36 %)

200 µg/g: -2.2 % (95 % CI: -4.9-4.7 %)

100 µg/g: 2.1 % (95 % CI: -2.7-14 %)

300 µg/g: -3.6 % (95 % CI: -7.0-2.3 %)

The method comparison study was performed according to the CLSI guideline EP09-A3. 186 extracted clinical stool samples were measured according to the instruction for use with the Quantum Blue® fCAL extended test and with the BÜHLMANN fCAL® ELISA. Measurements were performed over three days using three Quantum Blue® fCAL extended test cassette lots (figure 2).

Recovery: 102 %-121 %

Six stool specimen extracts were spiked with 150 µg/g calprotectin in calibrator material of human serum origin. The baseline extract was spiked with the corresponding amount of extraction buffer. Baseline and spiked samples were measured in eight replicates. One test cassette lot was used. The results are summarized in table 6.

Repeatability: 15.3-19.1 % CV

Within-laboratory precision: 18.0-23.0 % CV

Repeatability and within-laboratory precision were determined according to CLSI guideline EP05-A2. Four extracted stool samples, including extracts with calprotectin levels close to clinical decision points, were measured over ten days, in two independent runs each day, with two replicates per run. A single reagent lot was used (table 7).

Inter-lot precision: 16.5-20.6 % CV

Inter-lot precision was determined according to CLSI guideline EP05-A2. Four extracted stool samples, including extracts with calprotectin levels close to clinical decision points, were measured using three different reagent lots. The measurements were performed over five days, in a single run each day, with two replicates per run (table 8).

Limit of Blank (LoB): 6.7 µg/g calprotectin.

The LoB was established according to CLSI guideline EP17-A using extraction buffer to obtain 60 blank values. The study was performed with two different test cassette lots.

Limit of Detection (LoD): 18 µg/g calprotectin.

The LoD was established according to CLSI guideline EP17-A. Two extracted stool samples with calprotectin concentrations of 58 and 62 µg/g were diluted in extraction buffer to obtain a total of six samples with concentrations ranging from 1 x LoB (6.3 µg/g) to 4 x LoB (25.1 µg/g). The samples were measured in ten replicates to obtain 60 values. The study was performed with two different test cassette lots.

Limit of Quantification (LoQ):

Lower LoQ: ≤30 (28.2) µg/g calprotectin.

Upper LoQ: ≥1000 (1002) µg/g calprotectin.

The LoQ was established according to CLSI guideline EP17-A. To determine the lower LoQ, four extracted stool samples with calprotectin concentrations ranging from 19.1 to 37.3 µg/g were measured in ten replicates to produce 40 values (table 9). To determine the upper LoQ, four extracted stool samples with calprotectin concentrations ranging from 628 to 1001.7 µg/g were measured in ten replicates to produce 40 values (table 10). The study was performed with two different test cassette lots. Reference calprotectin values of extracted stool samples were determined with the BÜHLMANN fCAL® ELISA. The relative total error was calculated using the RMS model from estimates of precision and bias to the reference values, for each sample. The LoQ was defined as the lowest and highest sample concentration, for the lower LoQ and upper LoQ, respectively, which met the acceptance criterion of 30 % relative total error.

Linearity: 30-850 µg/g

The linear range of the Quantum Blue® fCAL extended assay was determined according to CLSI guideline EP06-A. Two extracted stool samples with low and high calprotectin concentrations were blended to obtain a total of 14 concentration levels covering and exceeding the expected measuring range. The blends were assayed in ten replicates on two test cassette lots. Mean calprotectin concentration values of each blend were plotted against the dilution factor used to obtain the blend. Linear as well as second and third order polynomial fitting was applied. Where coefficients of the second and third order polynomial fits were determined to be significant, the linear range was defined as the interval of calprotectin concentration in which deviation from the linear fit did not exceed 20 % relative concentration value or 20 µg/g (figure 3).

High dose hook effect

No high dose hook effect was observed for calprotectin concentrations up to 1500 µg/g. Decrease in mean signal below the 850 µg/g upper linear range limit was estimated for calprotectin concentrations above 4000 µg/g. No value below the highest clinical decision point of 300 µg/g was observed for any of the single replicate results for all high samples tested. Seven to eight extracted stool samples with calprotectin concentrations ranging from 1361 µg/g to 13'817 µg/g were measured in five replicates on three different test cassette lots.

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
Figure 1A Figure 1B Figure 1C Figure 1D

Figure 1

Clinical study – distinguishing organic disease from functional gastrointestinal disease

Calprotectin ¹ vs. clinical diagnosis	
Patient number and demographics	344 Age Median (IQR): 63 (53-72), 41.6 % male
Cut-off	50 µg/g
AUC	0.962
Sensitivity	88.8 %
Specificity	94.3 %
NPV	96.5 %
PPV	82.6 %
LR+	15.6
LR-	0.12

Table 3

¹ Results were obtained with the BÜHLMANN fCAL® ELISA.

Clinical studies – IBD monitoring

Calprotectin ¹ vs IBD activity determined by endoscopic findings	Study 1 Spain (ref. 10)	Study 2 Spain (ref. 11)	Study 3 Australia, New Zealand (ref.12)
Patient number and demographics	89 (CD ²) Ages: 32-58 44 % male	123 (UC ³) Ages: 18-85 66.4 % male	99 (CD ² after resection) Ages: 29-47 ³ 46.5 % male
Cut-off	272 µg/g	280 µg/g	100 µg/g
NPV	98 %	86 %	91 %
PPV	79 %	80.3 %	53 %

Table 4

¹ Results for studies 1 and 2 were obtained with the BÜHLMANN lateral flow assays (Quantum Blue® fCAL and Quantum Blue® fCAL high range). Results in study 3 were obtained with the BÜHLMANN fCAL® ELISA.

² CD = Crohn's disease patients

³ UC = Ulcerative Colitis patients

Clinical studies – IBD monitoring

Calprotectin ¹ vs future clinical remission or relapse	Study 4 UK (ref. 13)	Study 5 Spain (ref. 14)	Study 6 Spain (ref. 15)
Patient number and demographics	92 (CD ²) 38 % male	30 (CD ³) adalimumab therapy Ages: 24-64 43.3 % male	33 (CD ²) 20 (UC ³) infliximab therapy Ages: 18-68 47.2 % male
Follow-up time after calprotectin measurement	12 months	4 months	12 months
Patients in clinical relapse after follow-up	11 %	30 %	23 %
Cut-off	240 µg/g	204 µg/g	160 µg/g
NPV	96.8 %	100 %	96.1 %
PPV	27.6 %	75 %	68.7 %

Table 5

¹ Results for study 4 were obtained with the BÜHLMANN fCAL® ELISA. Results for studies 5 & 6 were obtained with the BÜHLMANN lateral flow assays (Quantum Blue® fCAL and Quantum Blue® fCAL high range).

² CD = Crohn's disease patients

³ UC = Ulcerative Colitis patients

Method comparison


Figure 2

APPENDIX I

TABLES AND FIGURES / TABELLEN UND ABBILDUNGEN / TABLES ET FIGURES / TABELLE E FIGURE / TABLAS E FIGURAS / TABELAS E FIGURAS

Recovery

Sample	1	2	3	4	5	6
Base value [$\mu\text{g/g}$]	77	111	141	240	232	349
Spike value [$\mu\text{g/g}$]	150	150	150	150	150	150
Expected value (base + spike) [$\mu\text{g/g}$]	227	261	291	390	382	499
Observed value (base + spike) [$\mu\text{g/g}$]	260	267	345	409	441	605
% Recovery (observed/expected)	115	102	118	105	115	121

Table 6

Within-laboratory precision

Mean Conc. [$\mu\text{g/g}$]	Repeatability CV [%]	Between -run Precision CV [%]	Between -day Precision CV [%]	Within-lab Precision CV [%]
51.7	19.1	10.2	7.6	23.0
118.6	18.3	0.0	3.3	18.6
295.0	15.3	5.3	7.9	18.0
647.8	16.8	11.1	5.9	21.0

Table 7

Reproducibility

Mean Conc. [$\mu\text{g/g}$]	Inter-lot precision CV [%]
49.2	16.5
112.2	17.0
278.4	20.1
682.1	20.6

Table 8

Limit of Quantification

Lower LoQ	Reference value (R) [$\mu\text{g/g}$]	Mean value obtained (O) [$\mu\text{g/g}$]	Bias [$\mu\text{g/g}$ (R-O)]	Precision [% CV]	Total Error [%]
Lot M0527	37.3	29.2	-8.1	17.5	25.7
	28.2	21.3	-6.9	16.7	27.8
	23.6	17.6	-6.0	25.6	31.9
	19.1	13.6	-5.5	20.6	32.2
Lot M2128	37.3	37.4	0.2	24.1	24.2
	28.2	25.1	-3.1	15.8	17.8
	23.6	21.0	-2.6	29.1	28.2
	19.1	18.5	-0.6	20.7	20.3

Table 9

Upper LoQ	Reference value (R) [$\mu\text{g/g}$]	Mean value obtained (O) [$\mu\text{g/g}$]	Bias [$\mu\text{g/g}$ (R-O)]	Precision [% CV]	Total Error [%]
Lot M0527	1001.7	752.6	-249.1	18.4	28.4
	746.0	706.9	-39.1	16.2	16.2
	678.6	704.2	25.6	14.0	15.1
	628.0	668.4	40.4	21.3	23.5
Lot M2128	1001.7	783.0	-218.7	22.6	28.1
	746.0	670.4	-75.6	13.9	16.1
	678.6	629.9	-48.7	17.2	17.5
	628.0	636.8	8.8	20.9	21.1

Table 10

Linearity


Figure 3

APPENDIX II

REFERENCES / REFERENZEN / RÉFÉRENCES / RIFERIMENTI / REFERENCIAS / REFERÊNCIAS

1. Fagerhol MK: *Calprotectin, a faecal marker of organic gastrointestinal abnormality*. Lancet 356, 1783-4 (2000)
2. Tibble JA et al.: *A simple method for assessing intestinal inflammation in Crohn's disease*. Gut 47,506-513 (2000).
3. Tibble JA et al.: *Use of surrogate markers of inflammation and Rome criteria to distinguish organic from nonorganic intestinal disease*. Gastroenterol 123, 450-460 (2002)
4. Jahnsen J, Røseth AG, Aadland E. *Measurement of calprotectin in faeces*. Tidsskr Nor Legeforen 128, 743-5 (2008)
5. Manz M et al. *Value of fecal calprotectin in the evaluation of patients with abdominal discomfort: an observational study*. BMC Gastroenterology 12, 5 (2012)
6. Pavlidis P. et al. *Diagnostic accuracy and clinical application of faecal calprotectin in adult patients presenting with gastrointestinal symptoms in primary care*. Scand J Gastroenterol. 48, 1048-54 (2013)
7. Fagerberg UL et al.: *Colorectal inflammation is well predicted by fecal calprotectin in children with gastrointestinal symptoms*. J Pediatr Gastroenterol Nutr 40, 450-5 (2005)
8. Qingling Zhu et al. *Fecal Calprotectin in Healthy Children Aged 1-4 Years*. PLoS ONE 11 (3) (2016)
9. Konikoff MR and Denson LA: *Role of fecal calprotectin as a biomarker of intestinal inflammation in inflammatory bowel disease*. Inflamm Bowel Dis 12(6), 524-34 (2006)
10. Lin JF et al. *Meta-analysis: fecal calprotectin for assessment of inflammatory bowel disease activity*. Inflamm Bowel Dis. Aug;20(8), 1407-15 (2014)
11. Lobatón T et al. *A new rapid test for fecal calprotectin predicts endoscopic remission and postoperative recurrence in Crohn's disease*. J Crohns Coliti, 7(12), 641-51 (2013)
12. Lobatón T et al. *A new rapid quantitative test for fecal calprotectin predicts endoscopic activity in ulcerative colitis*. Inflamm Bowel Dis. 19(5), 1034-42 (2013)
13. Wright EK et al., *Measurement of fecal calprotectin improves monitoring and detection of recurrence of Crohn's disease after surgery*. Gastroenterology. 148(5), 938-947 (2015)
14. Naismith GD et al. *A prospective evaluation of the predictive value of faecal calprotectin in quiescent Crohn's disease*. J Crohns Colitis. 8, 1022-9 (2014)
15. Ferreiro-Iglesias R et al. *Usefulness of a rapid faecal calprotectin test to predict relapse in Crohn's disease patients on maintenance treatment with adalimumab*. Scand J Gastroenterol. 23, 1-6 (2015)
16. Ferreiro-Iglesias R1 et al. *Fecal calprotectin as Predictor of Relapse in Patients with Inflammatory Bowel Disease Under Maintenance Infliximab Therapy*. J Clin Gastroenterol. 50(2), 147-51 (2015)
17. Guardiola J. et al. *Fecal Level of calprotectin Identifies Histologic Inflammation in Patients with Ulcerative Colitis In Clinical And Endoscopic Remission*. Clinical Gastroenterology and Hepatology. 12(11), 1865-70 (2014)
18. Lasson A et al. *Pharmacological intervention based on fecal calprotectin levels in patients with ulcerative colitis at high risk of a relapse: A prospective, randomized, controlled study*. United European Gastroenterol J. 3(1), 72-9 (2015)
19. Bressler B et al. *Clinicians' guide to the use of fecal calprotectin to identify and monitor disease activity in inflammatory bowel disease*. Can J Gastroenterol Hepatol. 29(7), 369-72 (2015)
20. Peyrin-BL et al. *Selecting Therapeutic Targets in Inflammatory Bowel Disease (STRIDE): Determining Therapeutic Goals for Treat-to-Target*. Am J Gastroenterol. 110, 1324-38 (2015)
21. Molander P et al. *Does Fecal calprotectin Predict Short-Term Relapse After Stopping Tnfalpha-Blocking Agents In Inflammatory Bowel Disease Patients In Deep Remission?* Journal of Crohn's and Colitis, 33-40 (2015)
22. De Vos M et al. *Consecutive fecal calprotectin measurements to predict relapse in patients with ulcerative colitis receiving infliximab maintenance therapy*. Inflamm Bowel Dis. 9, 2111-2117 (2013)

APPENDIX IV

SYMBOLS / SYMBOLE / SYMBOLES / SIMBOLI / SIMBOLOS

Symbol	Explanation	Symbol	Explanation
	Use by Verwendbar bis Utiliser jusqu'au Utilizzare entro Fecha de caducidad Validade até		Temperature limitation Zulässiger Temperaturbereich Limites de température Limiti di temperatura Limite de temperatura Limite de temperatura
REF	Catalogue number Bestellnummer Référence du catalogue Numero di catalogo Número de catálogo Número de catálogo	TC	Test cassette Test Kassette Cassette de testi Card di rilevazione Cartucho de prueba Cassete de teste
LOT	Batch code Chargenbezeichnung No. du lot No. del lotto No. de lote No. do lote	BUF EX	Extraction buffer Extraktionspuffer Tampon d'extraction Tampone di estrazione Tampón de extracción Tampão de extração
IVD	<i>In Vitro</i> Diagnostic Medical Device <i>In Vitro</i> Diagnostikum Dispositif médical de diagnostic <i>in vitro</i> Dispositivo medico-diagnostico <i>in vitro</i> Producto sanitario para diagnóstico <i>in vitro</i> Produto para Diagnóstico <i>in vitro</i>	CONTROL L	Control low Kontrolle Tief Contrôle Bas Controllo Basso Control Bajo Controle Baixo
	Content sufficient for <n> tests Ausreichend für "n" Ansätze Contenu suffisant pour „n“ tests Contenuto sufficiente per „n“ saggi Contenido suficiente para <n> ensayos Conteúdo suficiente para <n> testes	CONTROL H	Control high Kontrolle Hoch Contrôle Élevé Controllo Alto Control Alto Controle Alto
	Consult instructions for use Gebrauchsanweisung beachten Consulter le mode d'emploi Consultare le istruzioni per l'uso Consulte las instrucciones de uso Consultar as intruções de uso	RCC	RFID chip card RFID Chipkarte RFID Carte à puce RFID Carta chip RFID Tarjeta chip RFID Cartão Chip