Quantum Blue[®] C-Reactive Protein (CRP)

**Quantitative
Lateral Flow Assay**

LF-CRP25

25 tests

Revision date: 2013-07-29

ENGLISH

INTENDED USE

The Quantum Blue[®] C-Reactive Protein (CRP) Lateral Flow Assay is a sandwich-type immunoassay designed for the quantitative determination of CRP in human serum in combination with the BÜHLMANN Quantum Blue[®] Reader.

PRINCIPLE OF THE ASSAY

The test is designed for the selective measurement of CRP antigen by sandwich immunoassay. A monoclonal capture antibody (mAb) which is highly specific for CRP is lined onto the test membrane. A second monoclonal detection antibody conjugated to gold colloids is deposited in the conjugate release pad and released into the reaction system after addition of the test sample. The CRP/anti-CRP gold conjugate binds to the anti-CRP antibody coated onto the test membrane (test line; test band) and the remaining free anti-CRP gold conjugate binds to the goat anti-mouse antibody coated onto the test membrane (control line; control band). The signal intensities of the test line and the control line, respectively, are measured quantitatively by the BÜHLMANN Quantum Blue[®] Reader.

REAGENTS SUPPLIED AND PREPARATION

Reagents	Quantity	Code	Comments
Test Cartridge vacuum-sealed in a foil bag	25 pieces	B-LFCRP-TC	
Chase Buffer	1 vial 10 ml	B-LFCRP-CB	Ready to use
RFID Chip Card red	1 piece	B-LFCRP-RCCR	To be used with Firmware 1.0.4.x
RFID Chip Card white	1 piece	B-LFCRP-RCC	To be used with Firmware 2.2.17

Table 1

STORAGE AND SHELF LIFE OF REAGENTS

All kit components are stable at 2-8°C until the expiration date printed on the labels.

EQUIPMENT

- Precision pipettes with disposable tips: 10-200 µl
- Tubes for dilution of samples
- Timer (optional)
- Quantum Blue[®] Reader available from BÜHLMANN (order code: BI-POCTR-ABS)
- Gloves and Lab coat

PRECAUTIONS

SAFETY PRECAUTIONS

- Patient samples should be handled as potentially bio-hazardous (eg Hepatitis B virus and HIV). Please take appropriate precautions.
- Ensure that cuts, abrasions and skin lesions are properly protected and take precautions to avoid self-inoculation or splashing of mucous membranes.
- Wear appropriate personal protective equipment to avoid contact with eyes and skin.
- Unused solution should be disposed of according to local State and Federal regulations.

TECHNICAL PRECAUTIONS

- Components must not be used after the expiry date printed on the labels. Do not mix different lots of reagents.
- Test performance will be adversely affected if reagents are incorrectly diluted, modified or stored under conditions other than those as detailed in: storage and shelf life of reagents.
- Avoid contamination of reagents.

- Let the reagents adjust to reach room temperature. Mix well (vortex) the reagents before use.
- Cartridges cannot be re-used.
- Test results should be interpreted in conjunction with information available from clinical assessment of the patient and other diagnostic procedures.

SPECIMEN COLLECTION AND STORAGE

Less than 50 µl of serum is required.

Collect blood into plain venipuncture tubes (for Serum), avoid hemolysis, mix by inverting sample tube several times and leave to clot for 45 minutes at room temperature (18-28°C) protected from light. Centrifuge at 1800 x g for 15 minutes at room temperature (18-28°C) and collect the serum. Do not heat-inactivate samples.

Serum samples can be stored refrigerated at 2-8°C for up to 7 days. For longer storage, keep serum samples at -20°C. The samples are stable for at least 4 months.

Note: The use of lipaemic, hemolytic and icteric samples has to be validated for this assay. Lipaemic samples can be avoided by asking patients to fast for at least 12 hours prior to the sample being taken.

ASSAY PROCEDURE

The Assay Procedure consists of two individual parts:

1. Dilution of serum samples with chase buffer
Prior measurement dilute the serum sample 1:10 with Chase Buffer.
For example: Mix 15 µl serum sample with 135 µl Chase Buffer (B-LFCRP-CB).
2. Lateral Flow Assay/Reading
Add 100 µl of diluted serum sample onto the sample loading port of the test cartridge.
3. Lateral Flow Assay Procedure and Readout

There are two alternative methods available on the Quantum Blue[®] Reader: <CRP_22> and <CRP_0>. Select one of these methods before starting the experiments.

Load the lot specific parameters from the RFID Chip Card.

1. Method <CRP_22> with internal timer
 - Load the test cartridge onto the test cartridge holder of the Reader
 - Add 100 µl of diluted serum sample onto the sample loading port of the test cartridge
 - Close the cartridge holder and start the measurement by pressing the start button
 - The scan starts automatically after 12 minutes (1320 seconds).
2. Method <CRP_0> without internal timer
 - Add 100 µl of diluted serum sample onto the sample loading port of the test cartridge
 - Incubate for 22 minutes +/- 1 minute (set a timer manually)
 - Load the test cartridge onto the test cartridge holder of the Reader
 - Scan the cartridge with the Quantum Blue[®] Reader by pressing the start (<ENTER>) button immediately.

Remark: Please refer to the Quantum Blue[®] Reader Manual to learn about the basic functions and how to initialize and operate the Reader, especially how to select test methods, and how to load lot specific parameters from the RFID Chip Card in order to get the samples measured.

QUALITY CONTROL

- A thorough understanding of this instruction for use is necessary for the successful use of the product. Reliable results will be obtained only by following this instruction for use accurately.
- Patient samples that are not properly handled may cause inaccurate results.
- If the precision of the assay does not correlate with the established limits and repetition excludes errors in technique, check the following issues: i) pipetting, temperature controlling and timing devices, ii) expiration dates of reagents and iii) storage and incubation conditions.

VALIDATION OF RESULTS

- For a valid test result, the Control Line (C) must be visible in any case (see Figure 1A and 1B). It is used as functional test control only and cannot be used for the interpretation of the Test Line (T). If the Test Line (T) is not detectable after 22 minutes of incubation time (Figure 1A), the concentration of CRP present in the blood sample is below the detection limit. If a Test Line (T) is detectable after 22 minutes of incubation time (Figure 1B), the CRP concentration present in the blood sample is calculated by the Quantum Blue[®] Reader.
- If only the Test Line (T) is detectable after 22 minutes of incubation time (Figure 1C), the test result is invalid and the CRP assay has to be repeated using another Test Cartridge.
- If neither the Control Line (C) nor the Test Line (T) is detectable after 22 minutes of incubation time (Figure 1D), the test result is invalid and the CRP assay has to be repeated using another Test Cartridge.
- As the Quantum Blue[®] Reader allows a quantitative evaluation of the Test (T) and Control (C) Lines, an additional validity check of the Control Line (C) is undertaken. If the signal intensity of the Control Line (C) is below a specific threshold after 22 minutes of incubation time, the test result is also invalid and the CRP assay has to be repeated using another Test Cartridge.

STANDARDIZATION AND INTERPRETATION OF RESULTS

- The Lateral Flow Assay is calibrated to the NIBSC 85/506 reference material.
- The BÜHLMANN Quantum Blue[®] Reader uses a lot-specific standard curve to calculate the CRP concentration. This lot-specific standard curve is generated with accuracy pool samples with assigned CRP concentrations.
- For quantitative measurements of samples reading above 20 mg/L, dilute the serum samples 1:250 with chase buffer and assay them again according to the procedure. The measured concentration must then be multiplied by the reciprocal volume factor x 25 to obtain the final result.
- CRP values should be used as supplementary data available to the physician in establishing a diagnosis.

PERFORMANCE CHARACTERISTICS

- **Limit of Blank (LoB): ≤0.34 mg/L CRP.** The LoB has been established in accordance with CLSI protocol EP17-A in three independent runs using three different lots of Test Cartridges with 60 blank values in total by using CRP free serum as a sample.
- **Limit of Detection (LoD): ≤0.61 mg/L CRP.** The LoD has been established with CRP low serum. The samples

were measured with three different lots of Test Cartridges in three independent runs of 20 replicates each. The averaged SD values were determined and the LoD has there from been calculated in accordance with CLSI protocol EP17-A.

- **Limit of Quantification, LoQ: Lower LoQ: ≤1.0 mg/LCRP; Upper LoQ: ≥25 mg/L CRP.** The LoQ has been established with the precision profile. The limit of quantification corresponds to the concentration of CRP with an imprecision below 25% CV allowing a quantitative measurement between 1.0 (lower LoQ) and 25 mg/L (upper LoQ). Refer to Table 6.
- **Dynamic range:** The assay range is between 1.0 and 25 mg/L.
- **Precision: Repeatability: ≤19.3 % CV. Total precision: ≤24.0 % CV.** The repeatability of the Quantum Blue[®] CRP assay was calculated from 6 samples supplemented with purified CRP covering the assay range according to CLSI EP5-A2. Each sample was tested according to the assay procedure in one run using three lots with 20 replicates each. The repeatability was below 20 % CV. Refer to Table 7.
- **Inter-lot precision: ≤24.3 % CV; interlab precision: ≤25.2 % CV.** The inter-lot and interlab precision of the Quantum Blue[®] CRP assay was calculated from 6 samples supplemented with purified CRP covering the assay range. Each sample was tested according to the assay procedure in one run using three lots with 20 replicates each. The mean values of three different lots of Test Cartridges are presented in Table 8.
- **Linearity:** A sample containing 65 mg/L CRP was produced using CRP free serum and CRP stock solution. The sample was diluted in 16 steps using CRP free serum and subsequently assayed on four test cartridges. The mean was calculated for each dilution according to CLSI EP6-A. The results showed linearity within 0.80 and 35 mg/L. Refer to Table 9.
- **Recovery: 80-120%.** A sample containing 65 mg/L CRP was produced using CRP free serum and CRP stock solution. The sample was diluted in 16 steps using CRP free serum as subsequently assayed on four test cartridges. The mean was calculated for each for each concentration. The recovery was acceptable between 0.8 and 35 mg/L and with an additional dilution of 1:25 up to 100 mg/L. The recovery varied from 80-120 %.
- **High Dose Hook effect:** No High dose hook effect was observed up to a concentration of 200 mg/L.
- **Method Comparison:** $R^2 = 0.96; y=1.063x$
22 samples were compared with KoneLab CRP High Sensitivity. The correlation data are illustrated in Figure 2.

INTERFERING SUBSTANCES

Interference substances were evaluated in accordance with CLSI protocol EP7-A2.

Strongly hemolytic, lipemic or icteric samples with CRP concentration of up to 2 mg/L might interfere with the test leading to up to 40 % falsely elevated CRP results: The following substances were tested triglycerides (Intralipid[®] 13.2 mg/ml), conjugated bilirubin (0.2 mg/ml), unconjugated bilirubin (0.29 mg/ml) or haemoglobin (2 mg/ml).

Other substances and/or factors have not been investigated in this study. Interferences cannot be excluded.

EXPECTED VALUES

The following reference interval for serum (97.5th percentile) has been determined from asymptomatic blood donors:

97.5 th percentile	4.76 mg/L
Median:	1.12 mg/L
n = 99	

PERFORMANCE LIMITATIONS

Temperature

The assay has to be performed at a temperature between 18-22 °C.

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Figure 1


Figure 2:

Method comparison LF-CRP with Konelab „CRP High Sensitivity


Table 6:

LOB, LOD and LOQ

Lot	LOB (mg/L)	LOD (mg/L)	LOQ (mg/L)
Lot 1	0.21	0.61	1.0 to 30
Lot 2	0.20	0.52	0.6 to 25
Lot 3	0.34	0.58	0.6 to 25

Table 7

Summary of within-lab precision

Sample No.	Concentration [mean mg/L]	Repeatability [CV]	Total Precision [CV]
S1	1.16	17.6%	19.7%
S2	2.30	18.9%	21.8%
S3	3.25	19.3%	22.1%
S4	10.0	14.9%	19.4%
S5	16.5	19.1%	22.7%
S6	23.3	18.5%	24.0%

Table 8:

Inter-lot and inter-lab precision

Sample No.	Concentration [mean mg/L]	Inter-Lot Precision [CV]	Inter-Lab Precision [CV]
S1	1.04	21.3%	25.2%
S2	1.9	21.8%	17.9%
S3	2.9	15.4%	15.3%
S4	9.5	13.1%	18.1%
S5	15.9	24.3%	25.1%
S6	20.9	19.0%	23.3%

Table 9:

Linearity

Lot	Repeatability	Linear Range	Recovery (mean)
Lot 1	22.6%	0.78 – 43.3 mg/L	108.4%
Lot 2	23.5%	0.83 – 36.7 mg/L	102.0%
Lot 3	24.8%	0.81 – 36.1 mg/L	83.4%

REFERENCES/ LITERATURREFERENZEN/ REFERENCES/ RIFERIMENTI/ REFERENCIAS

- Powell LJ. C-reactive protein--a review. Am J Med Technol,45(2): 38-42 (1979)
- Gewurz H, Mold C, Siegel J, Fiedel B. C-reactive protein and the acute phase response. Adv Intern Med. 27:345-72 (1982)
- Patel VB, Robbins MA, Topol EJ. C-reactive protein: a 'golden marker' for inflammation and coronary artery disease. Cleve Clin J Med. 68(6):521-524, 527-34 (2001)
- Vermeire S, Van Assche G, Rutgeerts P. The role of C-reactive protein as an inflammatory marker in gastrointestinal diseases. Nat Clin Pract Gastroenterol Hepatol. 2(12):580-6 (2005)
- Vermeire S, Van Assche G, Rutgeerts P. Laboratory markers in IBD: useful, magic, or unnecessary toys?Gut.55(3):426-31 (2006)

SYMBOLS/ SYMBOLE/ SYMBOLES/ SIMBOLI/ SIMBOLOS

Symbol	Explanation	Symbol	Explanation
	Use By Verwendbar bis Utiliser jusqu'au Utilizzare entro Fecha de caducidad		Temperature limitation Zulässiger Temperaturbereich Limites de température Limiti di temperatura Limite de temperatura
REF	Catalogue number Bestellnummer Référence du catalogue Numero di catalogo Número de catálogo	TC	Test Cartridge Testkassette Cartouche de testi Card di rilevazione Cartucho de prueba
LOT	Batch code Chargenbezeichnung Code du lot Codice del lotto Codigo de lote	BUFCHASE	Chase Buffer Laufpuffer Tampon de chase Tampone di chase Tampón de chase
IVD	In Vitro Diagnostic Medical Device In Vitro Diagnostikum Dispositif médical de diagnostic in vitro Dispositivo medico-diagnostico in vitro Producto sanitario para diagnóstico in vitro	RCC	RFID Chip card RFID Chipkarte RFID Carte à puce RFID Carta chip RFID Tarjeta chip
	Contains sufficient for <n> tests Ausreichend für „n“ Ansätze Contenu suffisant pour „n“ tests Contenuto sufficiente per „n“ saggi Contenido suficiente para <n> ensayos		Consult Instructions for Use- Gebrauchsanweisung beachten Consulter le mode d'emploi Consultare le istruzioni per l'uso Consulte las instrucciones de uso