Direct Saliva MELATONIN

RIA

RK-DSM2

200 tests

Revision date: 2012-11-20

ENGLISH

INTENDED USE

The BÜHLMANN Direct Saliva Melatonin Radio Immunoassay (RIA) test kit is intended for the direct, quantitative determination of melatonin in human saliva (1-4).

PRINCIPLE OF THE ASSAY

The Direct Saliva Melatonin RIA kit measures melatonin by a double-antibody radio immunoassay based on the Kennaway G280 anti-melatonin antibody (5,6). Human saliva samples and reconstituted standards and controls are incubated with the anti-melatonin antibody and ¹²⁵I-melatonin. ¹²⁵I-melatonin competes with melatonin present in samples, standards and controls. After 20 hours of incubation, solid-phase second antibody is added to the mixture in order to precipitate the antibody-bound fraction. After aspiration of the unbound fraction, the antibody bound fraction of ¹²⁵I-melatonin is counted.

REAGENTS SUPPLIED AND PREPARATION

| Reagents | Quantity | Code | Reconstitution |
|--|-------------------|------------------|--|
| Incubation Buffer | 2 vials 100 ml | B-DSM-IB | Ready to use |
| Antiserum anti-melatonin antibody | 2 vials 11 ml | B-DSM-AS | Ready to use |
| Tracer ¹²⁵ I-melatonin | 2 vials 11 ml | B-MEL-TR | Ready to use |
| Calibrator A - E¹⁾ lyophilized melatonin standards | 5 vials | B-MEL-CASET | Reconstitute with 5 ml of Incubation Buffer |
| Controls low/high²⁾ lyophilized melatonin | 2 vials | B-DSM- CONSET | Reconstitute with 5 ml of Incubation Buffer |
| 2nd Antibody solid phase bound second antibody | 2 vials 11 ml | B-AB2 | Ready to use |

Table 1

¹⁾ After reconstitution the Calibrators A to E contain 0.5, 1.5, 5, 15 and 50 pg/ml of melatonin, respectively. Reconstitute each vial with 5 ml of Incubation Buffer. **Leave for at least 30 minutes at 2-8°C and vortex.**

²⁾ Lot specific amounts of melatonin. Reconstitute each vial with 5 ml of Incubation Buffer. **Leave for at least 30 minutes at 2-8°C and vortex.**

STORAGE AND SHELF LIFE OF REAGENTS

| Unopened reagents | |
|--|---|
| All unopened kit components are stable at 2-8°C until the expiration date. | |
| Opened / reconstituted reagents | |
| Incubation Buffer | Stable at 2-8°C until expiration date printed on the labels. |
| Antiserum | |
| Tracer | |
| Calibrators | Stable for at least 4 months after reconstitution at 2-8°C. |
| Controls | |
| 2 nd Antibody | Store refrigerated (Do not freeze!) Stable at 2-8°C until expiration date printed on the label. |

Table 2

PRECAUTIONS

SAFETY PRECAUTIONS

- Radioactive Material: This kit contains radioactive material which does not exceed 74 kBq of ¹²⁵Iodine.
- Receipt, acquisition, possession, use, and transfer are subject to the local regulations. Unused solutions and radioactive waste should be disposed of according to local State and Federal regulations.

TECHNICAL PRECAUTIONS

- Read carefully the instructions prior to carrying out the test. Test performance will be adversely affected if reagents are incorrectly diluted, modified or stored under

conditions other than those as detailed in this instruction for use.

- Incorrect results for standard curve, controls or samples may be obtained, if the 2nd antibody was not mixed sufficiently before pipetting. Do not freeze the 2nd antibody!
- Components must not be used after the expiry date printed on the labels.
- Do not mix different lots of reagents.
- Let the reagents adjust to reach room temperature. Reconstitute the lyophilized reagents as indicated. Mix well (vortex) the reagents before use.
- Counting time should be selected in order to keep statistical counting error small: e.g., at 2000 cpm the counting error is at 5%; at 10000 cpm it is only 1%.
- If the initial concentration of an unknown sample reads above the highest calibrator, the sample should be further diluted with incubation buffer and tested again according to the assay procedure.
- The procedure was tested and validated for human saliva samples. If other saliva specimen have to be used, it is recommended to either validate possible matrix effects with melatonin-free saliva specimen or extract the saliva samples by C18 reversed-phase column extraction (B-MEC).

MATERIALS REQUIRED BUT NOT PROVIDED

- BÜHLMANN Saliva Collection Devices, B-SLEEPCHECK16 or B-SVC (50).
- 100 µl, 400 µl, 1000 µl and 5000 µl precision pipettes (or preferably a 100-1000 µl adjustable multipipette) with disposable tips.
- Disposable polystyrene tubes for the RIA (preferably conical tubes; e.g. Sarstedt # 57.477).
- Deionized double distilled water
- Refrigerated centrifuge.
- Vortex mixer.
- Stir bar and magnetic stirrer.
- Aspiration device.
- Gamma-counter.

SALIVA COLLECTION AND STORAGE

Collection of Saliva

- Collect saliva using the BÜHLMANN Saliva Collection Device. One of these cotton swabs can absorb up to 3 ml of saliva. The procedure calls for 1 ml of saliva for duplicate determinations. **Do not use cotton swabs containing citric acid.**
- Subjects ideally refrain from eating during the sampling period. If this is not possible, subjects should be permitted to eat immediately after a collection only and rinse their mouths with water 15 minutes before the next collection.
- Avoid beverages containing artificial colorants as well as alcoholic beverages during the sampling period. Coffee is not permitted - it is counter indicated in sleep studies anyway. In view of some recent reports on melatonin in some foods, it is suggested that bananas should not be eaten during the sampling period.
- Subjects should avoid brushing their teeth, with or without toothpaste, during sampling periods. It is likely that subjects with gingivitis will contaminate the saliva with plasma or even whole blood leading to unknown consequences.
- Do not stimulate saliva flow by chewing gums or eating lemons.
- When collecting saliva at night, a dim flash light or a ≤100 lux yellow light should be used in order to avoid a possible light influence on the melatonin profile.

Storage:

Saliva samples absorbed in cotton swabs can be stored in the saliva collection device for up to 7 days at 2-8°C. If not assayed within one week after collection, samples should be frozen and can be stored for at least 6 months at ≤ -20°C. Repeated freeze-thaw cycles should be avoided. Thaw (if frozen) and centrifuge the saliva samples for 5 minutes at 1000 x g and 2-8°C prior to testing. Transfer the clear supernatant into a new vial.

ASSAY PROCEDURE

1. Label conical polystyrene tubes in duplicates (8x2): A to E for the calibrators, NSB for the blank tubes, MB for the maximum binding tubes and T for the total activity tubes. Label additional tubes in duplicate for controls and saliva samples.
- 2a. Pipet 500 µl of incubation buffer into the NSB tubes and 400 µl of incubation buffer into the MB tubes.
- 2b. Pipet 400 µl of the calibrators A to E into the corresponding tubes.
- 2c. Pipet 400 µl of the controls and saliva samples into each of the correspondingly marked tubes.
3. Add 100 µl of the melatonin antibody to all tubes **except** the NSB and T tubes.
4. Add 100 µl of the ¹²⁵I-melatonin tracer to all tubes. Vortex. Remove the T tubes, they will need no further processing until counting at step 10.
5. Incubate all tubes for 20± 4 hours at 2-8°C.
- 6a. Invert the bottle containing the solid phase second antibody several times, add a stir bar and resolve the sediment using a magnetic stirrer.
- 6b. While stirring the second antibody suspension continuously, add 100 µl of the suspension to all assay tubes (except the T tubes). Vortex.
7. Incubate for 15± 1 minutes at 2-8°C.
8. Add 1 ml of cold, deionized water (2-8°C) to all assay tubes (except the T tubes).
9. Centrifuge for 2 minutes at 2000 x g and 2-8°C. Aspirate the supernatants (except the T tubes) and retain the precipitates for counting.
10. Count the tubes for 2 minutes in a gamma-counter.

RESULTS & STANDARDIZATION

1. Record the cpm for all tubes (T, NSB, MB, Calibrators A, B, C, D, E, samples and controls) and calculate the mean cpm for each pair of tubes.
2. Subtract the mean assay blank (NSB tubes) from the respective mean of each pair of tubes:
$$\text{Net cpm} = \text{cpm}_{\text{Average}} - \text{cpm}_{\text{Average NSB}}$$
3. Calculate the binding of each pair of tubes as a percent of maximum binding (MB tubes), with the NSB-corrected cpm of the MB tubes taken as 100%:
$$\text{percentbound} = \frac{\text{netcpm}}{\text{netMBcpm}} \times 100$$
4. Prepare a lin/log graph paper and plot the percent bound on the vertical axis against the melatonin concentration (pg/ml) on the horizontal axis for each of the calibrators. Draw the best fitting curve or calculate the standard curve using a four parameter algorithm or equivalent.
5. Determine the melatonin concentrations for the patient controls and samples from this standard curve. Alternative data reduction methods are equally acceptable.

See Table 11 and Figure 1 for examples of a standard curve. This standard curve is for demonstration purposes

only. A standard curve must be generated for each set of samples to be assayed.

Standardization: Direct Saliva MELATONIN RIA is calibrated with UV/VIS: $\epsilon_{278} = 6300 \text{ M}^{-1}\text{cm}^{-1}$ in methanol.

QUALITY CONTROL

A thorough understanding of this instruction for use (IFU) is necessary for the successful use of the product. Reliable results will be obtained only by using precise laboratory techniques (current GLP guidelines) and accurately following this instruction for use.

The reproducibility of standard curve parameters and control values should be within established limits of laboratory acceptability. The confidence limits for the Controls are lot-specific and printed on the QC Data Sheet added to the kit. If the precision of the assay does not correlate with the established limits and repetition excludes errors in technique, check the following issues: i) pipetting, temperature controlling and timing devices ii) expiration dates of reagents iii) storage and incubation conditions iv) purity of water.

LIMITATIONS

Melatonin results should be interpreted in conjunction with information available from clinical assessment of the patient and other diagnostic procedures.

PERFORMANCE CHARACTERISTICS

The assay performance characteristics have been validated in duplicates.

Intra-Assay Precision (Within-Run): 2.6-20.1%. The intra-assay precision was calculated from results of 20 pairs of values from each saliva sample from different daytime in a single run. The values are presented in Table 12.

Inter-Assay Precision (Run-to-Run): 6.6-16.7%. The inter-assay precision was calculated from results of 20 pairs of values from each saliva sample from different daytime in 20 different runs. The values are presented in Table 13.

Detection Limit (LoB): 0.2 pg/ml (0.9 pmol/l). Twenty zero Calibrator (MB) replicates were assayed in a single run. The minimum detectable concentration of melatonin in 400 µl incubation buffer was calculated by subtracting two standard deviations of averaged zero Calibrator duplicates from the counts at maximum binding and intersecting this value with the standard curve obtained in the same run.

Detection Limit (LoQ): 0.9 pg/ml (4.0 pmol/l). The functional least detectable dose (FLDD) of the assay is the minimal melatonin concentration in saliva that can be measured with an intra-assay coefficient of variation (CV) of less than 10%. The limit of quantification was determined from 13 different saliva samples each measured at least 10 times in a single run.

Dilution Linearity/Parallelism: 93.1%. Two human saliva sample spiked with 40 pg/ml of melatonin were diluted with Incubation Buffer and subsequently assayed according to the assay procedure. The values are presented in Table 14.

Spiking Recovery: 106%. Three saliva samples were spiked with increasing amounts of melatonin and analyzed according to the assay procedure. The values are presented in Table 15.

Specificity: In Table 16 the following cross-reactivities of the melatonin antiserum were found at 50 % binding.

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ED₂₀: 23.7 pg/ml D₅₀: 6.6 pg/ml ED₈₀: 1.47 pg/ml

Table 13 Inter-Assay Precision

| Sample | Mean [pg/ml] | SD [pg/ml] | CV (%) |
|---------------|--------------|------------|--------|
| daytime | 0.82 | 0.137 | 16.7 |
| evening | 8.99 | 0.593 | 6.6 |
| nighttime | 25.58 | 2.148 | 8.4 |
| early morning | 3.39 | 0.256 | 7.5 |
| Mean | | | 9.8 |

Table 14 Dilution Linearity/Parallelism

| Sample | Dilution | Observed [pg/ml] | Expected [pg/ml] | O/E [%] |
|-----------------|------------------|------------------|------------------|---------|
| spiked Saliva 1 | spiked: 40 pg/ml | 44.4 | --- | --- |
| | 1 in 2 (50.0%) | 25.1 | 22.2 | 112.9 |
| | 1 in 4 (25.0%) | 11.3 | 11.1 | 101.6 |
| | 1 in 8 (12.5%) | 4.92 | 5.55 | 88.6 |
| | 1 in 16 (6.3%) | 2.60 | 2.78 | 93.5 |
| | 1 in 32 (3.1%) | 1.51 | 1.39 | 108.6 |
| | 1 in 64 (1.6%) | 0.60 | 0.69 | 87.0 |
| spiked Saliva 2 | spiked: 40 pg/ml | 47.2 | --- | --- |
| | 1 in 2 (50.0%) | 22.6 | 23.6 | 95.6 |
| | 1 in 4 (25.0%) | 10.8 | 11.8 | 91.2 |
| | 1 in 8 (12.5%) | 4.97 | 5.9 | 84.2 |
| | 1 in 16 (6.3%) | 2.65 | 2.95 | 89.9 |
| | 1 in 32 (3.1%) | 1.23 | 1.48 | 83.5 |
| | 1 in 64 (1.6%) | 0.60 | 0.74 | 81.1 |
| Mean | | | | 93.1 |

Figure 1 Example of Standard curve


Table 12 Intra-Assay Precision

| Sample | Mean [pg/ml] | SD [pg/ml] | CV (%) |
|---------------|--------------|------------|--------|
| daytime | 0.60 | 0.121 | 20.1 |
| evening | 3.56 | 0.145 | 4.1 |
| nighttime | 24.42 | 1.169 | 4.8 |
| early morning | 7.24 | 0.188 | 2.6 |
| Mean | | | 7.9 |

Table 15 Spiking Recovery

| Basic Value | Added [pg/ml] | Expected [pg/ml] | Observed [pg/ml] | Recovery [%] |
|-------------|---------------|------------------|------------------|--------------|
| 0.69 | 1 | 1.69 | 1.68 | 99 |
| | 2 | 2.69 | 2.50 | 93 |
| | 5 | 5.69 | 5.82 | 102 |
| | 10 | 10.69 | 11.94 | 112 |
| | 20 | 20.69 | 25.03 | 121 |
| | 40 | 40.69 | 42.95 | 106 |
| 1.35 | 1 | 2.35 | 2.51 | 107 |
| | 2 | 3.35 | 3.45 | 103 |
| | 5 | 6.35 | 7.39 | 116 |
| | 10 | 11.35 | 12.12 | 107 |
| | 20 | 21.35 | 25.11 | 118 |
| | 40 | 41.35 | 45.11 | 109 |
| 0.74 | 2.5 | 3.24 | 3.31 | 102 |
| | 5 | 2.74 | 7.02 | 122 |
| | 10 | 10.74 | 10.79 | 100 |
| | 20 | 20.74 | 21.20 | 102 |
| | 40 | 40.74 | 36.88 | 91 |
| | 80 | 80.74 | 82.60 | 102 |
| Mean | | | | 106 |

Table 16

Specificity

| Compound | Crossreactivity [%] |
|------------------------------|---------------------|
| Melatonin | 100 |
| 6-Sulfatoxymelatonin | 0.002 |
| Serotonin | <0.001 |
| 5-Hydroxy-Indole Acetic Acid | <0.001 |
| N-Acetylserotonin | 0.027 |
| 5-Methoxytryptamine | 0.003 |
| 5-Methoxytryptophane | 0.001 |
| 5-Methoxytryptophol | 0.001 |
| Methoxy psoralen | <0.001 |

Table description: cf. "Results" and "Performance Characteristics" (page 4).

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REFERENCES/ LITERATURREFERENZEN/ RÉFÉRENCES/ RIFERIMENTI/ REFERENCIAS




1. Koorengevel, KM, *et al.* *A forced desynchrony study of circadian pacemaker characteristics in seasonal affective disorder.* J Biol Rhythms **17**, 463-75 (2002).
2. Wirz-Justice, A, *et al.* *No evidence for a phase delay in human circadian rhythms after a single morning melatonin administration.* J Pineal Res **32**, 1-5 (2002).
3. Graw, P, *et al.* *Early morning melatonin administration impairs psychomotor vigilance.* Behav Brain Res **121**, 167-72 (2001).
4. Danilenko, KV, *et al.* *Is sleep per se a Zeitgeber in humans.* J Biol Rhythms **18**, 170-8 (2003).
5. Vaughan G M: *New sensitive serum melatonin radioimmunoassay employing the Kennaway G280 antibody: Syrian hamster morning adrenergic response.* J Pineal Res **15**, 88-103 (1993).
6. Voultsios, A, *et al.* *Salivary melatonin as a circadian phase marker: validation and comparison to plasma melatonin.* J Biol Rhythms **12**, 457-466 (1997)

| RADIOIMMUNOASSAY PROCEDURE | | | | | | | |
|---------------------------------------|-------------------------------|---------------------------------------|-----------------------|--------------------|--|-----------------------------|---|
| Polystyrene tubes in duplicate | Incubation Buffer (µl) | Standard, Control, Sample (µl) | Antiserum (µl) | Tracer (µl) | | Second Antibody (µl) | |
| Total | -- | -- | -- | 100 | | -- | Vortex and incubate for 15±2 minutes at 2-8°C |
| NSB | 500 | -- | -- | 100 | | 100 | |
| MB | 400 | -- | 100 | 100 | | 100 | |
| Std A 0.5 pg/ml | -- | 400 | 100 | 100 | | 100 | add 1 ml of deionized water (except T tubes) and centrifuge for 2 minutes at 2-8°C and 2000 x g |
| Std B 1.5 pg/ml | -- | 400 | 100 | 100 | | 100 | |
| Std C 5.0 pg/ml | -- | 400 | 100 | 100 | | 100 | |
| Std D 15.0 pg/ml | -- | 400 | 100 | 100 | | 100 | |
| Std E 50.0 pg/ml | -- | 400 | 100 | 100 | | 100 | |
| Control LOW | -- | 400 | 100 | 100 | | 100 | |
| Control HIGH | -- | 400 | 100 | 100 | | 100 | Vortex and incubate at 2-8°C for 20±4 hours |
| Samples | -- | 400 | 100 | 100 | | 100 | aspirate supernatant (except T tubes) and count for 2 minutes |

Table 17

SYMBOLS/SYMBOLE/ SYMBOLES/SIMBOLI/ SIMBOLOS

| Symbol | Explanation |
|---|--|
|  | Use By Verwendbar bis Utiliser jusqu'au Utilizzare entro Fecha de caducidad |
|  | Catalogue number Bestellnummer Référence du catalogue Numero di catalogo Número de catálogo |
|  | Batch code Chargenbezeichnung Code du lot Codice del lotto Codigo de lote |
|  | <i>In Vitro</i> Diagnostic Medical Device <i>In Vitro</i> Diagnostikum Dispositif médical de diagnostic <i>in vitro</i> Dispositivo medico-diagnostico <i>in vitro</i> Producto sanitario para diagnóstico <i>in vitro</i> |
|  | Consult Instructions for Use- Gebrauchsanweisung beachten Consulter le mode d'emploi Consultare le istruzioni per l'uso Consulte las instrucciones de uso |
|  | Contains sufficient for <n> tests Ausreichend für „n“ Ansätze Contenu suffisant pour „n“ tests Contenuto sufficiente per „n“ saggi Contenido suficiente para <n> ensayos |
|  | Temperature limitation Zulässiger Temperaturbereich Limites de température Limiti di temperatura Limite de temperatura |
|  | Radioactive Material Radioaktives Material matériel radioactif materiale radioattivo material radiactivo |

| Symbol | Explanation |
|--|--|
|  | Incubation Buffer Inkubations-Puffer Tampon d'incubation Tampone d'incubazione Tampón de incubación |
|  | Antiserum Antiserum Antisérum Antisiero Antisuero |
|  | Tracer Tracer Traceur Elemento tracciante Trazador |
|  | Calibrator A - E Kalibrator A - E Calibrateur A - E Calibratore A - E Calibrador A - E |
|  | Control Low Kontrolle tief Contrôle bas Controllo basso Control bajo |
|  | Control High Kontrolle hoch Contrôle élevé Controllo alto Control alto |
|  | 2 nd Antibody 2. Antikörper 2 ^{ème} Anticorps Secondo anticorpo Segundo anticuerpo |

CE

IVD

Printing Date
2013-01-14