# Quantum Blue<sup>®</sup> Calprotectin Ascites

**Quantitative  
Lateral Flow Assay**

LF-ASC25      25 tests

Revision date: 2013-03-28

## ENGLISH

### INTENDED USE

The Quantum Blue® Calprotectin Ascites, LF-ASC25, is a quantitative test designed for the determination of elevated Calprotectin concentrations in human ascites in combination with the Quantum Blue® Reader. The test has been designed for point of care testing for patients suffering from liver cirrhosis or other diseases accompanied by ascites accumulation in the peritoneum. It is used as a marker for elevated polymorphonuclear cell count (PMN).

It is used as an aid to diagnose the development of spontaneous bacterial peritonitis (SBP) (1).

For professional use only\*.

\*Canada, Taiwan: For laboratory use only.

### PRINCIPLE OF THE ASSAY

The test is designed for the selective measurement of Calprotectin (MRP8/14) antigen by sandwich immunoassay. A monoclonal capture antibody (mAb) highly specific for Calprotectin is coated onto the test membrane. A second monoclonal detection antibody conjugated to gold colloids is deposited onto the conjugate release pad and released into the reaction system after addition of the diluted ascites sample. The Calprotectin/anti-Calprotectin gold conjugate binds to the anti-Calprotectin antibody coated on the test membrane (test line; test band) and the remaining free anti-Calprotectin gold conjugate binds to the goat anti-mouse antibody coated on the test membrane (control line; control band). The signal intensities of the test line and the control line are measured quantitatively by the BÜHLMANN Quantum Blue® Reader.

### REAGENTS SUPPLIED AND PREPARATION

Reagents	Quantity	Code	Comments
<b>Test Cartridge</b> vacuum-sealed in a foil bag	25 pieces	B-CAL-TC	
<b>Chase Buffer</b>	1 bottle 10 ml	B-ASC-CB	Ready to use
<b>RFID Chip Card</b>	1 piece	B-ASC-RCC	White plastic card

Table 1

### STORAGE AND SHELF LIFE OF REAGENTS

All kit components are stable at 2-8°C until the expiration date printed on the labels.

### PRECAUTIONS

#### SAFETY PRECAUTIONS

- None of the reagents of this test contains components of human origin.
- Patient specimens should be handled as if capable of transmitting infections and should be handled in accordance with good laboratory practices using appropriate precautions.
- Unused solution should be disposed of according to local State and Federal regulations.

#### TECHNICAL PRECAUTIONS

##### Kit components

- All reagents and test samples must be equilibrated at room temperature (18-28°C) before starting the assay. Mix well (vortex) the reagents before use.
- Components must not be used after the expiry date printed on the labels.
- Do not mix different lots of reagents.
- Test Cassettes cannot be re-used.

### Test Procedure

- Read carefully the instructions prior to carrying out the test. Test performance will be adversely affected, if reagents are incorrectly diluted, handled or stored under conditions other than those as detailed in this instruction for use.
- Use the white RFID Chip Card in order to change lot-specific test parameters.
- The Quantum Blue® Reader must be switched on and programmed for the Quantum Blue® Calprotectin Ascites assay before starting the assay (see Quantum Blue® Reader Manual).
- Patient samples that are not properly handled may cause inaccurate results.
- Diluted samples should be used within several hours and cannot be stored for a longer time period.
- Undiluted ascites samples can be stored at < -20°C for at least 3 months.

### MATERIALS REQUIRED BUT NOT PROVIDED

- Quantum Blue® Reader available from BÜHLMANN (order code: BI-POCTR-ABS)
- Precision pipettes with disposable tips: 10-200 µl
- Centrifuge for microtubes
- Eppendorf tubes for dilution of the samples
- Timer
- Gloves and Lab coat

### SPECIMEN COLLECTION AND STORAGE

Collect ascites samples into clean tubes and store refrigerated at 2-8°C for up to 7 days.

For longer storage freeze the samples at <-20°C. The samples are stable for at least 3 months.

**Important:** The sample must be collected without any chemical or biological additions in the collection device.

### ASSAY PROCEDURE

#### Sample Dilution

- Dilute ascites samples prior to using them for analysis 1:5 with Chase Buffer (e.g. 20 µl sample and 80 µl Chase Buffer).
- Vortex.

#### Lateral Flow Assay Procedure and Readout

There are two alternative methods available on the Quantum Blue® Reader: <ASC\_720> and <ASC\_0>. Select one of these methods before starting the experiments.

Load the lot specific parameters from the RFID Chip Card.

##### 1. Method <ASC\_720> with internal timer

- Load the test cassette onto the test cassette holder of the Reader
- Add 60 µl of diluted ascites sample onto the sample loading port of the test cassette
- Close the cassette holder and start the measurement by pressing the start button
- The scan starts automatically after 12 minutes (720 seconds).

##### 2. Method <ASC\_0> without internal timer

- Add 60 µl of diluted ascites sample onto the sample loading port of the test cassette
- Incubate for 12 minutes +/- 1 minute (set a timer manually)
- Load the test cassette onto the test cassette holder of the Reader

- Scan the cassette with the Quantum Blue® Reader by pressing the start (<ENTER>) button immediately.

**Remark:** Please refer to the Quantum Blue® Reader Manual to learn about the basic functions and how to initialize and operate the Reader, especially how to select test methods, and how to load lot specific parameters from the RFID Chip Card in order to get the samples measured.

#### QUALITY CONTROL

- If the precision of the assay does not correlate with the established limits and repetition excludes errors in technique, check the following issues: i) pipetting, temperature controlling and timing devices, ii) expiration dates of reagents and iii) storage and incubation conditions.
- Result of the self-test of the Quantum Blue® Reader performed at startup of the instrument has to be valid.

#### VALIDATION OF RESULTS

- For a valid test result, the Control Line (C) must be visible in any case (see Appendix I, Figures 1A and 1B). It is used as functional test control only and cannot be used for the interpretation of the Test Line (T). If the Test Line (T) is not detectable after 12 minutes of incubation time (Figure 1A), the concentration of Calprotectin present in the ascites sample is below the detection limit. If a Test Line (T) is detectable after 12 minutes of incubation time (Figure 1B), the Calprotectin concentration present in the ascites sample is calculated by the Quantum Blue® Reader.
- If only the Test Line (T) is detectable after 12 minutes of incubation time (Figure 1C), the test result is invalid and the Calprotectin assay has to be repeated using another test cassette.
- If neither the Control Line (C) nor the Test Line (T) are detectable after 12 minutes of incubation time (Figure 1D), the test result is invalid and the assay has to be repeated using another test cassette.
- As the Quantum Blue® Reader allows for a quantitative evaluation of the Test (T) and Control (C) Lines, an additional validity check of the Control Line (C) is undertaken. If the signal intensity of the Control Line (C) is below a lot-specific threshold after 12 minutes of incubation time, the test result is also invalid and the Calprotectin Ascites assay has to be repeated using another test cassette.

#### STANDARDIZATION AND INTERPRETATION OF RESULTS

- The Lateral Flow Assay is calibrated with the BÜHLMANN Calprotectin ELISA (order code: EK-CAL).
- The BÜHLMANN Quantum Blue® Reader uses a lot-specific standard curve to calculate the Calprotectin Ascites concentration. This lot-specific standard curve is generated with the median values ( $n \geq 10$  measurements each) from  $\geq 10$  calibration points obtained from different samples with known Calprotectin Ascites concentrations. The assay range is between 0.18 and 1.8 µg/mL.
- For quantitative measurements of samples reading above 1.8 µg/mL dilute the ascites samples 1:50 (instead of 1:5) with chase buffer and test them again according to the procedure. The measured

concentration must then be multiplied by the volume factor x 10 to obtain the final result.

#### LIMITATIONS

- The reagents supplied with this kit are optimized to measure human Calprotectin in ascites samples.
- The BÜHLMANN Quantum Blue® Calprotectin Ascites rapid test can be used to obtain results earlier than PMN count. It should not replace PMN count but can aid in establishing the diagnosis. However, a negative result does not exclude disease progression to spontaneous bacterial peritonitis, SBP.
- Results obtained with the BÜHLMANN Quantum Blue® Calprotectin Ascites rapid assay should be confirmed by measuring other biomarkers (such as counting polymorphonuclear cells PMN in ascites) and the anamnesis of the patient, before any decision concerning therapy or medical intervention is being made.
- Ascites samples contaminated with erythrocytes or hemoglobin might cause false positive results and should not be used.

#### RESULTS & INTERPRETATION

The Calprotectin ascites rapid test helps quickly establish a first estimate if an elevated PMN count is to be expected and thus could be helpful to the clinician.

A clinical pilot study (1), using the BÜHLMANN Quantum Blue® rapid test showed that it is useful to predict levels.

According to data shown by Burri et al. (1) we recommend the following reference intervals for interpretation of calprotectin results:

PMN count  $\leq 250$  cells ( $n=111$ ):

Median: 0.38 µg/mL,  
IQR 0.38-0.562.

PMN count  $\geq 250$  cells ( $n=19$ ):

Median: 2.78 µg/mL,  
IQR 2.05-5.37.

At a cut of 0.51 µg/mL calprotectin, a sensitivity of 100% and a specificity of 84.7% was reached, with 6.53 LR+ and 0.0 LR-. The overall accuracy des Testes was 87.7%. NPV of calprotectin testing in ascites was 100%. In this study, none of the patients with elevated PMN count would have been missed by the rapid test (refer to Figure 3).

Recommended interpretation of calprotectin results:

<0.51 µg/mL: PMN counts of > 250 are rarely possible.

$\geq 0.51$  µg/mL: PMN counts of > 250 potentially occur.

#### PERFORMANCE CHARACTERISTICS

**Limit of Blank (LoB):** <0.044 µg/mL calprotectin.

**Limit of Detection (LoD):** <0.09 µg/mL calprotectin.

**Repeatability: 19% CV.** The repeatability of the Quantum Blue® Calprotectin assay was calculated from 6 ascites samples containing 0.18 and 1.6 µg/g calprotectin. Each sample was tested according to the assay procedure in one run of 20 replicates set up within 10 minutes. The mean values of three different lots of Test Cassettes are presented in (Table 6). The repeatability varied between 13.7-34.5 % CV.

**Linearity: 0.18 to 1.8 µg/mL.** Four ascites samples with elevated Calprotectin concentrations were diluted with an ascites sample with a low level of calprotectin. Each dilution was subsequently assayed according to the assay procedure. The results showed linearity within the indicated measuring range of 0.18 to 1.8 µg/mL of the Quantum Blue® Calprotectin Ascites assay for all 4 samples (R2 = 0.986 - 0.995).

**Method Comparison: R<sup>2</sup> = 0.82; y=1.001x-0.090 µg/mL**  
49 ascites samples within the indicated measuring range of the Quantum Blue® Calprotectin Ascites assay were analyzed according to the assay procedure and compared with the values obtained by the BÜHLMANN Calprotectin ELISA. The correlation data are illustrated in Figure 4.

**Limit of Quantification (LoQ): Lower LoQ: ≤0.18 µg/mL calprotectin; Upper LoQ: ≥1.8 µg/mL calprotectin.** The LoQ has been established with twelve ascites samples at concentrations between 0.12 and 2.46 µg/mL. The samples were measured with 20 replicates each. The resulting precision profile is shown in Figure 2. The limit of quantification corresponds to the concentration of Calprotectin with an imprecision below 25% CV allowing a quantitative measurement within the range from 0.18 (lower LoQ) to 1.8 µg/mL (upper LoQ).

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
Figure 1D

Figure 3

Calprotectin in ascites and PMN count

Box Plot: Calprotectin levels in ascites, classified according to the patients PMN count: >250 or < 250 cells/ml.

Table 6

Precision (Repeatability)

Sample	Concentration on mean [µg/mL]	Repeatability [CV]	Total Precision [CV]
ASC1	0.19	16.7%	18.8%
ASC2	0.32	12.9%	13.7%
ASC3	0.53	11.1%	18.8%
ASC4	0.75	21.0%	23.7%
ASC5	1.28	26.2%	26.2%
ASC6	1.60	32.7%	34.5%
<b>Mean</b>		<b>20.1%</b>	<b>22.6%</b>

Figure 2

Precision Profile



Figure 4

Method Comparison

## REFERENCES/ LITERATURREFERENZEN/ REFERENCES/ RIFERIMENTI/ REFERENCIAS

1. Burri E. et al: *Measurement of calprotectin in ascitic fluid to identify elevated polymorphonuclear cell count*; World J. Gastroenterol., 2013 April 7; 19(13):2028-2036
2. Homann C. et al.: *Ascites fluid and plasma calprotectin concentrations in liver disease*. Scand J Gastroenterol 2003; **38**(4): 415-420.
3. Erwin Biecker et al.: *Diagnosis and therapy of ascites in liver cirrhosis*; World J Gastroenterol 2011 March 14; 17(10): 1237-1248.
4. Garcia-Tsao G. *Current management of the complications of cirrhosis and portal hypertension: variceal hemorrhage, ascites, and spontaneous bacterial peritonitis*. Gastroenterology 2001; 120(3): 726-748
5. D'Amico G, Garcia-Tsao G, Pagliaro L. *Natural history and prognostic indicators of survival in cirrhosis: a systematic review of 118 studies*. J Hepatol 2006; 44(1): 217-231
6. Rimola A, et al. *Diagnosis, treatment and prophylaxis of spontaneous bacterial peritonitis: a consensus document*. International Ascites Club. J Hepatol 2000; 32(1): 142-153.
7. Runyon BA. *Management of adult patients with ascites due to cirrhosis: an update*. Hepatology 2009; **49**(6): 2087-2107.
8. Gines P et al.: *EASL clinical practice guidelines on the management of ascites, spontaneous bacterial peritonitis, and hepatorenal syndrome in cirrhosis*. J Hepatol 2010; **53**(3): 397-417
9. Ebanna et al.: *Fluid level of Calprotectin in chronic liver disease - Malignant and Non Malignant*; Bull. Alex. Fac. Med. 44 No.3, 2008

## APPENDIX III

## SYMBOLS/ SYMBOLE/ SYMBOLES/ SIMBOLI/ SIMBOLOS

Symbol	Explanation
	Use By Verwendbar bis Utiliser jusqu'au Utilizzare entro Fecha de caducidad
	Catalogue number Bestellnummer Référence du catalogue Numero di catalogo Número de catálogo
	Batch code Chargenbezeichnung Code du lot Codice del lotto Codigo de lote
	<i>In Vitro</i> Diagnostic Medical Device <i>In Vitro</i> Diagnostikum Dispositif médical de diagnostic <i>in vitro</i> Dispositivo medico-diagnostico <i>in vitro</i> Producto sanitario para diagnóstico <i>in vitro</i>
	Contains sufficient for <n> tests Ausreichend für "n" Ansätze Contenu suffisant pour „n“ tests Contenuto sufficiente per „n“ saggi Contenido suficiente para <n> ensayos
	Consult Instructions for Use- Gebrauchsanweisung beachten Consulter le mode d'emploi Consultare le istruzioni per l'uso Consulte las instrucciones de uso

Symbol	Explanation
	Temperature limitation Zulässiger Temperaturbereich Limites de température Limiti di temperatura Limite de temperatura
	Test Cassette Test Kassette Cartouche de testi Card di rilevazione Cartucho de prueba
	Chase Buffer Lauf-Puffer Tampon de chase Tampone di chase Tampón de chase
	RFID Chip card RFID Chipkarte RFID Carte à puce RFID Carta chip RFID Tarjeta chip