# VASOPRESSIN

## RIA

RK-AR1 100 tests

Revision date: 2012-09-12

**INTENDED USE**

This double antibody radioimmuno-assay is designed for the quantitative *in vitro* diagnostic determination of immunoreactive **arginine vasopressin** (anti-diuretic hormone, ADH) in EDTA plasma and urine after extraction (1-3).

**PRINCIPLE OF THE ASSAY**

Immunoreactive [Arg8] vasopressin is measured by a double-antibody radioimmunoassay using a modification of the method of Glick and Kagan (4). Column (or ethanol) extracted plasma samples and calibrators are pre-incubated for 24 hours with the anti-vasopressin antibody; <sup>125</sup>I-vasopressin competes then with vasopressin present in samples and calibrators for the same antibody binding sites. After a second 24 hour incubation, the solid phase second antibody is added to the mixture, the antibody-bound fraction is finally precipitated and counted.

The procedure recommends a solid-phase extraction of the plasma samples with reversed-phase columns (phenylsilylsilica). Alternatively, an extraction of plasma samples with ethanol may also be used. However, the values using the ethanol procedure are somewhat higher due to unspecific matrix effects (see Page 3).

**REAGENTS SUPPLIED AND PREPARATION**

Reagents	Quantity	Code	Reconstitution
<b>Extraction Columns (1ml)</b> Reversed phase columns containing 100 mg phenylsilylsilica	10 pieces	B-AEC	For column preparation and conditioning see chapter "Extraction"
<b>Phosphate Buffer</b> lyophilized Buffer	1 vial	B-AR1-PB	Reconstitute with 100 ml of deionized water
<b>Antiserum</b> lyophilized rabbit anti-vasopressin antibody	1 vial	B-AR1-AS	Reconstitute with 10 ml of deionized water
<b>Tracer</b> lyophilized <sup>125</sup> I-Vasopressin	1 vial	B-ADH-TR	Reconstitute with 11 ml of Phosphate Buffer
<b>Calibrator<sup>1)</sup></b> Lyophilized synthetic arginine vasopressin in Phosphate Buffer	1 vial	B-AR1-CA	Reconstitute with 5 ml of deionized water
<b>Controls Normal / High<sup>2)</sup></b> Arginine Vasopressin in a buffer matrix	2 vials	B-AR1-CONSET	Reconstitute with 5 ml of Phosphate Buffer
<b>Second Antibody</b> Cellulose coated anti-rabbit antibody	1 vial 11 ml	B-AB2	Ready to use

Table 1

<sup>1)</sup> Reconstitution of the Calibrator results in a stock solution of 80 pg/ml arginine vasopressin.

<sup>2)</sup> Lot specific amounts of arginine vasopressin in buffer matrix. Refer to the additional QC Data Sheets for exact concentrations.

**Note:** DO NOT EXTRACT NEITHER CONTROLS NOR CALIBRATORS.

**STORAGE AND SHELF LIFE OF REAGENTS**

Unopened Reagents	
The Extraction columns are stable at 18-28°C. The other kit components are stable at 2-8°C. Do not use past kit expiration date printed on the labels. Do not freeze the Second Antibody.	
Opened / Reconstituted Reagents	
Extraction columns	Used columns should be stored at 18-28°C and protected from light and dust.
Phosphate Buffer	Stable for 2 months at 2-8°C
Antiserum	Stable for 2 months at -20°C
Tracer	Store at -20°C. Aliquot if repeated use is expected
Calibrator	Stable for 2 months at -20°C
Controls	Aliquot if repeated use is expected
Second Antibody	Store at 2-8°C (do not freeze)

Table 2

**SAFETY PRECAUTIONS**

- **Radioactive Material:** This kit contains radioactive material which does not exceed 56 kBq of <sup>125</sup>Iodine.
- Receipt, acquisition, possession, use, and transfer are subject to the local regulations. Unused solutions and radioactive waste should be disposed of according to local State and Federal regulations.
- All kit reagents except of the extraction columns (B-AEC) and the second antibody (B-AB2) contain components of human origin. Although tested and found negative for HBV surface antigen, HCV and HIV1/2 antibodies, the reagents should be handled as if capable of transmitting infections and should be handled in accordance with good laboratory practices using appropriate precautions.

**TECHNICAL PRECAUTIONS**

- Read carefully the instructions prior to carrying out the test. Test performance will be adversely affected if reagents are incorrectly diluted, modified or stored under conditions other than those as detailed in this instruction for use.
- Components must not be used after the expiry date printed on the labels.
- Do not mix different lots of reagents.
- Let the reagents adjust to reach room temperature. Reconstitute the lyophilized reagents as indicated. Mix well (vortex) the reagents before use.
- Counting time should be selected in order to keep statistical counting error small: e.g., at 2000 cpm the counting error is at 5%; at 10000 cpm it is only 1%.
- If the initial concentration of an unknown sample reads above the highest calibrator, the sample should be diluted with phosphate buffer and tested again according to the assay procedure.

**MATERIALS REQUIRED BUT NOT PROVIDED**

- 100 µl, 400 µl, 500 µl, 1000 µl, and 5000 µl precision pipettes (or preferably a 100 - 1000 µl adjustable multipipette) with disposable tips.
- 10 ml volumetric pipette.
- Disposable polystyrene tubes for the preparation of the calibrator dilutions.
- Disposable conical polystyrene tubes to run the assay (e.g. Sarstedt # 57.477).
- Disposable polypropylene, polystyrene or glass tubes for the preparation of plasma extracts.
- Cylinder for the preparation of phosphate buffer.
- Deionized or distilled water.
- Methanol p.a. (e.g. Merck # 6009).
- 1 N HCl (hydrochloric acid) (e.g. Merck # 9057).
- Extraction vacuum manifold for applying the extraction columns (optional).
- Vacuum concentrator or supply of particle free nitrogen.
- Refrigerated centrifuge.
- Vortex rotator.
- Stir bar and magnetic stirrer.
- Aspiration device.
- Gamma counter.

**SPECIMEN COLLECTION AND STORAGE**

Appropriate sample collection is essential to ensure accurate results of the vasopressin analysis. If the procedure calls for true basal levels, the patient must be fasting for at least 12 hours and must stay recumbent, without any stress and in a quiet environment, for at least 1 hour prior to blood collection.

Collect at least 5 ml blood into an **EDTA venipuncture tube** and immediately place the sample on ice. Centrifuge at 2-8°C at 2000 x g for 15 minutes within 10 minutes of blood collection. Separate the plasma from the cells and freeze the specimen immediately at -20°C or precede to the extraction procedures (see below).

Vasopressin **in urine** remains stable for at least 2 months at -20°C, if such samples have been acidified with 10% (v/v) of 1 N HCl (hydrochloric acid) immediately after collection (6). After a short centrifugation step, these samples may be added directly to the conditioned extraction columns (see below). The procedure calls for 1 ml of EDTA plasma or urine (duplicate determinations).

### EXTRACTION

We recommend an extraction procedure using a reversed-phase column extraction which is highly specific for the adsorption and the subsequent elution of arginine vasopressin. Alternatively, an ethanol extraction method may also be used.

- The extraction columns provided with this kit can **each be utilized up to five times** if used according to the extraction procedures described in this protocol. Used columns should be stored at 18-28°C and protected from light and dust.
- The column extraction methods described below result in **recoveries** of more than 90% with <sup>125</sup>I-vasopressin or <sup>3</sup>H-vasopressin, respectively, spiked into human EDTA plasma samples.
- If samples have to be measured containing very low concentrations of arginine vasopressin the sample application volume may be increased up to 4 ml without a notable change in extraction recoveries (see Extractive Concentration in the assay performance below).
- The column extraction procedures were **tested and validated for human EDTA plasma** and **urine** samples. If other specimens such as **animal plasma** are used, it is recommended to validate the extraction recovery using a spike of <sup>125</sup>I-vasopressin in the specimen. The high ionic strength in **urine** due to high salt concentrations might interfere with the radioimmunoassay for vasopressin. Therefore, column extraction of urine samples after acidification with HCl is highly recommended (6). However, Panzali *et al.* (8) have assayed 17 urine samples without column extraction and found results close to those obtained with extracted samples.

### Extraction Procedure using Centrifugation

#### Sample Pretreatment

- Mark one polypropylene, polystyrene or glass tube for each sample to be extracted.
- Add 1 ml of corresponding plasma sample.
- Add 100 µl of 1 N HCl, vortex and centrifuge for 5 minutes at 2000 x g.
- Use acidified plasma supernatant in step Sample application

#### Column Preparation and Conditioning

- Mark one extraction column for each sample to be extracted and place into polypropylene, polystyrene or glass centrifugation tubes.
- Add 1 ml of methanol to columns and centrifuge for 1 minute at 200 x g. Repeat this step once.
- Add 1 ml of distilled or deionized water to columns, centrifuge for 1 minute at 200 x g. Repeat this step once.
- Empty tubes to avoid tips of extraction columns from contacting eluates.

#### Sample Application

- Load 1 ml of acidified sample onto the correspondingly marked column and centrifuge for 1 minute at 200 x g.

#### Washing

- Add 1 ml of distilled or deionized water to columns, centrifuge for 1 minute at 500 x g. Repeat this step once.

#### Elution of Extract

- Place each extraction column into a clean correspondingly marked polypropylene, polystyrene or glass tube.
- Add 1 ml of methanol containing 0.5 % (v/v) of 1 N HCl solution to columns and centrifuge for 1 minute at 200 x g.
- Use column for extracting the next sample (up to 5 times) or store column at 18-28°C and protected from light and dust.

#### Evaporation and Reconstitution of Extract

- Evaporate the methanol solution to dryness using a vacuum concentrator with a cold trap. Alternatively, use a 37°C heating block or water bath and evaporate the methanol to dryness with a stream of particle free nitrogen.
- Reconstitute the samples with 1 ml of phosphate buffer, vortex.
- Equilibrate the extracts for 30 minutes at 2-8°C and vortex them from time to time. Store reconstituted extracts capped and frozen at -20°C if not assayed immediately.

### Extraction Procedure Using Vacuum Manifold

With the help of negative pressure the fluids will be passed through the column. The procedure is the same as described in the Extraction procedure using centrifugation.

The following flow rates should be used:

- Sample application and elution should be done with a flow rate of 2 ml/min
- All the other fluids can pass the column with a flow rate of 5 ml/min.

### Ethanol Extraction of Samples

Alternatively, samples may be extracted using "ethanol precipitation".

- Mix 1 ml of each patient sample with 5 ml of chilled 98 % (v/v) ethanol. Vortex for 1 minute.
- Centrifuge for 20 minutes at 2-8°C and 1000 x g. Decant into another clean polypropylene, polystyrene or glass tube.
- Dry the supernatant in a vacuum concentrator. Alternatively, use a 37°C heating block or water bath and evaporate the ethanol to dryness with a stream of compressed and particle free nitrogen or air.
- Reconstitute the samples with 1 ml of phosphate buffer, vortex and store them frozen at -20°C if not assayed immediately.

With the ethanol extraction, samples show usual recoveries of 65-80%. Therefore, we recommend using as an internal recovery control 1 ml of EDTA plasma spiked with 50 µl of <sup>125</sup>I-Vasopressin. Compared to the total activity this factor can be used to calculate the recovery rate.

## ASSAY PROCEDURE

### Preparation of the standard curve:

In order to obtain an entire standard curve, serial dilutions of the Calibrator are prepared as follows:

- Label eight tubes A through H and pipet 1.0 ml of phosphate buffer into tubes B through H.
- Pipet 1.0 ml of the reconstituted Calibrator stock solution (80 pg/ml) into tubes A and B, vortex.
- Transfer 1.0 ml from tube B to tube C, vortex.
- Continue to transfer 1.0 ml from each tube until dilution series is completed. The corresponding concentrations of arginine vasopressin will be:

A 80 pg/ml (73.6 pmol/L)	E 5.0 pg/ml (4.6 pmol/L)
B 40 pg/ml (36.8 pmol/L)	F 2.5 pg/ml (2.3 pmol/L)
C 20 pg/ml (18.4 pmol/L)	G 1.25 pg/ml (1.1 pmol/L)
D 10 pg/ml (9.2 pmol/L)	H 0.63 pg/ml (0.6 pmol/L)

**Note: Allow all reagents for steps 1 to 4 to come to room temperature (18-28°C) prior to use.**

1. Label 11 polystyrene tubes in duplicate: A to H (Calibrators), NSB (blank), MB (maximum binding) and T (total activity). Label additional polystyrene tubes in duplicate for patient samples and controls.
- 2a. Pipet 500 µl of phosphate buffer into the NSB tubes and 400 µl into the MB tubes.
- 2b. Pipet 400 µl of each of the A to H Calibrators into the corresponding tubes.
- 2c. Pipet 400 µl of the extracted patient samples and controls (without extraction) into each of the corresponding tubes.
3. Add 100 µl of the vasopressin antiserum to all tubes except the NSB and T tubes. Vortex.
4. Incubate for 24 hours ± 3 hours at 2-8°C.
5. Add 100 µl of vasopressin tracer to all tubes. Vortex. Cap and remove the T tubes, they will require no further processing until counting step 12.
6. Incubate for 24 hours ± 3 hours at 2-8°C.
- 7a. Invert several times the bottle containing the solid phase second antibody, add a stir bar and place the bottle on a magnetic stirrer.
- 7b. While stirring continuously, add 100 µl of the suspension to each assay tube (except T tubes). Vortex.
8. Incubate for 20 minutes ± 1 minute at 18-28°C.
9. Add 1 ml of deionized water to each tube (except T tubes) **without resuspending the precipitates.**
10. Centrifuge for 5 minutes at 2-8°C and 1000 x g.
11. Aspirate all supernatants (**except the T tubes**) and retain the precipitates for counting.
12. Count all tubes for 1 minute in a gamma counter.

## INTERPRETATION OF RESULTS & STANDARDIZATION

Record the cpm for all tubes (T, NSB, MB, Calibrators A-H samples and Controls) and calculate the mean cpm for each pair of tubes. Subtract the mean assay blank (NSB tubes) from the respective mean of each pair of tubes:

$$\text{Net cpm} = \text{cpm}_{\text{Average}} - \text{cpm}_{\text{Average NSB}}$$

Calculate the binding of each pair of tubes as a percent of maximum binding (MB tubes), with the NSB-corrected cpm of the MB tubes taken as 100%.

$$\text{percentbound} = \frac{\text{netcpm}}{\text{netMBcpm}} \times 100$$

Prepare a lin/log graph paper and plot the percent bound on the vertical axis against the Vasopressin concentration (pg/ml) on the horizontal axis for each of the Calibrators. Draw the best fitting curve or calculate the standard curve using a four parameter algorithm. Determine Vasopressin concentrations

for the patient samples and Controls from this standard curve. Alternative data reduction methods are equally acceptable.

In order to calculate the concentration in pmol/l, multiply the pg/ml values by a factor of 0.92.

Refer to Table 11 and Figure 1 for examples of results and standard curves. *These results and standard curves are for demonstration purposes only. A standard curve must be generated for each set of samples to be assayed.*

**Standardization:** The BÜHLMANN arginine vasopressin Calibrator is calibrated against WHO reference preparation 77/501.

## QUALITY CONTROL

A thorough understanding of this instruction for use is necessary for the successful use of the product. Reliable results will be obtained only by using precise laboratory techniques (current GLP guidelines) and accurately following this instruction for use.

The accuracy of each actual calibrator lot is controlled every three months by comparing with WHO reference preparation 77/501.

The reproducibility of standard curve parameters and control values should be within established limits of laboratory acceptability. The confidence limits for the Controls are lot-specific and printed on the QC Data Sheet added to the kit.

If the precision of the assay does not correlate with the established limits and repetition excludes errors in technique, check the following issues: i) pipetting, temperature controlling and timing devices ii) expiration dates of reagents iii) storage and incubation conditions iv) purity of water.

## PERFORMANCE LIMITATIONS

- It is mandatory using EDTA plasma **ONLY** in order to inhibit metalloprotease activities.
- Use of conical polystyrene tubes is strongly recommended. During step 11 of the assay procedure, a more solid pellet will be achieved and the following aspiration of the supernatant can be done much easier.
- Patient samples that are not properly collected and handled may cause inaccurate arginine vasopressin results (see Specimen Collection). Patient samples must be frozen immediately in order to ensure correct values at the time of measurement.
- Test results should be interpreted in conjunction with information available from clinical assessment of the patient and other diagnostic procedures.
- Vasopressin is extremely unstable. It is crucial to freeze the samples immediately. Transportation has to be at -20 °C.

## PERFORMANCE CHARACTERISTICS

The assay performance characteristics have been validated in duplicates.

**Intra-Assay Precision (Within-Run): 7.6%.** The intra-assay precision was calculated from the results of 20 pairs of values from 3 column-extracted plasma samples spiked with different amount of arginine vasopressin in a single run. The values are presented in Table 12.

**Inter-Assay Precision (Run-to-Run): 10.0%.** The inter-assay precision was determined from the results of 20 pairs of values from 2 column-extracted and 2 non-extracted plasma samples assayed in 20 consecutive runs. The values are presented in Table 13.

**Detection limit (LoB): 0.39 pg/ml (0.36 pmol/l).** Twenty Zero Calibrator duplicates (phosphate buffer) were assayed in a single run. The Limit of Blank of vasopressin was calculated by subtracting two standard deviations from the mean cpm of

the maximum binding and intersecting this value with the standard curve obtained in the same run.

**Detection limit (LoQ): 1.9 pg/ml** (1.75 pmol/l). Five different serum samples with values between 1.3 and 4.1 pg/ml were tested 10 times each in duplicates in one assay. CV and mean value were calculated for each sample. The least detectable dose (limit of quantification) was observed at 10% CV.

**Dilution Linearity/Parallelism: 95.9%.** A human plasma sample was diluted serially with phosphate buffer before and after column extraction and subsequently assayed according to the assay procedure. The values are presented in Table 14.

**Spiking Recovery: 93.4%.** Three plasma samples were spiked with increasing amounts of synthetic arginine vasopressin, extracted with the column extraction method and analyzed according to the assay procedure. The values are presented in Table 15.

**Extractive Concentration: 96.7%.** Increasing volumes of a human plasma sample containing a low concentration of arginine vasopressin were applied onto extraction columns and extracted according to the protocol (see pages 3). Each of the resulting extract was reconstituted in 1 ml of phosphate buffer and subsequently assayed according to the assay procedure. The values are presented in Table 16.

**Specificity:** Cross-reactions of the vasopressin antiserum were determined at 50% binding and are presented in Table 17.

**Method comparison:** The Vasopressin RIA has been compared with the VASOPRESSIN direct RIA (RK-VPD, without column extraction). Results obtained with 28 blood donors (EDTA plasma) show an excellent correlation (see Figure 2).

#### Reference Intervals

In an evaluation including EDTA plasma samples from 53 apparently healthy blood donors, we have established the **normal value** (95<sup>th</sup> percentile) of: 3.6 pg/ml (see Table 18).

This referencel interval should be used as guideline only. It is recommended that each laboratory establishes its own expected range for its patient population.

Values for blood donors may be slightly elevated due to unknown fasting conditions and unknown position before donation. Therefore, correctly collected basal values are expected to be in the lower range or even undetectable. Most procedures need dynamic testing by stimulation or suppression of vasopressin release.

The clinical validity of the Bühlmann Vasopressin-RIA (RK-AR1) was demonstrated in many clinical and research publications (5-7).

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
Table 11 **Example of Results**  
(Incubation at 28°C)

	cpm	B/T [%]	B/B <sub>0</sub> [%]	Conc [pg/ml]	CV [%]
Total	18212	100.0			
Total	18187	100.0			
<b>Total Avg.</b>	18200	100.0			
NSB	216	1.2			
NSB	210	1.2			
<b>NSB Avg.</b>	213	1.2			
MB	5743	31.6	100.0	0	
MB	6054	33.3	100.0	0	
<b>MB Avg.</b>	5898	32.4	100.0	0	
Calibrator H	5524	30.4	93.4	0.63	
Calibrator H	5490	30.2	92.8	0.63	
<b>Calibrator H Avg.</b>	5507	30.3	93.1	0.63	4.8
Calibrator G	5081	27.9	85.6	1.25	
Calibrator G	5208	28.6	87.8	1.25	
<b>Calibrator G Avg.</b>	5144	28.3	86.7	1.25	11.2
Calibrator F	4365	23.9	74.0	2.5	
Calibrator F	4359	23.9	73.9	2.5	
<b>Calibrator F Avg.</b>	4362	23.9	73.9	2.5	0.5
Calibrator E	3548	19.5	58.7	5	
Calibrator E	3520	19.3	58.1	5	
<b>Calibrator E Avg.</b>	3533	19.4	58.4	5	3.9
Calibrator D	2165	11.9	34.3	10	
Calibrator D	2082	11.4	32.9	10	
<b>Calibrator D Avg.</b>	2124	11.7	33.6	10	3.9
Calibrator C	1509	8.3	22.8	20	
Calibrator C	1348	7.4	20.0	20	
<b>Calibrator C Avg.</b>	1428	7.9	21.4	20	11.4
Calibrator B	826	4.5	10.8	40	
Calibrator B	819	4.5	10.7	40	
<b>Calibrator B Avg.</b>	823	4.5	10.7	40	1.1
Calibrator A	642	3.5	7.5	80	
Calibrator A	596	3.3	6.7	80	
<b>Calibrator A Avg.</b>	619	3.4	7.1	80	16.1
Control Normal	3340		55.0	5.1	
Control Normal	3335		54.9	5.1	
<b>Control Norm Avg.</b>	3337		55.0	5.1	0.2
Control High	1633		25.0	15.4	
Control High	1517		22.9	17.2	
<b>Control High Avg.</b>	1575		24.0	16.3	7.6

ED-20 = 20.3 pg/ml

ED-50 = 6.1 ng/ml

ED-80 = 1.8 pg/ml

Figure 1 **Example of a Standard Curve**Table 12 **Intra-Assay Precision (Within-Run)**

Sample	Mean Value [pg/ml]	SD [pg/ml]	CV [%]
Plasma 1	1.9	0.2	8.5
Plasma 2	4.6	0.3	7.6
Plasma 3	10.1	0.7	6.6
Mean			7.6

Table 13 **Inter-Assay Precision (Run-to-Run)**

Sample	Mean Value [pg/ml]	SD [pg/ml]	CV [%]
Plasma 4	2.3	0.2	11.0
Plasma 5	10.2	0.9	9.1
Plasma 6*	2.6	0.3	11.4
Plasma 7*	16.5	1.4	8.6
Mean			10.0

\*non-extracted samples

Table 14 **Dilution Linearity/Parallelism**

Sample	Dilution	Observed [pg/ml]	Expected [pg/ml]	O/E [%]
Plasma 6 before extraction	1:1	18.9	--	--
	1:2	9.2	9.5	97
	1:4	4.6	4.7	98
	1:8	2.2	2.4	92
	1:16	1.2	1.2	98
Plasma 6 after extraction	1:1	17.6	--	--
	1:2	8.9	8.8	101
	1:4	4.3	4.4	98
	1:8	2.1	2.2	95
	1:16	1.0	1.1	88
Mean				95.9

Table 15 **Spiking Recovery**

Sample	Spiked with Vasopressin	Observed Value [pg/ml]	Expected Value [pg/ml]	Recovery [%]
EDTA Pool1	5 pg/ml	4.9	5.8	84.9
	25 pg/ml	26.0	25.8	100.5
	50 pg/ml	51.7	50.8	101.7
EDTA Pool 2	5 pg/ml	5.3	6.0	86.9
	25 pg/ml	25.8	26.0	98.9
	50 pg/ml	48.6	51.0	95.3
EDTA Pool 2	5 pg/ml	5.1	6.0	84.0
	25 pg/ml	25.0	26.0	96.0
	50 pg/ml	47.2	51.0	92.4
Mean				93.4

Table 16 **Extractive Concentration**

Sample	Sample Load	Observed [pg/ml]	Expected [pg/ml]	O/E [%]
Plasma 9	1 ml	1.8	--	--
	2 ml	3.6	3.6	100
	3 ml	5.4	5.4	100
	4 ml	6.5	7.2	90
Mean				96.7

Table 17 **Specificity**

Peptide	Cross-Reaction [%]
Arginine vasopressin	100
Arg8-Vasotocin	<0.001
<b>Lysine vasopressin*</b>	<b>0.016</b>
Desmopressin (DDAVP)	0.056
Oxytocin	<0.001

\*In Pigs and Hippopotamus the arginine at position eight is replaced by lysine

Figure 2 Method Comparison RK-VPD vs. RK-AR1


Table 18 Expected Values

	[pg/ml]
n	53
Mean	2.2
Median	2.1
SD	0.88
5th - 95th percentile	1.0 - 3.6
Min / Max	0.9 / 4.7
Mean + 2SD	3.96

**Table description:** cf. "Interpretation of RESULTS & Standardization" (page 4), "Performance Characteristics" and "REFERENCE INTERVALS" (page 5).

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**RADIOIMMUNOASSAY PROCEDURE**

Polystyrene tubes in duplicate	Phosph. Buffer (µl)	Standard, Control, Sample (µl)	Antiserum (µl)		Tracer (µl)		Second Antibody (µl)	
Total	--	--	--		100		--	Vortex and incubate for 20 minutes (± 1 min) at 18-28°C
NSB	500	--	--		100		100	
MB	400	--	100		100		100	
Std A 80 pg/ml	--	400	100		100		100	
Std B 40 pg/ml	--	400	100	Vortex and incubate at 2-8°C for 24 hours (± 3hrs)	100		100	Add 1ml of deionized water (except T tubes) and centrifuge for 5 minutes at 2-8°C and 1000 x g
Std C 20 pg/ml	--	400	100		100		100	
Std D 10 pg/ml	--	400	100		100		100	
Std E 5 pg/ml	--	400	100		100		100	
Std F 2.5 pg/ml	--	400	100		100		100	
Std G 1.25 pg/ml	--	400	100		100		100	
Std H 0.63 pg/ml	--	400	100		100		100	
Control NORMAL	--	400	100		100		100	Aspirate supernatant (except t tubes) and count for 1 minute
Control ELEVATED	--	400	100		100		100	
Sample	--	400	100		100		100	

SYMBOLS/ SYMBOLE/ SYMBOLES/ SIMBOLI/ SIMBOLOS

Symbol	Explanation
	Use By Verwendbar bis Utiliser jusqu'au Utilizzare entro Fecha de caducidad
<b>REF</b>	Catalogue number Bestellnummer Référéce du catalogue Numero di catalogo Número de catálogo
<b>LOT</b>	Batch code Chargenbezeichnung Code du lot Codice del lotto Codigo de lote
<b>IVD</b>	<i>In Vitro</i> Diagnostic Medical Device <i>In Vitro</i> Diagnostikum Dispositif médical de diagnostic <i>in vitro</i> Dispositivo medico-diagnostico <i>in vitro</i> Producto sanitario para diagnóstico <i>in vitro</i>
	Temperature limitation Zulässiger Temperaturbereich Limites de température Limiti di temperatura Limite de temperatura
	Consult Instructions for Use- Gebrauchsanweisung beachten Consulter le mode d'emploi Consultare le istruzioni per l'uso Consulte las instrucciones de uso
	Contains sufficient for <n> tests Ausreichend für „n“ Ansätze Contenu suffisant pour „n“ tests Contenuto sufficiente per „n“ saggi Contenido suficiente para <n> ensayos
	Radioactive Material Radioaktives Material Matériel radioactif Materiale radioattivo Material radiactivo

Symbol	Explanation
<b>AEC</b>	Extraction Columns Extraktionssäulen Colonne d'extraction Colonne d' estrazione Columnas de extracción
<b>BUF H3PO4</b>	Phosphate Buffer Phosphat-Puffer Tampon phosphate Tampone fosfato Tampón fosfato
<b>Ab</b>	Antiserum Antiserum Antisérum Antisiero Antisuero
<b>TR</b>	Tracer Tracer Traceur Tracciante Trazador
<b>CAL</b>	Calibrator Kalibrator Calibreur Calibratore Calibrador
<b>CONTROL N</b>	Control Normal Normalkontrolle Contrôle normal Controllo normale Control normal
<b>CONTROL H</b>	Control High Kontrolle hoch Contrôle élevé Controllo alto Control alto
<b>Ab2</b>	2 <sup>nd</sup> Antibody 2. Antikörper 2 <sup>ème</sup> Anticorps Secondo anticorpo Segundo anticuerpo

Printing Date  
2013-01-14