# anti-MAG Antibodies ELISA

MAG = mijelin udruženi glikoprotein

za *in vitro* dijagnostičku upotrebu

EK-MAG 96 testova

Datum izdavanja: 2024-07-22  
Verzija A3

# SRPSKI

## NAMENA

Anti-MAG Antibodies ELISA je *in vitro* dijagnostički test za semikvantitativno određivanje anti-MAG IgM antitela u uzorcima humanog seruma. Test je namenjen kao pomoć prilikom postavljanja dijagnoze neuropatije udružene sa anti-MAG antitelima, zajedno sa drugom kliničkim i laboratorijskim nalazima.

Samo za laboratorijsku upotrebu.

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Anti-MAG Antibodies ELISA test omogućava merenje IgM antitela na mijelin udruženi glikoprotein (MAG) u serumu, ELISA testom po principu sendviča. Mikrotitar ploča je obložena prečišćenim MAG humanog mozga. Serumi pacijenata, kontrole i kalibratori se dodaju u udubljenja mikrotitar ploče. Posle 2 sata inkubacije na 2 - 8°C i koraka ispiranja, detekciono antitelo konjugovano sa peroksidazom rena (HRP) detektuje anti-MAG antitela vezana za humani MAG na ploči. Posle još 2 sata inkubacije i daljih koraka ispiranja, dodaje se hromogeni HRP supstrat, tetrametilbenzidin (TMB), (nastaje plava boja), a zatim se reakcija zaustavlja (boja se menja u žutu). Apsorbanca se meri na 450 nm.

Nivo anti-MAG antitela se određuje korišćenjem kalibracione krive generisane na osnovu izmerenih vrednosti kalibratora i izražava se kao BÜHLMANN jedinica titra (BTU).

## OBEZBEĐENI REAGENSI I PRIPREMA

Reagensi	Količina	Kod	Rekonstitucija
<b>Mikrotitar ploča</b> 96 udubljenja obložena humanim MAG	12 stripova sa po 8 udubljenja u ramu	B-MAG-MP	Spremna za upotrebu
<b>Folija za ploču</b>	3 komada	-	Spremna za upotrebu
<b>Koncentrat pufera za ispiranje (10x)</b>	1 boca x 100 mL	B-MAG-WB	Razblažiti sa 900 mL dejonizovane vode
<b>Pufer za inkubaciju</b> Sa konzervansima	1 boca x 100 mL	B-MAG-IB	Spreman za upotrebu
<b>Kalibratori A do D<sup>1</sup></b> Liofilizat sa konzervansima	4 bočice	B-MAG-CASET	Dodati 1 mL Pufera za inkubaciju
<b>Kontrola Niska i Visoka<sup>2</sup></b> Liofilizat sa konzervansima	2 bočice	B-MAG-CONSET	Dodati 1 mL Pufera za inkubaciju

Reagensi	Količina	Kod	Rekonstitucija
<b>IgM obeležen enzimom</b> anti-humano IgM antitelo konjugovano sa HRP-om u puferskom matriksu sa konzervansima	1 bočica x 11 mL	B-MAG-ELM	Spreman za upotrebu Rastvor plave boje
<b>TMB supstrat</b> TMB u citratnom puferu	1 bočica x 11 mL	B-TMB	Spreman za upotrebu
<b>Rastvor za zaustavljanje reakcije</b> 0.25 M sumporna kisleina	1 bočica x 11 mL	B-STC	Spreman za upotrebu <b>Korozivan agens</b>

Tabela 1

<sup>1</sup> Nakon rekonstitucije, kalibratori A, B, C i D sadrže 70000, 15000, 3000 i 1000 BÜHLMANN jedinica titra (BTU) anti-MAG antitela, tim redosledom.

<sup>2</sup> Sadržaj anti-MAG antitela u kontrolama je specifičan za lot. Za tačnu koncentraciju pogledati dodatni dokument kontrole kvaliteta.

## ČUVANJE I ROK UPOTREBE REAGENASA

Zatvoreni / neotvoreni reagensi	
Čuvati na 2-8 °C. Ne koristiti reagense nakon isteka roka upotrebe datog na nalepnicama.	
Otvoreni / Rekonstituisani reagensi	
Mikrotitar ploča	Stripove koji se neće koristiti vratiti odmah u kesicu zajedno sa desikantom i ponovo zatvoriti celom dužinom zip zatvarača. Čuvati do 1 meseca na 2-8 °C.
Razblaženi pufer za ispiranje	Čuvati do 1 meseca na 2-8 °C.
Pufer za inkubaciju	
Enzimski obeležen IgM	
TMB supstrat	
Kontrole	Alikvotirati nakon rekonstitucije i čuvati na ≤-20 °C. Čuvati do 1 meseca na ≤-20 °C. <sup>1</sup>
Kalibratori	
Rastvor za zaustavljanje reakcije	Čuvati do 1 meseca na 18-28 °C.

Tabela 2

<sup>1</sup> Rekonstituisani kalibratori i kontrole mogu se zamrznuti i odmrznuti najviše tri puta tokom 1 meseca.

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Tabela 5

CI – interval pouzdanosti

ROC AUC – area under receiver operating characteristic curve

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
---

## REFERENCE

1. Blirup-Jensen, S., Johnson, A. M. & Larsen, M. Protein standardization V: Value transfer. A practical protocol for the assignment of serum protein values from a Reference Material to a Target Material. *Clin. Chem. Lab. Med.* **46**, 1470–1479 (2008).
2. CLSI guidelines EP30-A - Characterization and Qualification of CommuTabela Reference Materials for Laboratory Medicine (2010).
3. Kuijff, M. L. et al. Detection of anti-MAG antibodies in polyneuropathy associated with IgM monoclonal gammopathy. *Neurology* **73**, 688–695 (2009).
4. Matà, S. et al. IgM monoclonal gammopathy-associated neuropathies with different IgM specificity. *Eur. J. Neurol.* **18**, 1067–1073 (2011).
5. Stork, A. C. J. et al. Classical and lectin complement pathway activity in polyneuropathy associated with IgM monoclonal gammopathy. *J. Neuroimmunol.* **290**, 76–79 (2016).
6. Stork, A. C. J. et al. Fcγ receptor IIIA genotype is associated with rituximab response in antimyelin-associated glycoprotein neuropathy. *J Neurol Neurosurg Psychiatry.* **85**, 916–918 (2014).
7. Liberatore, G. et al. Sensitivity and specificity of a commercial ELISA test for anti-MAG antibodies in patients with neuropathy. *J. Neuroimmunol.* **345**, (2020).
8. Taams, N. E. et al. Clinical relevance of serum antibodies to GD1b in immune-mediated neuropathies. *J. Peripher. Nerv. Syst.* **23**, 227–234 (2018).
9. Campagnolo, M. et al. Polyneuropathy with anti-sulfatide and anti-MAG antibodies: Clinical, neurophysiological, pathological features and response to treatment. *J. Neuroimmunol.* **281**, 1–4 (2015).
10. Vallat, J-M. et al. The Wide Spectrum of Pathophysiologic Mechanisms of Paraproteinemic Neuropathy. *Neurology* **96**, 214-225 (2021).

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