anti-MAG Antibodies ELISA

MAG = mijelin udruženi glikoprotein

za *in vitro* dijagnostičku upotrebu

EK-MAG 96 testova

Datum izdavanja: 2024-07-22
Verzija A3

NAMENA

Anti-MAG Antibodies ELISA je *in vitro* dijagnostički test za semikvantitativno određivanje anti-MAG IgM antitela u uzorcima humanog seruma. Test je namenjen kao pomoć prilikom postavljanja dijagnoze neuropatije udružene sa anti-MAG antitelima, zajedno sa drugom kliničkim i laboratorijskim nalazima.

Samo za laboratorijsku upotrebu.

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Anti-MAG Antibodies ELISA test omogućava merenje IgM antitela na mijelin udruženi glikoprotein (MAG) u serumu, ELISA testom po principu sendviča. Mikrotitar ploča je obložena prečišćenim MAG humanog mozga. Serumi pacijenata, kontrole i kalibratori se dodaju u udubljenja mikrotitar ploče. Posle 2 sata inkubacije na 2 - 8°C i koraka ispiranja, detekciono antitelo konjugovano sa peroksidazom rena (HRP) detektuje anti-MAG antitela vezana za humani MAG na ploči. Posle još 2 sata inkubacije i daljih koraka ispiranja, dodaje se hromogeni HRP supstrat, tetrametilbenzidin (TMB), (nastaje plava boja), a zatim se reakcija zaustavlja (boja se menja u žutu). Apsorbanca se meri na 450 nm.

Nivo anti-MAG antitela se određuje korišćenjem kalibracione krive generisane na osnovu izmerenih vrednosti kalibratora i izražava se kao BÜHLMANN jedinica titra (BTU).

OBEZBEĐENI REAGENSI I PRIPREMA

Reagensi	Količina	Kod	Rekonstitucija
Mikrotitar ploča 96 udubljenja obložena humanim MAG	12 stripova sa po 8 udubljenja u ramu	B-MAG-MP	Spremna za upotrebu
Folija za ploču	3 komada	-	Spremna za upotrebu
Koncentrat pufera za ispiranje (10x)	1 boca x 100 mL	B-MAG-WB	Razblažiti sa 900 mL deionizovane vode
Pufer za inkubaciju Sa konzervansima	1 boca x 100 mL	B-MAG-IB	Spremna za upotrebu
Kalibratori A do D¹ Liofilizat sa konzervansima	4 bočice	B-MAG-CASET	Dodati 1 mL Pufera za inkubaciju
Kontrola Niska i Visoka² Liofilizat sa konzervansima	2 bočice	B-MAG-CONSET	Dodati 1 mL Pufera za inkubaciju

Reagensi	Količina	Kod	Rekonstitucija
IgM obeležen enzimom anti-humano IgM antitelo konjugovano sa HRP-om u puferskom matriksu sa konzervansima	1 bočica x 11 mL	B-MAG-ELM	Spreman za upotrebu Rastvor plave boje
TMB supstrat TMB u citratnom puferu	1 bočica x 11 mL	B-TMB	Spreman za upotrebu
Rastvor za zaustavljanje reakcije 0.25 M sumporna kisleina	1 bočica x 11 mL	B-STC	Spreman za upotrebu Korozivan agens

Tabela 1

¹ Nakon rekonstitucije, kalibratori A, B, C i D sadrže 70000, 15000, 3000 i 1000 BÜHLMANN jedinica titra (BTU) anti-MAG antitela, tim redosledom.

² Sadržaj anti-MAG antitela u kontrolama je specifičan za lot. Za tačnu koncentraciju pogledati dodatni dokument kontrole kvaliteta.

ČUVANJE I ROK UPOTREBE REAGENASA

Zatvoreni / neotvoreni reagensi	
Čuvati na 2-8 °C. Ne koristiti reagense nakon isteka roka upotrebe datog na nalepnicama.	
Otvoreni / Rekonstituisani reagensi	
Mikrotitar ploča	Stripove koji se neće koristiti vratiti odmah u kesicu zajedno sa desikantom i ponovo zatvoriti celom dužinom zip zatvarača. Čuvati do 1 meseca na 2-8 °C.
Razblaženi pufer za ispiranje	Čuvati do 1 meseca na 2-8 °C.
Pufer za inkubaciju	
Enzimski obeležen IgM	
TMB supstrat	
Kontrole	Alikvotirati nakon rekonstitucije i čuvati na ≤-20 °C. Čuvati do 1 meseca na ≤-20 °C. ¹
Kalibratori	
Rastvor za zaustavljanje reakcije	Čuvati do 1 meseca na 18-28 °C.

Tabela 2

¹ Rekonstituisani kalibratori i kontrole mogu se zamrznuti i odmrznuti najviše tri puta tokom 1 meseca.

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Tabela 5

CI – interval pouzdanosti

ROC AUC – area under receiver operating characteristic curve

INTERFERIRAJUĆE SUPSTANCE

Senzitivnost Anti-MAG Antibodies ELISA testa na oralne i injekcione farmaceutske preparate, kao i na endogene supstance procenjena je prema CLSI smernici EP07-A3. Bias rezultata veći od 20% se smatra interferencijom.

Nisu otkrivene interferencije kod sledećih supstanci do navedenih koncentracija: intravenski imunoglobulin (20 mg/mL), cladribine (273 ng/mL), Interferon alfa-2a (49,5 ng/mL), ibuprofen (0,22 mg/mL), reumatoidni faktor (680 IU/mL), hemoglobin (10 mg/mL), hemoliza (10 mg/mL), trigliceridi (20 mg/mL), konjugovani bilirubin (0,4 mg/mL), nekonjugovani bilirubin (0,4 mg/mL).

invalid


TABELE I SLIKE

Primer rezultata

	Nivo [BTU]	Apsorbanca [OD]	Izrač. nivo [BTU]	CV [%]
Slepa proba 1		0,046		
Slepa proba 2		0,049		
Prosek		0,048		
Kalibrator A	70000	2,195	70497	
Kalibrator A	70000	2,188	69508	
Prosek	70000	2,191	70000	0,2
Kalibrator B	15000	1,272	15313	
Kalibrator B	15000	1,245	14693	
Prosek	15000	1,258	15000	1,5
Kalibrator C	3000	0,417	3070	
Kalibrator C	3000	0,400	2931	
Prosek	3000	0,408	3000	2,9
Kalibrator D	1000	0,135	1009	
Kalibrator D	1000	0,132	991	
Prosek	1000	0,134	1000	1,5
Kontrola NISKA		0,360	2602	
Kontrola NISKA		0,376	2731	
Prosek		0,368	2666	3,1
Kontrola VISOKA		1,395	18433	
Kontrola VISOKA		1,383	18090	
Prosek		1,389	18261	0,6
Uzorak 1		0,001	255	
Uzorak 1		0,009	297	
Prosek		0,005	276	116,5
Uzorak 2		1,092	11599	
Uzorak 2		0,969	9511	
Prosek		1,030	10555	8,5

Tabela 6

Primer standardne krive (OD₄₅₀)

Slika 1

Preciznost unutar laboratorije

ID	Srednji nivo, BTU	n	Unutar serije		Između serija		Između dana		Unutar laboratorije	
			SD	%CV	SD	%CV	SD	%CV	SD	%CV
S2	2251	80	267	11,8	199	8,8	130	5,8	357	15,9
S3	8849	80	349	3,9	314	3,6	122	1,4	485	5,5
S4	19683	80	622	3,2	1492	7,6	908	4,6	1855	9,4
S5	37185	80	1684	4,5	3083	8,3	1466	3,9	3806	10,2

Tabela 7

ID	Opis	n	Srednji nivo, BTU	% Unutar kategorije
S1	< 1'000 BTU (negativan)	80	213	98,8

Tabela 8

Reproducibilnost

ID	Srednji nivo, BTU	n	Unutar serije		Između serija		Između dana		Unutar laboratorije	
			SD	%CV	SD	%CV	SD	%CV	SD	%CV
S2	2802	75	181	6,5	517	18,4	261	9,3	606	21,6
S3	9052	75	258	2,9	821	9,1	279	3,1	904	10,0
S4	18241	75	531	2,9	1146	6,3	1475	8,1	1942	10,6
S5	34713	75	893	2,6	2740	7,9	2023	5,8	3521	10,1

Tabela 9

ID	Opis	n	Srednji nivo, BTU	% Unutar kategorije
S1	< 1'000 BTU (negativan)	75	55	100,0

Tabela 10

Ukrštena reaktivnost

Antitela	#	
Antitela na citoplazmu neutrofila (ANCA)	13	
Anti-nukleusna antitela (ANA)	23	
Anti-tireoglobulinska antitela (anti-Tg)	5	
Anti-ribonukleoproteinska antitela	1	
Anti-gangliozidna antitela	anti-GD1b	4
	anti-GM1	3
	anti-GQ1b	5
Anititela na acetilholinske receptore i antitela na tirozin kinazu specifičnu za mišice	7	

Tabela 11

Bolesti sa simptomima koji su slični perifernim neuropatijama (značaj za diferencijalnu dijagnozu)	#
Alkoholizam	1
Dijabetes	5
Amiotrofna lateralna skleroza	14
Šagasova bolest	5
Idiopatska periferna neuropatija	1
Sarkoidoza	4
Waldenstrom-ova bolest	4
Wegenerova granulomatoza	1

Tabela 12

Kliničke performanse

Studija	Pozitivne kontrole	Negativne kontrole	Cut-off	Sens.	Spec.
Kuijf et al., 2009	DPN +IgM MG (n = 68)	OPN + HC (n = 139)	1500 BTU	0,72	0,97
Mata et al., 2011	MGUS PN (n = 46)		3200 BTU	0,37	
Campagnolo et al., 2015	DPN +IgM MG (n = 20)	HDC + HC (n = 3)	1000 BTU	0,94	1,00
Stork et al., 2014	DPN +IgM MG (n = 26)		1000 BTU	0,69	
Stork et al., 2016	DPN +IgM MG (n = 83)	HC (n=83)	1000 BTU	0,59	1,00
Taams et al., 2018	MGUS PN (n = 101)		1500 BTU	0,51	
Liberatore et al., 2020		OPN + HC (n = 222)	7000 BTU		1,00

Tabela 13

DPN+IgM MG, Demijelinizirajuća polineuropatija sa IgM monoklonskom gamopatijom; MGUS PN, Polineuropatija udružena sa monoklonskom gamopatijom nepoznatog značaja; OPN, druga polineuropatija; HC, kontrola poreklom od zdravih osoba; HDC, kontrola poreklom od osoba obolelih od hematoloških bolesti

invalid

Priprema uzoraka / kontrola / kalibratora

Razblažiti uzorke seruma 1:1000 sa (hladnim)! puferom za inkubaciju i dobro promešati vorteksiranjem

Rekonstituisati ili odmrznuti (aliquote) rekonstituisanih kontrola i kalibratora

ostaviti 30 minuta na 2-8 °C

anti-MAG Antibodies ELISA

obložena mikrotitar ploča

isprati 4 x sa $\geq 300 \mu\text{L}$ (hladnog)! pufera za ispiranje

100 μL pufera za inkubaciju, kalibratora, kontrola ili uzoraka seruma (1:1000)

inkubirati 2 sata ($\pm 5 \text{ min}$) na 2-8 °C

isprati 4 x sa $\geq 300 \mu\text{L}$ (hladnog)! pufera za ispiranje

dodati 100 μL enzimskog obeleživača

inkubirati 2 sata ($\pm 5 \text{ min}$) na 2-8 °C

isprati 4 x sa $\geq 300 \mu\text{L}$ (hladnog)! pufera za ispiranje

dodati 100 μL TMB supstrata

inkubirati 30 minuta ($\pm 2 \text{ min}$) na 18-28 °C na šejkeru za ploče ~400-600 rpm

dodati 100 μL rastvora za zaustavljanje reakcije

→ Očitati apsorbancu na 450 nm (unutar 30 minuta)

VREME DOBIJANJA REZULTATA: 5 SATI

REFERENCE

1. Blirup-Jensen, S., Johnson, A. M. & Larsen, M. Protein standardization V: Value transfer. A practical protocol for the assignment of serum protein values from a Reference Material to a Target Material. *Clin. Chem. Lab. Med.* **46**, 1470–1479 (2008).
2. CLSI guidelines EP30-A - Characterization and Qualification of CommuTabela Reference Materials for Laboratory Medicine (2010).
3. Kuijff, M. L. et al. Detection of anti-MAG antibodies in polyneuropathy associated with IgM monoclonal gammopathy. *Neurology* **73**, 688–695 (2009).
4. Matà, S. et al. IgM monoclonal gammopathy-associated neuropathies with different IgM specificity. *Eur. J. Neurol.* **18**, 1067–1073 (2011).
5. Stork, A. C. J. et al. Classical and lectin complement pathway activity in polyneuropathy associated with IgM monoclonal gammopathy. *J. Neuroimmunol.* **290**, 76–79 (2016).
6. Stork, A. C. J. et al. Fcγ receptor IIIA genotype is associated with rituximab response in antimyelin-associated glycoprotein neuropathy. *J Neurol Neurosurg Psychiatry.* **85**, 916–918 (2014).
7. Liberatore, G. et al. Sensitivity and specificity of a commercial ELISA test for anti-MAG antibodies in patients with neuropathy. *J. Neuroimmunol.* **345**, (2020).
8. Taams, N. E. et al. Clinical relevance of serum antibodies to GD1b in immune-mediated neuropathies. *J. Peripher. Nerv. Syst.* **23**, 227–234 (2018).
9. Campagnolo, M. et al. Polyneuropathy with anti-sulfatide and anti-MAG antibodies: Clinical, neurophysiological, pathological features and response to treatment. *J. Neuroimmunol.* **281**, 1–4 (2015).
10. Vallat, J-M. et al. The Wide Spectrum of Pathophysiologic Mechanisms of Paraproteinemic Neuropathy. *Neurology* **96**, 214-225 (2021).

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