

Application Note

Products	BÜHLMANN fCAL® turbo: BÜHLMANN fCAL® turbo Reagent Kit (B-KCAL-RSET) BÜHLMANN fCAL® turbo Control Kit (B-KCAL-CONSET) BÜHLMANN fCAL® turbo Calibrator Kit (B-KCAL-CASET)	CE
Analyzer	Beckman Coulter DxC 500 AU	
Version	20240119	

Before installation, please read the appropriate assay instructions for use. Additionally, refer to the analyzer manual for additional instructions. The reagents supplied are ready to use. Equilibrate reagents at room temperature before loading. Mix gently before loading onto the instrument. Load according to the instrument manual. Use the designated bottles provided by the instrument manufacturer. Avoid bubble formation.

Instrument Settings DxC 500 AU

Reagent Name Fecal Calprotectin REF C34811 (OUS) D02562 (US) DxC 500 AU Settings
Calibrator Name Fecal Calprotectin Calibrator REF C34810

Reagent ID 256

TEST CONFIGURATION & CHEMISTRY DETAILS

Assay Name	Test	Rev	Discipline	Chemistry
Test ID	fCAL		Calculated Result	<input type="checkbox"/>
LIS Code	fCAL		Result Type	Quantitative ▼

UNITS AND RANGE SETTINGS

Use Settings from	None ▼	Units	µg/g ▼	Decimal Places	x.X ▼	Other
Test Kind	General ▼	Revision	01	<input checked="" type="checkbox"/> Multi Reagent Switch		
Reagent Name	fCAL	Reagent ID	256	<input type="checkbox"/> FSE Test		
ABB Name	FCA1G	Parameter Long Name	Fecal Calprotectin FCA1G Fecal			
Region	<input checked="" type="checkbox"/> US	<input checked="" type="checkbox"/> OUS	<input checked="" type="checkbox"/> AP	<input type="checkbox"/> JP	<input checked="" type="checkbox"/> EU	<input type="checkbox"/> Other

GENERAL PARAMETERS

SAMPLE VOLUME	Sample Volume	8.0 µL	Dilution	0 ▼ µL	REACTION OD LIMIT	Low	-2.0000	High	3.0000
	Predilution Rate	1 ▼			REACTION BLANK OD LIMIT	First: Low	-2.0000	High	3.0000
REAGENT VOLUME	R1-1	100 µL	Dilution	0 µL		Last: Low	-2.0000	High	3.0000
	R2-1	20 µL	Dilution	0 µL	ANALYTICAL MEASURING RANGE	Low	30.0	High	2000.0
WAVELENGTH	Primary	540 nm	Secondary	None nm	MANUFACTURER FACTOR	A	1	B	0
METHOD	Fixed 1 ▼				REAGENT ONBOARD STABILITY		90 Days		0 Hours
REACTION SLOPE	+				LIH INFLUENCE CHECK	<input type="checkbox"/> Perform LIH check			
MEASURING POINT	Point 1: First	12	Last	16	Lipemia	+ ▼			
	Point 2: First		Last		Icterus	+ ▼			
Linearity Limit	0 %				Hemolysis	+ ▼			
Lag Time Check	<input type="checkbox"/> Perform Lag Time Check								

Product	BÜHLMANN fCAL® turbo
Analyzer	Beckman Coulter DxC 500 AU
Version	20240119



CALIBRATION PARAMETERS									
Base Unit	Decimal Place	Unit 1	Factor 1	Unit 2	Factor 2	Unit 3	Factor 3	Unit 4	Factor 4
µg/g	1	None	0	None	0	None	0	None	0

CALIBRATOR SPECIFIC

Calibration Type Counts

Formula MB Factor

Calibrator Name Positive Cutoff

SLOPE CHECK Number of Levels

Slope Check

CALIBRATION OD AND CONCENTRATION PARAMETERS

Use highest calibrator for Upper AMR

	Calibrator Name	Conc*	Factor/OD Range Low	Factor/OD Range High
Point 1	fCAL-1		-2.0000	3.0000
Point 2	fCAL-2		-2.0000	3.0000
Point 3	fCAL-3		-2.0000	3.0000
Point 4	fCAL-4		-2.0000	3.0000
Point 5	fCAL-5		-2.0000	3.0000
Point 6	fCAL-6		-2.0000	3.0000
Point 7				

STABILITY AND INTERVAL

Reagent Blank Stability Days Hours Interval

Calibration Stability Days Hours Interval

OD DELTA CHECK

Reagent Blank

Calibration

*Lot Dependent

PROZONE CHECK PARAMETERS

Logic Check 1

Check Points: Point 1 Point 2 Point 3

Decision Values: Value 1 Value 2 Value 3

Limit Points: Limit 1 Limit 2

Check Pattern:

Logic Check 2

Check Points: Point 1 Interval

Decision Values: Value 1 Value 2

Limit Points: Limit 1 Limit 2

Logic Check 3

Check Points: Point 1 Interval

Decision Values: Value 1 Value 2

Limit Points: Limit 1 Limit 2

Product	BÜHLMANN fCAL® turbo
Analyzer	Beckman Coulter DxC 500 AU
Version	20240119



Performance Data

Parameter	Acceptance Criteria	Performance
Method comparison	±15 % bias at clinical decision points of 80 µg/g and 160 µg/g	80 µg/g: -3.3% 160 µg/g: 3.5% (see Table 1)
Precision	≤ 15% CV for samples ≥ 50 µg/g	Total Precision: 4.4 to 10.2% (see Table 2)
Analytical sensitivity	LoB ≤ LoD LoD ≤ LoQ LoQ ≤ 30 µg/g	LoB: 10.5 µg/g LoD: 18.1 µg/g LoQ: 30.0 µg/g
Analytical measuring interval (AMI)		30 to 2000 µg/g
Extended measuring interval (EMI)		30 to 14160.7 µg/g
Linearity	R ² ≥ 0.95 Allowable nonlinearity: samples < 75µg/g: 7.5 µg/g ; sample ≥75µg/g: 10%	R ² = 0.999 29.5 to 14160.7 µg/g
Sample carry-over	Mean carry-over ≤ 0.32% Otherwise a technical precaution must be included in the instrument-specific application note.	No significant sample carry-over
Calibration curve stability	Drift analysis with deviation limit from baseline of ± 20%	32 days
On-board stability	up to 30 days at 2-15°C	up to 90 days at 2-15°C

Table 1 Detailed method comparison performance.

N	Reference range	Passing-Bablok Regression Analysis					Bland-Altman Analysis		
		Slope	Intercept (µg/g)	Bias % at 80 µg/g	Bias % at 160 µg/g	r	Mean bias %	Lower LoA %	Upper LoA %
45	39.1 to 9075.4	1.103	-10.9	-3.3	3.5	0.998	6.3	-22.6	35.1

Table 2 Detailed precision performance

ID	Mean µg/g	Within-run (repeatability)	Between-day	Between-run	Total Precision
		CV%			
P1	53.0	6.6	6.7	3.9	10.2
P2	81.6	5.7	3.3	3.3	7.4
P3	196.6	1.9	8.0	2.5	8.6
P4	800.6	0.9	7.5	2.5	7.9
P5	1786.5	0.8	3.5	2.4	4.4
P6	6840.4	1.2	6.5	3.0	7.3