

## Application Note

<b>Products</b>	BÜHLMANN fPELA® turbo: BÜHLMANN fPELA® turbo Reagent Kit (B-KPELA-RSET) BÜHLMANN fPELA® turbo Control Kit (B-KPELA-CONSET) BÜHLMANN fPELA® turbo Calibrator Kit (B-KPELA-CASET)	<b>CE</b>
<b>Analyzer</b>	Beckman Coulter AU480 / AU680 / AU5800 / DxC 700 AU	
<b>Version</b>	20230420	

Before installation, please read the appropriate assay instructions for use. Additionally, refer to the analyzer manual for additional instructions.

The reagents supplied are ready to use. Equilibrate reagents at room temperature before loading. Mix gently before loading onto the instrument. Load according to the instrument manual. Use the designated bottles provided by the instrument manufacturer. Avoid bubble formation.

## Instrument Settings AU480

Reagent ID: 257

Specific Test Parameters			
General	LIH	ISE	Range
Test Name:	FPE1G ▾	< >	Type: Other-1 ▾ Operation: Yes ▾
Sample Volume	3 μL	Dilution 0 μL	OD Limit
Pre-Dilution Rate	1		Min. OD -2.000 Max. OD 3.000
R1(R1-1)	160 μL	Dilution 0 μL	Reagent OD limit:
			First Low -2.000 High 3.000
			Last Low -2.000 High 3.000
R2 (R2-1)	30 μL	Dilution 0 μL	Dynamic Range Low 9.8 High 500.0
Wavelength:	Pri. 540 nm	Sec. None nm	Correlation Factor A 1 B 0
Method:	END ▾		Factor for Maker A 1 B 0
Reaction slope:	+ ▾		Onboard Stability 91 Day 0 Hour
Measuring Point 1:	First 11	Last 16	LIH Influence Check No ▾
Measuring Point 2:	First	Last	Lipemia + ▾
Linearity:	%		Icterus + ▾
No Lag Time:	No ▾		Hemolysis + ▾

Specific Test Parameters			
General	ISE	Range	
Test Name:	FPE1G ▾	< >	Type: Other-1 ▾
Value/Flag:	# ▾	Level L: #	Level H: #
Specific Ranges:			Panic Value
	From	To	Low High
<input type="checkbox"/> 1.	Sex Year Month	Year Month	Low High
<input type="checkbox"/> 2.	# # #	# # #	# #
<input type="checkbox"/> 3.	# # #	# # #	# #
<input type="checkbox"/> 4.	# # #	# # #	# #
<input type="checkbox"/> 5.	# # #	# # #	# #
<input type="checkbox"/> 6.	# # #	# # #	# #
<input type="checkbox"/> 7.	No demographics		# #
<input type="checkbox"/> 8.	Not within expected values		# #
Unit	μg/g	Decimal Places	1

<b>Product</b>	BÜHLMANN fPELA® turbo
<b>Analyzer</b>	Beckman Coulter AU480 / AU680 / AU5800 / DxC 700 AU
<b>Version</b>	20230420



**Calibration Specific**

General | ISE

Test Name:  < > Type   Use Serum Cal.

Calibration Type:  Formula:  Counts:

<Calibrator Parameters>

	Calibrator*	OD	Conc**	Factor Range		Slope Check
				Low	High	
Point 1:	<input type="text" value="*"/>	<input type="text"/>	<input type="text" value="**"/>	<input type="text" value="-9999999"/>	<input type="text" value="9999999"/>	<input type="text" value="None"/>
Point 2:	<input type="text" value="*"/>	<input type="text"/>	<input type="text" value="**"/>	<input type="text" value="-9999999"/>	<input type="text" value="9999999"/>	<input type="checkbox"/> Reagent Blank
Point 3:	<input type="text" value="*"/>	<input type="text"/>	<input type="text" value="**"/>	<input type="text" value="-9999999"/>	<input type="text" value="9999999"/>	<input type="checkbox"/> Calibration
Point 4:	<input type="text" value="*"/>	<input type="text"/>	<input type="text" value="**"/>	<input type="text" value="-9999999"/>	<input type="text" value="9999999"/>	Advanced Calibration
Point 5:	<input type="text" value="*"/>	<input type="text"/>	<input type="text" value="**"/>	<input type="text" value="-9999999"/>	<input type="text" value="9999999"/>	Operation <input type="text" value="Yes"/>
Point 6:	<input type="text" value="*"/>	<input type="text"/>	<input type="text" value="**"/>	<input type="text" value="-9999999"/>	<input type="text" value="9999999"/>	Interval (RB/ACAL) <input type="text" value="Lot/Lot"/>
Point 7:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 8:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 9:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 10:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	

<Point Cal. For Master Curve> No. of Correction Points  Use Master Curve   Lot Calibration

	Calibrator	OD	Conc	OD Range		Stability
				Low	High	
Point 1:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Reagent Blanks Calibration <input type="text" value="30"/> Day <input type="text"/> Hour
Point 2:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Calibration <input type="text" value="30"/> Day <input type="text"/> Hour

MB Type Factor:  1-Point Calibration Point   With CONC-0

\*user defined \*\*lot specific

<b>Product</b>	BÜHLMANN fPELA® turbo
<b>Analyzer</b>	Beckman Coulter AU480 / AU680 / AU5800 / DxC 700 AU
<b>Version</b>	20230420



## AU680

Reagent ID: 257

Specific Test Parameters	
General	Range
Test Name: <input type="text" value="FPE1G"/> < >	Type: <input type="text" value="Other-1"/> Operation: <input type="text" value="Yes"/>
Sample Volume: <input type="text" value="3"/> $\mu\text{L}$ Dilution: <input type="text" value="0"/> $\mu\text{L}$ OD Limit	Min. OD: <input type="text" value="-2.000"/> Max. OD: <input type="text" value="3.000"/>
Pre-Dilution Rate: <input type="text" value="1"/>	Reagent OD limit: First Low: <input type="text" value="-2.000"/> High: <input type="text" value="3.000"/>
Reagents Volume: R1(R1-1) <input type="text" value="160"/> $\mu\text{L}$ Dilution: <input type="text" value="0"/> $\mu\text{L}$	Last Low: <input type="text" value="-2.000"/> High: <input type="text" value="3.000"/>
R2 Volume: <input type="text" value="30"/> $\mu\text{L}$ Dilution: <input type="text" value="0"/> $\mu\text{L}$	Dynamic Range Low: <input type="text" value="9.8"/> High: <input type="text" value="500.0"/>
Common Reagent Type: <input type="text" value="None"/> Name: <input type="text"/>	Correlation Factor A: <input type="text" value="1"/> B: <input type="text" value="0"/>
Wavelength: Pri. <input type="text" value="540"/> nm Sec. <input type="text" value="None"/> nm	Factor for Maker A: <input type="text" value="1"/> B: <input type="text" value="0"/>
Method: <input type="text" value="END"/>	Onboard Stability: <input type="text" value="91"/> Days <input type="text" value="0"/> Hour
Reaction slope: <input type="text" value="+"/> $\nabla$	LIH Influence Check: <input type="text" value="No"/> $\nabla$
Measuring Point 1: First <input type="text" value="11"/> Last <input type="text" value="16"/>	Lipemia: <input type="text" value="+"/> $\nabla$
Measuring Point 2: First <input type="text"/> Last <input type="text"/>	Icterus: <input type="text" value="+"/> $\nabla$
Linearity: <input type="text"/> %	Hemolysis: <input type="text" value="+"/> $\nabla$
No Lag Time: <input type="text" value="No"/> $\nabla$	

Specific Test Parameters																																																		
General	Range																																																	
Test Name: <input type="text" value="FPE1G"/> < >	Type: <input type="text" value="Other-1"/>																																																	
Value/Flag: <input type="text" value="#"/> $\nabla$ Level L: <input type="text" value="#"/> Level H: <input type="text" value="#"/>																																																		
Specific Ranges:	Panic Value																																																	
<table border="1"> <thead> <tr> <th></th> <th>From</th> <th>To</th> <th>Low</th> <th>High</th> </tr> </thead> <tbody> <tr> <td><input type="checkbox"/> 1. Sex</td> <td><input type="text" value="#"/></td> <td><input type="text" value="#"/></td> <td><input type="text" value="#"/></td> <td><input type="text" value="#"/></td> </tr> <tr> <td><input type="checkbox"/> 2. Year</td> <td><input type="text" value="#"/></td> <td><input type="text" value="#"/></td> <td><input type="text" value="#"/></td> <td><input type="text" value="#"/></td> </tr> <tr> <td><input type="checkbox"/> 3. Month</td> <td><input type="text" value="#"/></td> <td><input type="text" value="#"/></td> <td><input type="text" value="#"/></td> <td><input type="text" value="#"/></td> </tr> <tr> <td><input type="checkbox"/> 4. Year</td> <td><input type="text" value="#"/></td> <td><input type="text" value="#"/></td> <td><input type="text" value="#"/></td> <td><input type="text" value="#"/></td> </tr> <tr> <td><input type="checkbox"/> 5. Month</td> <td><input type="text" value="#"/></td> <td><input type="text" value="#"/></td> <td><input type="text" value="#"/></td> <td><input type="text" value="#"/></td> </tr> <tr> <td><input type="checkbox"/> 6. Year</td> <td><input type="text" value="#"/></td> <td><input type="text" value="#"/></td> <td><input type="text" value="#"/></td> <td><input type="text" value="#"/></td> </tr> <tr> <td><input type="checkbox"/> 7. No demographics</td> <td></td> <td></td> <td><input type="text" value="#"/></td> <td><input type="text" value="#"/></td> </tr> <tr> <td><input type="checkbox"/> 8. Not within expected values</td> <td></td> <td></td> <td><input type="text" value="#"/></td> <td><input type="text" value="#"/></td> </tr> </tbody> </table>		From	To	Low	High	<input type="checkbox"/> 1. Sex	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="checkbox"/> 2. Year	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="checkbox"/> 3. Month	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="checkbox"/> 4. Year	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="checkbox"/> 5. Month	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="checkbox"/> 6. Year	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="checkbox"/> 7. No demographics			<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="checkbox"/> 8. Not within expected values			<input type="text" value="#"/>	<input type="text" value="#"/>	<table border="1"> <thead> <tr> <th>Low</th> <th>High</th> </tr> </thead> <tbody> <tr> <td><input type="text" value="#"/></td> <td><input type="text" value="#"/></td> </tr> </tbody> </table>	Low	High	<input type="text" value="#"/>	<input type="text" value="#"/>
	From	To	Low	High																																														
<input type="checkbox"/> 1. Sex	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>																																														
<input type="checkbox"/> 2. Year	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>																																														
<input type="checkbox"/> 3. Month	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>																																														
<input type="checkbox"/> 4. Year	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>																																														
<input type="checkbox"/> 5. Month	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>																																														
<input type="checkbox"/> 6. Year	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>																																														
<input type="checkbox"/> 7. No demographics			<input type="text" value="#"/>	<input type="text" value="#"/>																																														
<input type="checkbox"/> 8. Not within expected values			<input type="text" value="#"/>	<input type="text" value="#"/>																																														
Low	High																																																	
<input type="text" value="#"/>	<input type="text" value="#"/>																																																	
Unit: <input type="text" value="µg/g"/> Decimal Places: <input type="text" value="1"/>																																																		

Calibration Specific																																																	
General	ISE																																																
Test Name: <input type="text" value="FPE1G"/> < >	Type: <input type="text" value="Other-1"/> <input type="checkbox"/> Use Serum Cal.																																																
Calibration Type: <input type="text" value="6 AB"/> Formula: <input type="text" value="SPLINE"/> Counts: <input type="text" value="2"/>	Slope Check: <input type="text" value="None"/>																																																
<Calibrator Parameters>	Allowable Range Check: <input type="checkbox"/> Reagent Blank <input type="checkbox"/> Calibration <input type="checkbox"/>																																																
<table border="1"> <thead> <tr> <th>Calibrator*</th> <th>OD</th> <th>Conc **</th> <th>Factor Range</th> </tr> <tr> <th></th> <th></th> <th></th> <th>Low High</th> </tr> </thead> <tbody> <tr> <td>Point 1: <input type="text" value="*"/></td> <td><input type="text"/></td> <td><input type="text" value="**"/></td> <td><input type="text" value="-9999999"/> <input type="text" value="9999999"/></td> </tr> <tr> <td>Point 2: <input type="text" value="*"/></td> <td><input type="text"/></td> <td><input type="text" value="**"/></td> <td><input type="text" value="-9999999"/> <input type="text" value="9999999"/></td> </tr> <tr> <td>Point 3: <input type="text" value="*"/></td> <td><input type="text"/></td> <td><input type="text" value="**"/></td> <td><input type="text" value="-9999999"/> <input type="text" value="9999999"/></td> </tr> <tr> <td>Point 4: <input type="text" value="*"/></td> <td><input type="text"/></td> <td><input type="text" value="**"/></td> <td><input type="text" value="-9999999"/> <input type="text" value="9999999"/></td> </tr> <tr> <td>Point 5: <input type="text" value="*"/></td> <td><input type="text"/></td> <td><input type="text" value="**"/></td> <td><input type="text" value="-9999999"/> <input type="text" value="9999999"/></td> </tr> <tr> <td>Point 6: <input type="text" value="*"/></td> <td><input type="text"/></td> <td><input type="text" value="**"/></td> <td><input type="text" value="-9999999"/> <input type="text" value="9999999"/></td> </tr> <tr> <td>Point 7: <input type="text"/></td> <td><input type="text"/></td> <td><input type="text"/></td> <td><input type="text"/></td> </tr> <tr> <td>Point 8: <input type="text"/></td> <td><input type="text"/></td> <td><input type="text"/></td> <td><input type="text"/></td> </tr> <tr> <td>Point 9: <input type="text"/></td> <td><input type="text"/></td> <td><input type="text"/></td> <td><input type="text"/></td> </tr> <tr> <td>Point 10: <input type="text"/></td> <td><input type="text"/></td> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>	Calibrator*	OD	Conc **	Factor Range				Low High	Point 1: <input type="text" value="*"/>	<input type="text"/>	<input type="text" value="**"/>	<input type="text" value="-9999999"/> <input type="text" value="9999999"/>	Point 2: <input type="text" value="*"/>	<input type="text"/>	<input type="text" value="**"/>	<input type="text" value="-9999999"/> <input type="text" value="9999999"/>	Point 3: <input type="text" value="*"/>	<input type="text"/>	<input type="text" value="**"/>	<input type="text" value="-9999999"/> <input type="text" value="9999999"/>	Point 4: <input type="text" value="*"/>	<input type="text"/>	<input type="text" value="**"/>	<input type="text" value="-9999999"/> <input type="text" value="9999999"/>	Point 5: <input type="text" value="*"/>	<input type="text"/>	<input type="text" value="**"/>	<input type="text" value="-9999999"/> <input type="text" value="9999999"/>	Point 6: <input type="text" value="*"/>	<input type="text"/>	<input type="text" value="**"/>	<input type="text" value="-9999999"/> <input type="text" value="9999999"/>	Point 7: <input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Point 8: <input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Point 9: <input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Point 10: <input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Advanced Calibration: <input type="checkbox"/> Operation: <input type="text" value="Yes"/> Interval (RB/ACAL): <input type="text" value="Lot / Lot"/>
Calibrator*	OD	Conc **	Factor Range																																														
			Low High																																														
Point 1: <input type="text" value="*"/>	<input type="text"/>	<input type="text" value="**"/>	<input type="text" value="-9999999"/> <input type="text" value="9999999"/>																																														
Point 2: <input type="text" value="*"/>	<input type="text"/>	<input type="text" value="**"/>	<input type="text" value="-9999999"/> <input type="text" value="9999999"/>																																														
Point 3: <input type="text" value="*"/>	<input type="text"/>	<input type="text" value="**"/>	<input type="text" value="-9999999"/> <input type="text" value="9999999"/>																																														
Point 4: <input type="text" value="*"/>	<input type="text"/>	<input type="text" value="**"/>	<input type="text" value="-9999999"/> <input type="text" value="9999999"/>																																														
Point 5: <input type="text" value="*"/>	<input type="text"/>	<input type="text" value="**"/>	<input type="text" value="-9999999"/> <input type="text" value="9999999"/>																																														
Point 6: <input type="text" value="*"/>	<input type="text"/>	<input type="text" value="**"/>	<input type="text" value="-9999999"/> <input type="text" value="9999999"/>																																														
Point 7: <input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>																																														
Point 8: <input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>																																														
Point 9: <input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>																																														
Point 10: <input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>																																														
<Point Cal. For Master Curve>	Use Master Curve: <input type="checkbox"/> Lot Calibration: <input type="checkbox"/>																																																
<table border="1"> <thead> <tr> <th>Calibrator</th> <th>OD</th> <th>Conc</th> <th>OD Range</th> <th>Stability</th> </tr> <tr> <th></th> <th></th> <th></th> <th>Low High</th> <th>Reagent Blanks Calibration</th> </tr> </thead> <tbody> <tr> <td>Point 1: <input type="text"/></td> <td><input type="text"/></td> <td><input type="text"/></td> <td><input type="text"/></td> <td><input type="text" value="30"/> Day <input type="text" value="0"/> Hour</td> </tr> <tr> <td>Point 2: <input type="text"/></td> <td><input type="text"/></td> <td><input type="text"/></td> <td><input type="text"/></td> <td><input type="text" value="30"/> Day <input type="text" value="0"/> Hour</td> </tr> </tbody> </table>	Calibrator	OD	Conc	OD Range	Stability				Low High	Reagent Blanks Calibration	Point 1: <input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text" value="30"/> Day <input type="text" value="0"/> Hour	Point 2: <input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text" value="30"/> Day <input type="text" value="0"/> Hour																													
Calibrator	OD	Conc	OD Range	Stability																																													
			Low High	Reagent Blanks Calibration																																													
Point 1: <input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text" value="30"/> Day <input type="text" value="0"/> Hour																																													
Point 2: <input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text" value="30"/> Day <input type="text" value="0"/> Hour																																													
MB Type Factor: <input type="text"/> 1-Point Calibration Point: <input type="text"/> $\nabla$ <input type="checkbox"/> With CONC-0																																																	

\*user defined \*\*lot specific

<b>Product</b>	BÜHLMANN fPELA® turbo
<b>Analyzer</b>	Beckman Coulter AU480 / AU680 / AU5800 / DxC 700 AU
<b>Version</b>	20230420



## AU5800

Reagent ID: 257

Parameters		Specific Test Parameters			
General	LIH	ISE	HbA1c	Calculated Test	Range
Test Name:	FPE1G	<	>	Type: Other-1	Operation: Yes
Sample Volume	3 $\mu$ L	Dilution	0 $\mu$ L	OD Limit	
Pre-Dilution Rate	1	Diluent Bottle	#	Min.OD	-2.000
Rgt. Volume	R1(R1-1) 160 $\mu$ L	Dilution	0 $\mu$ L	Max.OD	3.000
	R1-2	Dilution		Reagent OD Limit	
	R2(R2-1) 30 $\mu$ L	Dilution	0 $\mu$ L	1 <sup>st</sup> Low	-2.000
				High	3.000
				Last Low	-2.000
				High	3.000
Common Rgt. Type	None	Name	None	Dynamic Range Low	9.8
Wavelength	540 nm	Sec.	None	High	500.0
Method	END			Correlation Factor A	1
Reaction Slope	+			Factor for Maker A	1
Measuring Point1 1 <sup>st</sup>	11	Last	16	Onboard Stability Period	91 Day
Measuring Point2 1 <sup>st</sup>		Last		LIH Influence Check	No
Linearity Limit	%			Lipemia	+
Lag Time Check	No			Icterus	+
				Hemolysis	+

Parameters		Specific Test Parameters			
General	LIH	ISE	HbA1c	Calculated Test	Range
Test Name:	FPE1G	<	>	Type: Other-1	
Value/Flag:	#			Level	Low # High #
Specific Ranges:		From	To	Low	High
	Sex	Year	Month	Year	Month
<input type="checkbox"/> 1.	#	#	#	#	#
<input type="checkbox"/> 2.	#	#	#	#	#
<input type="checkbox"/> 3.	#	#	#	#	#
<input type="checkbox"/> 4.	#	#	#	#	#
<input type="checkbox"/> 5.	#	#	#	#	#
<input type="checkbox"/> 6.	#	#	#	#	#
7.	Standard demographics				#
8.	Not within expected values				#
Panic Value	Low	#	High	#	Unit: $\mu$ g/g
					Decimal Places: 1

Parameters		Calibration Parameters			
Calibrators	Calibration Specific				
General	ISE				
Test Name:	FPE1G	<	>	Type: Other-1	Cuvette .
		<input type="checkbox"/> Use Serum Cal.			
Calibration Type:	6 AB	Formula:	SPLINE	Counts:	2
<Calibrator Parameters>					
	Calibrator*	OD	Conc**	Range	
Point 1:	*		**	Low High	Slope Check
Point 2:	*		**	-9999999 9999999	
Point 3:	*		**	-9999999 9999999	Allowance Range Check
Point 4:	*		**	-9999999 9999999	<input type="checkbox"/> Reagent Blank
Point 5:	*		**	-9999999 9999999	<input type="checkbox"/> Calibration
Point 6:	*		**	-9999999 9999999	
Point 7:					Advanced Calibration
Point 8:					Operation
Point 9:					Yes
Point 10:					Interval (RB/ACAL)
					Lot/Lot
<Point Cal. For	No. of Correction Points		Use Master Curve		<input type="checkbox"/> Lot Calibration
Master Curve>			OD Range		
	Calibrator	OD	Conc	Low High	Stability
Point-1					Reagent Blank
Point-2					Calibration
					30 Day
					0 Hour
MB Type Factor:	1-Point Calibration Point				<input type="checkbox"/> with Conc-0

\*user defined \*\*lot specific

<b>Product</b>	BÜHLMANN fPELA® turbo
<b>Analyzer</b>	Beckman Coulter AU480 / AU680 / AU5800 / DxC 700 AU
<b>Version</b>	20230420



## DxC 700 AU

Reagent ID: 257

General	LIH	ISE	Calculated Test	Range
<b>Test Name:</b> <span style="border: 1px solid black; padding: 2px;">FPE1G</span> <b>Test No</b> <span style="border: 1px solid black; padding: 2px;"> </span> <b>Type:</b> <span style="border: 1px solid black; padding: 2px;">Other-1</span> <b>Operation</b> <span style="border: 1px solid black; padding: 2px;">Yes</span>				
Sample Volume <input style="width: 40px;" type="text" value="3"/> $\mu\text{L}$	Dilution <input style="width: 40px;" type="text" value="0"/> $\mu\text{L}$	OD Limit	Min. OD <input style="width: 60px;" type="text" value="-2.0000"/>	Max OD <input style="width: 60px;" type="text" value="3.0000"/>
Pre-Dilution Rate <input style="width: 40px;" type="text" value="1"/>			Reagent OD Limit	
Reagent Volume R1 (R1-1) <input style="width: 40px;" type="text" value="160"/> $\mu\text{L}$	Dilution <input style="width: 40px;" type="text" value="90"/> $\mu\text{L}$	1 <sup>st</sup> Low <input style="width: 60px;" type="text" value="-2.0000"/>	High <input style="width: 60px;" type="text" value="3.0000"/>	
R1-2 <input style="width: 40px;" type="text" value=""/> $\mu\text{L}$	Dilution <input style="width: 40px;" type="text" value=""/> $\mu\text{L}$	Last Low <input style="width: 60px;" type="text" value="-2.0000"/>	High <input style="width: 60px;" type="text" value="3.0000"/>	
R2 (R2-1) <input style="width: 40px;" type="text" value="30"/> $\mu\text{L}$	Dilution <input style="width: 40px;" type="text" value="15"/> $\mu\text{L}$	Analytical Measuring Range Low <input style="width: 60px;" type="text" value="9.8"/>	High <input style="width: 60px;" type="text" value="500"/>	
Common Reagent Type <span style="border: 1px solid black; padding: 2px;">None</span>	Name <span style="border: 1px solid black; padding: 2px;">None</span>	Correlation Factor A <input style="width: 40px;" type="text" value="1"/>	B <input style="width: 40px;" type="text" value="0"/>	
Wavelength Pri <input style="width: 40px;" type="text" value="540"/> nm	Sec <span style="border: 1px solid black; padding: 2px;">none</span> nm	Manufacturer Factor A <input style="width: 40px;" type="text" value="1"/>	B <input style="width: 40px;" type="text" value="0"/>	
Method <span style="border: 1px solid black; padding: 2px;">END</span>		Onboard Stability Period <input style="width: 40px;" type="text" value="91"/> Day	<input style="width: 40px;" type="text" value="0"/> Hour	
Reaction Slope <input style="width: 40px;" type="text" value="+"/>		LIH Influence Check <span style="border: 1px solid black; padding: 2px;">No</span>		
Measuring Point-1 1st <input style="width: 40px;" type="text" value="11"/>	Last <input style="width: 40px;" type="text" value="16"/>	Lipemia <input style="width: 40px;" type="text" value="+"/>		
Measuring Point-2 1st <input style="width: 40px;" type="text" value=""/>	Last <input style="width: 40px;" type="text" value=""/>	Icterus <input style="width: 40px;" type="text" value="+"/>		
Linearity Limit <input style="width: 40px;" type="text" value=""/>	%	Hemolysis <input style="width: 40px;" type="text" value="+"/>		
Lag Time Check <input style="width: 40px;" type="text" value=""/>				

General	LIH	ISE	Calculated Test	Range
<b>Test Name:</b> <span style="border: 1px solid black; padding: 2px;">FPE1G</span> <b>Test No</b> <span style="border: 1px solid black; padding: 2px;"> </span> <b>Type:</b> <span style="border: 1px solid black; padding: 2px;">Other-1</span>				
Value/Flag <span style="border: 1px solid black; padding: 2px;">#</span>	Level	Low <input style="width: 60px;" type="text" value="#"/>	High <input style="width: 60px;" type="text" value="#"/>	
<b>Specific Ranges</b>				
	Sex	Year	Month	Year
	From	To	Month	Other Type
<input type="checkbox"/> 1:	<input style="width: 40px;" type="text" value="#"/>	<input style="width: 40px;" type="text" value="#"/>	<input style="width: 40px;" type="text" value="#"/>	<input style="width: 40px;" type="text" value="#"/>
<input type="checkbox"/> 2:	<input style="width: 40px;" type="text" value="#"/>	<input style="width: 40px;" type="text" value="#"/>	<input style="width: 40px;" type="text" value="#"/>	<input style="width: 40px;" type="text" value="#"/>
<input type="checkbox"/> 3:	<input style="width: 40px;" type="text" value="#"/>	<input style="width: 40px;" type="text" value="#"/>	<input style="width: 40px;" type="text" value="#"/>	<input style="width: 40px;" type="text" value="#"/>
<input type="checkbox"/> 4:	<input style="width: 40px;" type="text" value="#"/>	<input style="width: 40px;" type="text" value="#"/>	<input style="width: 40px;" type="text" value="#"/>	<input style="width: 40px;" type="text" value="#"/>
<input type="checkbox"/> 5:	<input style="width: 40px;" type="text" value="#"/>	<input style="width: 40px;" type="text" value="#"/>	<input style="width: 40px;" type="text" value="#"/>	<input style="width: 40px;" type="text" value="#"/>
<input type="checkbox"/> 6:	<input style="width: 40px;" type="text" value="#"/>	<input style="width: 40px;" type="text" value="#"/>	<input style="width: 40px;" type="text" value="#"/>	<input style="width: 40px;" type="text" value="#"/>
7:	Standard demographics			<input style="width: 40px;" type="text" value="#"/>
8:	Not within expected values			<input style="width: 40px;" type="text" value="#"/>
Critical Limits	Low <input style="width: 40px;" type="text" value="#"/>	High <input style="width: 40px;" type="text" value="#"/>	Unit <span style="border: 1px solid black; padding: 2px;">µg/g</span>	Decimal Places <input style="width: 40px;" type="text" value="1"/>

Calibrators	General	ISE	
<b>Test Name:</b> <span style="border: 1px solid black; padding: 2px;">FPE1G</span> <b>Type:</b> <span style="border: 1px solid black; padding: 2px;">Other-1</span>			
<input type="checkbox"/> Use Serum Cal.			
Calibration Type: <span style="border: 1px solid black; padding: 2px;">6AB</span>	Formula: <span style="border: 1px solid black; padding: 2px;">Spline</span>	Counts: <input style="width: 40px;" type="text" value="2"/>	
<b>&lt;Calibrator Parameters&gt;</b>			Slope Check <span style="border: 1px solid black; padding: 2px;">none</span>
	Calibrator*	OD	Conc**
			Range
			Low
			High
<input type="checkbox"/> Point-1	*		**
<input type="checkbox"/> Point-2	*		**
<input type="checkbox"/> Point-3	*		**
<input type="checkbox"/> Point-4	*		**
<input type="checkbox"/> Point-5	*		**
<input type="checkbox"/> Point-6	*		**
<input type="checkbox"/> Point-7	*		**
MB Type Factor <input style="width: 40px;" type="text" value=""/>	1-Point Calibration Point <input style="width: 40px;" type="text" value=""/>	<input type="checkbox"/> with Conc-0	Interval (RB) <span style="border: 1px solid black; padding: 2px;">Lot</span>
			Interval (ACAL) <span style="border: 1px solid black; padding: 2px;">Lot</span>
		Stability	
		Reagent Blank <input style="width: 40px;" type="text" value="30"/> Day <input style="width: 40px;" type="text" value="0"/> Hour	
		Calibration <input style="width: 40px;" type="text" value="30"/> Day <input style="width: 40px;" type="text" value="0"/> Hour	

\*User Defined

\*\*Lot dependent

<b>Product</b>	BÜHLMANN fPELA® turbo
<b>Analyzer</b>	Beckman Coulter AU480 / AU680 / AU5800 / DxC 700 AU
<b>Version</b>	20230420



## Performance Data

Parameter	Acceptance Criteria	Performance
<b>Method comparison</b>	Linear regression fit with slope 0.8-1.2 Bland-Altman: Bias $\pm$ 20% Bias at cutoff: Passing-Bablok: $\pm$ 15%	Slope: 0.995 Mean bias: 2.3% Cutoff bias: 3.1% (see Table 1)
<b>Precision</b>	$\leq$ 10 % CV for samples $\geq$ 25 $\mu\text{g/g}$ and $\leq$ 20 % CV for samples $<$ 25 $\mu\text{g/g}$	Total Precision: 1.1% to 3.8% (see Table 2)
<b>Analytical sensitivity</b>	Goal Limit of Quantification (LoQ) $\leq$ 20 $\mu\text{g/g}$	9.8 $\mu\text{g/g}$
<b>Analytical measuring interval (AMI)</b>		9.8 to 500 $\mu\text{g/g}$
<b>Linearity</b>	$\leq$ 10 $\mu\text{g/g}$ or 10% for fPELA or recovery between 80 and 120% no CV criterion for samples $<$ LoQ	13.9 to 2391.6 $\mu\text{g/g}$
<b>Extended measuring interval (EMI)</b>		9.8 to 2391.6 $\mu\text{g/g}$
<b>High dose hook effect</b>	A sample with concentration of at least 2500 $\mu\text{g/g}$ measures a value $>$ 500 $\mu\text{g/g}$ in the first run.	No HDHE up to 19031 $\mu\text{g/g}$
<b>Sample carry-over</b>	Mean carry-over $<$ 0.75%	No significant sample carry-over
<b>Calibration curve stability</b>	Minimum stability of 14 days	30 days
<b>On-board stability</b>	Minimum stability of 30 days	up to 91 days at 2-15°C

**Table 1 Detailed method comparison performance.**

N	Reference range	Passing-Bablok Regression Analysis				Bland-Altman Analysis		
		Slope	Intercept	% Bias at 200 $\mu\text{g/g}$	r	Mean bias %	Lower LoA %	Upper LoA %
63	15.6 to 3852.0	0.995	7.1	3.1	0.998	2.3	-17.0	21.5

**Table 2 Detailed precision performance**

ID	Mean $\mu\text{g/g}$	Within-run (repeatability)	Between-day	Between-run	Total Precision
					CV%
Low1	81.5	3.2	1.5	1.1	3.8
Mid Low	142.1	2.3	0.6	0.0	2.3
Mid High	266.0	1.1	0.3	0.4	1.2
High1	427.7	0.8	0.7	0.3	1.1
High3	2010.1	3.4	1.7	0.0	3.8