

Application Note

Products	BÜHLMANN fPELA® turbo: BÜHLMANN fPELA® turbo Reagent Kit (B-KPELA-RSET) BÜHLMANN fPELA® turbo Control Kit (B-KPELA-CONSET) BÜHLMANN fPELA® turbo Calibrator Kit (B-KPELA-CASET)	CE
Analyzer	Beckman Coulter DxC 500 AU	
Version	20231219	

Before installation, please read the appropriate assay instructions for use. Additionally, refer to the analyzer manual for additional instructions. The reagents supplied are ready to use. Equilibrate reagents at room temperature before loading. Mix gently before loading onto the instrument. Load according to the instrument manual. Use the designated bottles provided by the instrument manufacturer. Avoid bubble formation.

Instrument Settings DxC 500 AU

Name Fecal Pancreatic Elastase REF C69666 DxC 500 AU Settings Name Fecal Pancreatic Elastase Calibrator REF C69668

Reagent ID 257

TEST CONFIGURATION & CHEMISTRY DETAILS

Assay Name	Test	Rev	Discipline	Chemistry
Test ID	fPELA		Calculated Result	<input type="checkbox"/>
LIS Code	fPELA		Result Type	Quantitative ▼

UNITS AND RANGE SETTINGS

Use Settings from	None ▼	Units	µg/g ▼	Decimal Places	x.x ▼	Other
Test Kind	General ▼	Revision	01	<input checked="" type="checkbox"/> Multi Reagent Switch		
Reagent Name	fPELA	Reagent ID	257	<input type="checkbox"/> FSE Test		
ABB Name	FPE1G	Parameter Long Name	Fecal Pancreatic Elastase C69666 FPE1G Fecal			
Region	<input checked="" type="checkbox"/> US	<input checked="" type="checkbox"/> OUS	<input checked="" type="checkbox"/> AP	<input type="checkbox"/> JP	<input checked="" type="checkbox"/> EU	<input type="checkbox"/> Other

GENERAL PARAMETERS

SAMPLE VOLUME	Sample Volume	3.0 µL	Dilution	0 µL	REACTION OD LIMIT	Low	-2.0000	High	3.0000	
	Predilution Rate	1			REACTION BLANK OD LIMIT	First: Low	-2.0000	High	3.0000	
REAGENT VOLUME	R1-1	160 µL	Dilution	0 µL		Last: Low	-2.0000	High	3.0000	
	R2-1	30 µL	Dilution	0 µL	ANALYTICAL MEASURING RANGE	Low	12.9	High	500.0	
WAVELENGTH	Primary	540 nm	Secondary	800 nm	MANUFACTURER FACTOR	A	1	B	0	
METHOD	END ▼				REAGENT ONBOARD STABILITY		91	Days	0	Hours
REACTION SLOPE	+				LIH INFLUENCE CHECK	<input type="checkbox"/> Perform LIH check				
MEASURING POINT	Point 1: First	11	Last	18	Lipemia	+ ▼				
	Point 2: First		Last		Icterus	+ ▼				
					Hemolysis	+ ▼				
Linearity Limit	0 %									
Lag Time Check			<input type="checkbox"/> Perform Lag Time Check							

Product	BÜHLMANN fPELA® turbo
Analyzer	Beckman Coulter DxC 500 AU
Version	20231219



CALIBRATION PARAMETERS									
Base Unit	Decimal Place	Unit 1	Factor 1	Unit 2	Factor 2	Unit 3	Factor 3	Unit 4	Factor 4
µg/g	1	None	0	None	0	None	0	None	0

CALIBRATOR SPECIFIC

Calibration Type Counts

Formula MB Factor

Calibrator Name Positive Cutoff

SLOPE CHECK Number of Levels

Slope Check

STABILITY AND INTERVAL

Reagent Blank Stability Days Hours Interval

Calibration Stability Days Hours Interval

CALIBRATION OD AND CONCENTRATION PARAMETERS

Use highest calibrator for Upper AMR

	Calibrator Name	Conc*	Factor/OD Range Low	Factor/OD Range High
Point 1	fPELA-1		-2.0000	3.0000
Point 2	fPELA-2		-2.0000	3.0000
Point 3	fPELA-3		-2.0000	3.0000
Point 4	fPELA-4		-2.0000	3.0000
Point 5	fPELA-5		-2.0000	3.0000
Point 6	fPELA-6		-2.0000	3.0000
Point 7				

*Lot dependent

OD DELTA CHECK

Reagent Blank

Calibration

PROZONE CHECK PARAMETERS

Logic Check 1

Check Points
 Point 1
 Point 2
 Point 3

Limit Points
 Limit 1
 Limit 2

Check Pattern
 Pattern

Logic Check 2

Decision Values
 Value 1
 Value 2
 Value 3

Check Points
 Point 1
 Interval

Limit Points
 Limit 1
 Limit 2

Logic Check 3

Decision Values
 Value 1
 Value 2

Check Points
 Point 1
 Interval

Limit Points
 Limit 1
 Limit 2

Decision Values
 Value 1
 Value 2

Product	BÜHLMANN fPELA® turbo
Analyzer	Beckman Coulter DxC 500 AU
Version	20231219



Performance Data

Parameter	Acceptance Criteria	Performance
Method comparison	Slope: 0.8- 1.2 Mean Bias: ≤20% Bias at cutoff: ± 15%	Slope: 1.001 Mean bias: 4.7% Bias at 200 µg/g: 4.0 (see Table 1)
Precision	≤ 10 % for samples ≥ 25 µg/g and ≤ 20 % for samples < 25 µg/g	Total Precision: 1.0% to 4.2% (see Table 2)
Analytical sensitivity	Limit of Quantification (LoQ): ≤ 20 µg/g	12.9 µg/g
Analytical measuring interval (AMI)		12.9 to 500 µg/g
Extended measuring interval (EMI)		12.9 to 2792.7 µg/g
Linearity	R2 ≥ 0.95 deviation from linearity of less than 10 µg/g / 10%	19.3 to 2792.7 µg/g
Sample carry-over	Mean carry-over ≤ 0.75% Otherwise a technical precaution must be included in the instrument-specific application note.	No significant sample carry-over
High Dose Hook Effect	A sample with concentration of 2500 µg/g must return a measured value > 500 µg/g in the first measurement.	No HDHE up to 7756 µg/g
Calibration curve stability	Time interval for re-calibration should be at least 14 days depending on the clinical chemistry analyzer.	41 days
Onboard stability		up to 91 days at 2-15°C

Table 1 Detailed method comparison performance.

N	Reference range	Passing-Bablok Regression Analysis				Bland-Altman Analysis		
		Slope	Intercept	% Bias at 200 µg/g	r	Mean bias %	Lower LoA	Upper LoA
48	7.0 to 2865.4	1.001	7.8	4.0	0.997	4.7	-14.0	23.5

Table 2 Detailed precision performance

ID	Mean µg/g	Within-run (repeatability)	Between- day	Between- run	Total Precision
Low	84.6	3.3	0.0	2.6	4.2
Medium Low	223.4	0.8	0.0	0.6	1.0
Medium High	303.4	1.1	0.3	0.6	1.3
High	481.5	0.5	0.6	0.6	1.0
High3	1630.9	0.9	0.7	0.5	1.2