BÜHLMANN kappa/lambda Reagent Set

**Determination of
kappa and lambda light chains of IgM type
by EK-MAG**

B-KLSET 100 Tests

Revision date: 2015-11-05

ENGLISH

INTENDED USE

The kappa/lambda (κ/λ) Reagent Set is intended as accessory set for BÜHLMANN anti-MAG ELISA, EK-MAG, in order to determine anti-MAG kappa and lambda light chains of IgM type in human serum.

REAGENTS SUPPLIED AND PREPARATION

Reagents	Quantity	Code	Reconstitution
Enzyme Label IgM kappa	1 vial 12 ml	B-MAG-ELMK	Ready to use
Enzyme Label IgM lambda	1 vial 12 ml	B-MAG-ELML	Ready to use

Table 1

STORAGE AND SHELF LIFE OF REAGENTS

Sealed/ Unopened Reagents
Store at 2-8°C. Do not use past expiration date.
Opened reagents
Stable at 2-8°C until expiration date.

Table 2

PRECAUTIONS

SAFETY PRECAUTIONS

- Serum samples may be potentially infectious and should be handled according to Good Laboratory Practice (GLP) using appropriate precautions.
- Clinical specimens may contain infectious agents e.g. Hepatitis B virus and HIV.
- Wear appropriate personal protective equipment to avoid contact with eyes and skin.
- Unused solution should be disposed of according to local State and Federal regulations.

TECHNICAL PRECAUTIONS

- Components must not be used after the expiry date printed on the labels.
- Do not mix different lots of reagents.
- Test performance will be adversely affected if reagents are incorrectly diluted, modified or stored under conditions other than those as detailed in Table 2.
- Avoid contamination of reagents.
- Every effort should be made to ensure that no cross contamination occurs between wells. It is essential that enzyme-conjugated monoclonal antibody is not allowed to contaminate reagents and equipment. Care should be taken to avoid touching or splashing the rim of the well. Conjugate allowed to dry on to the rim of the well may adversely affect the performance of the test.
- All equipment including washers and liquid handling devices must be free of microbial contamination and be correctly calibrated and maintained according to the manufacturer's instructions.
- Microwells cannot be re-used.
- Protect substrate from light.
- It is recommended that samples are run in duplicate until the laboratory is totally familiar with the assay.

- Test results should be interpreted in conjunction with information available from epidemiological studies, clinical assessment of the patient, and other diagnostic procedures.

REAGENTS REQUIRED BUT NOT PROVIDED

Reagents	Supplier	Code	Processing
Anti-MAG ELISA Kit	BÜHLMANN Laboratories AG	EK-MAG	Refer to instruction for use

Table 3

MATERIALS REQUIRED BUT NOT PROVIDED

- Precision pipettes with disposable tips: 2 μ l, 100 μ l and 1000 μ l pipettes.
- Disposable polystyrene or polypropylene tubes for the preparation of sample dilutions.
- 1000 ml cylinder for the dilution of the Wash Buffer Concentrate.
- Microtiter plate washer or squeeze bottle for Wash Buffer.
- Microtiter plate rotator.
- Microtiter plate reader for measurement of absorbance at 450 nm.

SPECIMEN COLLECTION AND STORAGE

- The procedure calls for <0.1 ml of blood or <50 μ l of serum.
- Lipemic, hemolytic and icteric samples should not be used in this assay. Lipemic samples can be avoided by asking patients to fast for at least 12 hours prior to the sample being taken.
- Collect blood into plain tubes, avoid hemolysis, leave to clot for one hour at RT (18-28°C), centrifuge for 15 minutes at approximately 1800 x g at RT and collect the serum.
- Store serum samples at \leq -20°C. Samples are stable for \geq 1 year if stored at \leq -20°C.
- Avoid repeated freeze-thaw cycles. Frozen samples should be thawed and mixed thoroughly by gentle swirling or inversion prior to use.

ASSAY PROCEDURE

1. Dilute all patient samples 1:1000 with cold Incubation Buffer (B-MAG-IB) (e.g. 2 μ l of serum + 2 ml of Incubation Buffer). Allow diluted samples to set for one hour at 18-28°C vortex from time to time. Put samples for 10 minutes on ice prior to pipetting in step 4.
2. Prepare sufficient strips (wells) of a MAG microtiter plate (B-MAG-MP) in order to test the desired number of samples and controls. Reseal excess strips in the foil pouch together with the desiccant packs **without delay**. Store them refrigerated.

Note: Use refrigerated reagent solutions in steps 3. to 9.

3. Wash the coated wells four times using at least 300 μ l of refrigerated Wash Buffer (B-MAG-WB) per well. Empty the wells and strike the plate firmly onto blotting paper.

- 4a. Pipet 100 µl of the High Control (B-MAG-CONSET) in quadruplicate into wells A1+A2 and B1 + B2.
 - 4b. Pipet 100 µl of each diluted sample in quadruplicate into the subsequent wells.
 5. Cover the plate with a Plate Sealer and incubate for 2 hours (± 5 min) at 2-8°C.
 6. Remove and discard the Plate Sealer. Empty the wells and wash four times using at least 300 µl of **refrigerated** Wash Buffer per well. Empty the wells and strike plate firmly onto blotting paper.
 - 7a. Add 100 µl of Enzyme Label IgM kappa to the first row (A1, B1, C1 ...) of wells for each sample/control and
 - 7b. add 100 µl of Enzyme Label IgM lambda to the second row of wells (A2, B2, C2 ...) for each sample/control
 8. Cover the plate with plate sealer and incubate for 2 hours (± 5 min) at 2-8°C.
 9. Remove and discard plate sealer. Empty the wells and wash four times using at least 300 µl of **refrigerated** Wash Buffer per well. Empty the wells and strike the plate firmly onto blotting paper.
- Note: Allow TMB substrate solution to reach 18-28°C.**
10. Add 100 µl of TMB Substrate Solution (B-MAG-TMB) to each well.
 11. Cover plate with a Plate Sealer, place it on a plate mixer set at 800-1000 rpm, protect plate from direct light and incubate for 30 minutes (± 5 min) at 18-28°C.
 12. Add 100 µl of Stop Solution (B-MAG-STS) to all wells. Remove air bubbles with a pipette tip. Proceed to step 13. within 30 minutes.
 13. Read absorbance at 450 nm in a Microtiter Plate reader.


Figure 1: Quitard A et al; 2011. Determination of subtypes IgM kappa and IgM lambda in patients with an anti-MAG antibody neuropathy; (2011).

RESULTS

Calculation and Cut-off

- Record the absorbance at 450 nm of each diluted patient sample.
- Average duplicate values.
- For calculation compare OD[kappa] with OD[lambda]

$$Ratio = \frac{OD[kappa]}{OD[lambda]}$$

The cut-off for monoclonality is set as follows (see Figure 1):

- kappa –light chain: Ratio > 4
- lambda –light chain: Ratio < 0.25

QUALITY CONTROL

A thorough understanding of this package insert is necessary for the successful use of the product. Refer also the information given in the package insert of the anti-MAG ELISA (order code EK-MAG).

Reliable results will be obtained only by using precise laboratory techniques (current GLP guidelines) and accurately following this package insert.

Since there is no control serum for anti-MAG antibodies commercially available, we recommend using a positive and negative serum pool for internal quality controls.

All controls must fall within established confidence limits.

The isotype of the control sample is printed on the additional data sheet.

PERFORMANCE LIMITATIONS

HAEMOLYSIS

Grossly haemolysed samples should not be used as lower assay values may be obtained.

LIPEMIC Serum

Samples should be taken from fasting individuals.

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