Quantum Blue[®] sCAL

Quantitative Lateral Flow Assay

LF-MRP25

25 tests

Revision date: 2016-11-23

ENGLISH

INTENDED USE

The Quantum Blue® sCAL lateral flow assay is an immunoassay designed for the quantitative determination of MRP8/14 in human serum in combination with the Quantum Blue® Reader.

MRP8/14 serum levels can be used to determine the inflammatory status of patients with chronic inflammatory diseases such as Rheumatoid Arthritis (RA) e.g. to monitor the treatment efficiency of drugs such as anti-TNF alpha antibodies.

For professional use only.

*Canada, Taiwan: For laboratory use only.

PRINCIPLE OF THE ASSAY

The test is designed for the selective measurement of MRP8/14 antigen by sandwich immunoassay. A monoclonal capture antibody which is highly specific for MRP8/14 is lined onto the membrane. A second monoclonal detection antibody conjugated to gold colloids is deposited onto the conjugate release pad and released into the reaction system after addition of the diluted serum. The antibody-gold-conjugate reacts with the MRP8/14 of the sample and then binds to the anti-MRP8/14 antibody which is lined onto the membrane (Test Line). The remaining free MRP8/14-gold-conjugate binds to the goat anti-mouse antibody which is lined onto the membrane (Control Line). The signal intensities of the Test Line and the Control Line, respectively, are measured quantitatively by the Quantum Blue® Reader.

REAGENTS SUPPLIED AND PREPARATION

Reagents	Quantity	Code	Comments
Test Cassette	25 pieces	B-LFCALUS-TC	In vacuum-sealed foil bag pouch
Chase Buffer	1 bottle 10 mL	B-LFMRP-CB	Ready to use
Controls Low* / High*	2 vials 0.5 mL/vials	B-LFMRP- CONSET	Ready to use
RFID Chip Card	1 piece	B-LFMRP-RCC	White plastic card

Table 1

* The controls contain lot specific amounts of MRP8/14. Refer to the additional QC data sheet for actual concentrations.

STORAGE AND SHELF LIFE OF REAGENTS

All kit components are stable at 2-8 °C until the expiration date printed on the labels.

PRECAUTIONS

Safety Precautions

- None of the reagents of this test contains components of human origin.
- Patient specimens should be handled as if capable of transmitting infections and should be handled in accordance with good laboratory practice (GLP) using appropriate precautions.

- Unused solution should be discarded according to local State and Federal regulations.

Technical Precautions

Kit Components

- All reagents and test samples must be equilibrated at room temperature (18-28 °C) before starting the assay. Mix (vortex) the reagents well before use.
- The chase buffer might show flocculation which does not have any influence on the quality of the test and the test result.
- Components must not be used after the expiry date printed on the labels.
- Do not mix different lots of reagents.
- Test cassettes cannot be re-used.

Assay Procedure

- Read carefully the instructions prior to carrying out the assay. Assay performance will be adversely affected, if reagents are incorrectly diluted, handled or stored under conditions other than those as detailed in this instruction for use.
- Use the white RFID chip card in order to change lot-specific test parameters.
- The Quantum Blue® Reader must be switched on and programmed for the Quantum Blue® sCAL assay before starting the assay (see Quantum Blue® Reader manual).
- Patient samples that are not properly handled may cause inaccurate results.
- Diluted samples should be used within several hours and cannot be stored for a longer time period.
- Undiluted serum samples can be stored at ≤ -20 °C for at least 2 years.

MATERIALS REQUIRED BUT NOT PROVIDED

- Quantum Blue® Reader available at BÜHLMANN (order code: BI-POCTR-ABS).
- Timer
- Precision pipettes with disposable tips: 1-10 μ L and 100-1000 μ L
- Centrifuge for micro tubes
- Disposable polystyrene or polypropylene tubes for dilution of samples
- Gloves and laboratory coat

SPECIMEN COLLECTION AND STORAGE

Collect blood into plain venipuncture tubes without any additives, avoid hemolysis, and let the serum clot at room temperature (18-28 °C) for at least 20 and for up to 60 minutes. Centrifuge at room temperature and at $\sim 2'000$ g for 15 minutes. Decant the serum. Serum samples can be stored refrigerated at 2-8 °C for up to 1 month. For longer storage, keep serum samples at ≤ -20 °C. The samples are stable for at least 2 years at ≤ -20 °C.

ASSAY PROCEDURE

1. Dilution of Serum Samples

- Prior to the measurement, dilute serum samples 1:10 with chase buffer (B-LFMRP-CB) (e.g. mix 20 µL sample with 180 µL chase buffer) in a test tube.
- Mix well (e.g. using a vortex mixer).

2. Lateral Flow Assay Procedure and Readout

- Load the lot specific parameters from the RFID chip card.
- Add 60 µL of diluted sample onto the sample loading port of the Test Cassette.
- Incubate for 12 minutes (set a timer manually).
- Load the Test Cassette onto the Test Cassette holder of the reader.
- Scan the Test Cassette with the Quantum Blue® Reader by pressing the start (“ENTER”) button immediately.
- For Low / High Controls: Repeat step 2 using 60 µL of Control instead of diluted serum.

Remark: Please refer to the Quantum Blue® Reader manual to learn about the basic functions and how to initialize and operate the reader, especially how to select test methods, and how to load lot-specific parameters from the RFID chip card in order to get the samples measured.

QUALITY CONTROL

- If the precision of the assay does not correlate with the established limits and repetition, exclude errors in technique, and check the following issues: i) pipets, thermometers and timers, ii) expiration date of reagents, iii) storage and incubation conditions.
- Result of the self-test of the Quantum Blue® Reader performed at startup of the instrument has to be valid.

VALIDATION OF RESULTS

- For a valid test result, the Control Line (C) must be visible (see appendix I, Figures 1A and 1B). It is used as functional test control only and cannot be used for the interpretation of the Test Line (T). If the Test Line (T) is not detectable after the incubation (Figure 1A), the concentration of MRP8/14 present in the sample is below the detection limit. If a Test Line (T) is detectable after the incubation (Figure 1B), the MRP8/14 concentration in the sample is calculated by the Quantum Blue® Reader.
- If only the Test Line (T) is detectable after the incubation (Figure 1C), the test result is invalid and the MRP8/14 assay has to be repeated using another Test Cassette.
- If neither the Control Line (C) nor the Test Line (T) are detectable after the incubation (Figure 1D), the test result is invalid and the assay has to be repeated using another Test Cassette.
- As the Quantum Blue® Reader allows for a quantitative evaluation of the test (T) and control (C) lines, an additional validity check of the Control Line (C) is undertaken. If the signal intensity of the Control Line (C) is below a lot-specific threshold after the incubation

time, the test result is also invalid and the MRP8/14 assay has to be repeated using another Test Cassette.

STANDARDIZATION AND INTERPRETATION OF RESULTS

- The lateral flow assay is calibrated against the BÜHLMANN sCAL ELISA (order code: EK-MRP8/14).
- The Quantum Blue® Reader uses a lot-specific standard curve to calculate the MRP8/14 concentration. The assay range is defined being between 0.5 and 10 µg/mL.
- For quantitative measurements of samples reading above 10 µg/mL dilute the serum samples additionally 1:10 (total 1:100) and test it again according to the assay procedure. The measured concentration must then be multiplied by the reciprocal volume factor to obtain the final result.

LIMITATIONS

- The reagents supplied with this kit are optimized to measure human MRP8/14 in serum samples.
- MRP8/14 results can be used supplementary to other information in order to establish the diagnosis.

RESULTS & INTERPRETATION

The Quantum Blue® sCAL test helps quickly establish a first estimate of the acute inflammatory status of the patient.

The following reference values might help characterize the inflammatory status of a patient. Patients with RA* (n = 20) show mean levels of serum MRP8/14 of 6.25 µg/mL (SD: 3.43 µg/mL) whereas patients responding to treatment (n = 20) showed mean levels of 2.92 µg/mL (SD: 1.45 µg/mL) (ref 1,2).

MRP8/14 normal values were determined using serum samples from apparently healthy and asymptomatic blood donors (adult men and woman at the age 18-70 years). 122 serum samples were analysed by EK-MRP8/14 due to the good correlation ($R^2 = 0.94$, slope = 0.93) with the LF-MRP25.

The MRP8/14 levels of patients with RA were significantly higher compared to healthy controls (n = 122) with a median serum MRP8/14 concentration of 1.16 µg/mL and the 95th percentile at a concentration of 2.92 µg/mL.

*RA (Rheumatoid arthritis)

PERFORMANCE CHARACTERISTICS

Limit of Blank (LoB): 0.25 µg/mL MRP8/14

Limit of Detection (LoD): 0.42 µg/mL MRP8/14

Repeatability: 17.2 % CV. The repeatability of the Quantum Blue® sCAL assay was calculated from 10 samples within 20 days. Each sample was tested in duplicates twice a day according to the assay procedure. The repeatability varied between 11.9 and 25.4 %.

Linearity: 0.5 to 10 µg/mL. Four serum samples with elevated MRP8/14 concentrations were diluted with chase buffer, B-LFMRP-CB. Each dilution was subsequently assayed according to the assay procedure. The results showed linearity within the indicated measuring range of 0.5 to 10 µg/mL of the Quantum Blue® sCAL assay.

High Dose Hook Effect: A high dose hook effect was not observed up to a concentration of 146 µg/mL MRP8/14.

Serum Indices: No interference is detected with the following substances up to the listed concentrations: Triglycerides (Intralipid® 1320 mg/dL; equivalent to 37 mmol/L triglyceride), conjugated bilirubin (342 µmol/L; 29 mg/dL), unconjugated bilirubin (342 µmol/L; 20 mg/dL) or haemoglobin (2 g/L (= 200 mg/dL)).

Method Comparison: $R^2 = 0.94$; $y = 0.93x + 0.51$ µg/mL

29 samples within the indicated measuring range of the Quantum Blue® sCAL assay were analyzed 10 x each according to the assay procedure and compared with the values of 1-4 x duplicates obtained with the BÜHLMANN sCAL ELISA. The correlation data are illustrated in Figure 3.

Limit of Quantification (LoQ):

Lower LoQ: 0.42 µg/mL MRP8/14

Upper LoQ: >10 µg/mL MRP8/14

The LoQ has been established with fifteen samples. The samples were measured with 10 replicates each. The resulting precision profile is shown in Figure 2. The limit of quantification corresponds to the concentration of MRP8/14 with an imprecision below 25 % CV allowing a quantitative measurement within the range from 0.42 µg/mL (lower LoQ) to >10 µg/mL (upper LoQ).

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