

CALEX[®] Cap



Fast

Easy

**Stable
fecal
extraction**


- Maximum workflow optimisation
- Setting the standard for easy fecal extraction
- Stable at room temperature

CALEX® Cap Unique Features

The prefilled tubes are ready to use in three simple steps:

- 1** Dip the dosing tip into the fecal sample and fill the grooves completely with the matrix.
- 2** Place the pin back into the tube, close tightly and extract by soaking and vortexing.
- 3** For turbidimetric assays only: Centrifuge the extract within the CALEX® tube for 10 min.

The resulting extract (1:500) is ready to use in the turbidimetric assay. The CALEX® extraction is highly correlating to the manual weighting method.



Pre-analytical side
Total laboratory automation and smooth rack filling with new sampling pinhead

Fully quantitative results

Analytical side
Cap ergonomics facilitates automated and manual decapping


- SIMPLICITY** Cap ergonomometry supports the simple and safe use of CALEX® Cap, significantly reducing hands-on time
- AUTOMATION** Designed for total laboratory automation (TLA) systems
- STABILITY** Extensive validation proves stability at RT for 7 days. The CALEX® Cap is suitable for air transportation according to IATA 650 (UN3373)
- EFFICIENCY** Optimal extract dilution of 1:500



From the Patient ...



...in turbo speed...



...to the Lab



BÜHLMANN Laboratories AG
Baselstrasse 55
4124 Schönenbuch
Switzerland

Phone +41 61 487 12 12
Fax orders +41 61 487 12 99
info@buhlmannlabs.ch
www.buhlmannlabs.ch



Ordering Codes:
B-CALEX-C50 50 devices
B-CALEX-C200 200 devices
B-CALEX-C500 500 devices

CE-marked products

CALEX® is a registered trademark of BÜHLMANN in many countries.
Parts of the CALEX® Cap are patent protected by:
EP2833795(B1); US9752967(B2); AU2016203121(B2); CA2997598(C);
JP6307132(B2); KR10-1875862(B1).