

# Quantum Blue<sup>®</sup> Infliximab trough level rapid test

Test Rapido  
per modulare  
tempestivamente  
la terapia




- Il primo test rapido per il monitoraggio dell'IFX
- Risultati in 15 minuti
- Range su siero 0,4 – 20 µg/mL
- Range esteso per il dosaggio d'induzione
- Marchio CE/IVD validato anche per bio-similari

# Quantum Blue® Infliximab test rapido

## Caratteristiche Principali

- Test rapido quantitativo per determinare la concentrazione di Infliximab nel siero
- Range di misura 0,4 - 20 µg/mL, estensibile per la fase d'induzione
- Copre il range tra 3 e 7 µg/mL per la definizione del trough level
- Utilizzato negli studi clinici in IBD e Artrite Reumatoide
- Validato per Remsima®, Inflectra® and Flixabi®
- Ottima correlazione con i test ELISA di riferimento in letteratura

## Trough levels Infliximab



\* Vande Casteele et al. 2015



BÜHLMANN Laboratories AG  
Switzerland  
info@buhlmannlabs.ch  
www.buhlmannlabs.ch  
Phone +41 61 487 12 12

BÜHLMANN Italia Srl  
Italy  
info@buhlmannlabs.it  
www.buhlmannlabs.it  
Fax ordini +41 61 487 12 99

Quantum Blue® is a registered trademark  
of BÜHLMANN in many countries.

**Codici ordini:**  
LF-TLIF10 10 card  
LF-TLIF25 25 card

**Prodotti correlati:**  
LF-TLAD10 10 card  
LF-TLAD25 25 card  
LF-CAL25 25 card  
LF-CALE25 25 card  
LF-CHR25 25 card

CE-marked products