

Quantum Blue® Infliximab

trough level rapid test

Test Rapido
per modulare
tempestivamente
la terapia

- Il Test Rapido Quantitativo leader per il monitoraggio dell'IFX
- Risultati in 15 minuti
- Range su Siero 0,4 – 20 µg/mL
- Range esteso per il dosaggio d'induzione
- Marchio CE validato anche per bio-similari



Quantum Blue® Infiximab trough level rapid test

Caratteristiche principali

- Test rapido quantitativo per determinare la concentrazione d'Infiximab
- Alta correlazione allo standard WHO per infiximab¹ (NIBSC 16/170) e a test ELISA accreditati di terza parte
- Range 0,4 -20 µg/mL, estesibile a 183,6 µg/mL per la fase d'induzione
- Clinicamente Validato^{2,3}
- Validato per Remsima®, Inflectra®, Zessly® e Flixabi®

Infiximab "Trough levels"

<3 µg/mL

Livello d'infiximab
sub-terapeutico
terapia da
adattare⁴

3-7 µg/mL

Dose Efficace
d'infiximab⁴

>7 µg/mL

Livello d'infiximab
sopra-terapeutico
terapia da
adattare⁴

Bibliografia:

¹Keller et al., 2020, First successful comparison of Quantum Blue® rapid TDM assay standardization with WHO international standard for infiximab

²Parra, S. et al., 2018, Infiximab Trough Levels and Quality of Life in Patients with Inflammatory Bowel Disease in Maintenance Therapy, *Gastroenterol Res and Pract*

³Magro, F. et al., 2017, Clinical performance of an infiximab rapid quantification assay, *Ther Adv Gastroenterol*

⁴Castele, N. V. et al., 2015, Trough Concentrations of Infiximab Guide Dosing for Patients With Inflammatory Bowel Disease, *Gastroenterology* (TAXIT trial)



BÜHLMANN Laboratories AG
Switzerland
info@buhlmannlabs.ch
www.buhlmannlabs.ch
Phone +41 61 487 12 12

BÜHLMANN Italia Srl
Italy
info@buhlmannlabs.it
www.buhlmannlabs.it
ordini: buhlmannitalia.ordini@legalmail.it

Quantum Blue® is a registered trademark
of BÜHLMANN in many countries.

Codici ordini:
LF-TLIF10 10 card
LF-TLIF25 25 card

Prodotti correlati:
LF-TLAD10 10 card
LF-TLAD25 25 card
LF-ADIF10 10 card
LF-ADIF25 25 card
LF-CALE25 25 card
LF-CHR25 25 card
CE-marked products