


CALcheck Blue™

CALPROTECTIN

qualitative
rapid test
easy to use

Differentiation between organic
and functional intestinal diseases



- Ideal Rapid test strip for IBD/IBS screen
- Cut off at 50 µg/g
- Results within 10 minutes
- CE mark

- No pipetting
- Rapid and simple test: reaction time 10 minutes
- Easy to read: one band one answer
- Control Band corresponds to cut off level for easy interpretation
- CALcheck Blue™ provides highly reproducible and consistent results
- CALcheck Blue™ has high correlation to quantitative Calprotectin ELISA and Quantum Blue® tests
- Overall agreement to ELISA: 97.6%; PPV: 100%; NPV: 96.6%

